


Surgical Din Baskets



Exceptionally robust and resistant to accidental damage

A comprehensive range of surgical instrument trays and baskets that provide a high degree of protection for storage, transportation and sterilising of costly instruments. Available in a variety of sizes & designs that accommodate a wide range of Surgical Instruments.

**Available
in different
Sizes & Designs**



Popular Sizes
Code: **JM/CR-4728**
240x120x50mm
Code: **JM/CR-4613**
240x250x50mm

Din Baskets

- Side Perforated
- Cross-hatched flat base mesh for durability
- Retracting handles
- Lids sold separately

- **Made from high quality Stainless Steel for durability and strength**
- **Burr free finish**
- **Electrolyte Polished**
- **Baskets are stackable and are designed with rounded corners and edges to avoid common problems of wraps being torn by sharp edges eliminating the possibility of user injury**
- **Top and bottom frame wire diameter 6mm**
- **Wire mesh with opening square 4.2mm (The wire diameter is 0.9mm)**


Fine Mesh Boxes

- Hinged removable Lids



Endoscope Baskets

- Removable lid with special twist fastener, protects against accidental opening
- Soft silicone clamping element on a stainless steel support for optimum fixing



Surgical Instrument Trays & Baskets

Din Basket - Side Perforated, Woven Base

CODE	SIZE
JM/MS-4728	240x120x50mm
JM/MS-4613	240x250x50mm
JM/MS-4615	240x250x70mm
JM/MS-4622	300x300x50mm
JM/MS-4623	350x350x60mm
JM/MS-4618	480x250x50mm
JM/MS-4619	480x250x60mm
JM/MS-4620	480x250x70mm
JM/MS-4625	480x480x50mm
JM/MS-4654	480x480x60mm
JM/MS-4628	500x300x50mm
JM/MS-4634	600x300x60mm

Din Basket - Side Perforated Flat Base

CODE	SIZE
JM/CR-4728	240x120x50mm
JM/CR-4613	240x250x50mm
JM/CR-4615	240x250x70mm
JM/CR-4622	300x300x50mm
JM/CR-4623	350x350x60mm
JM/CR-4618	480x250x50mm
JM/CR-4619	480x250x60mm
JM/CR-4620	480x250x70mm
JM/CR-4625	480x480x50mm
JM/CR-4654	480x480x60mm
JM/CR-4628	500x300x50mm
JM/CR-4634	600x300x60mm

POPULAR SIZES

Din Basket - Full Wire Mesh

CODE	SIZE
JM/MS-4645	240x120x50mm
JM/MS-4646	240x250x50mm
JM/MS-4648	240x250x70mm
JM/MS-4656	300x300x50mm
JM/MS-4657	350x350x60mm
JM/MS-4649	480x250x50mm
JM/MS-4650	480x250x60mm
JM/MS-4651	480x250x70mm
JM/MS-4653	480x480x50mm
JM/MS-4654	480x480x60mm
JM/MS-4657	500x300x50mm
JM/MS-4658	600x300x60mm

Din Basket Lids

CODE	SIZE
JM/MS-4635	240x120mm
JM/MS-4636	240x250mm
JM/MS-4639	300x300mm
JM/MS-4640	350x350mm
JM/MS-4637	480x250mm
JM/MS-4638	480x480mm
JM/MS-4641	500x300mm
JM/MS-4642	600x300mm

Endoscope Basket - Single Scope

CODE	SIZE
JM/MS-4601	290x80x52mm
JM/MS-4604	460x80x52mm
JM/MS-4607	670x80x52mm

Endoscope Basket - Double Scope

CODE	SIZE
JM/MS-4603	290x80x52mm
JM/MS-4605	460x80x52mm
JM/MS-4607	670x80x52mm
JM/MS-4609	640x160x52mm

Fine Mesh Box With Lid

CODE	SIZE
JM/CR-4678	80x30mm
JM/CR-4730	80x80x40mm
JM/CR-4736	100x100x40mm
JM/CR-4670	120x155x50mm
JM/CR-4676	178x127x25mm
JM/CR-4677	178x229x25mm

Perforated Trinket Box With Lid

CODE	SIZE
JM/MS-4678	80x30mm
JM/MS-4730	80x80x40mm
JM/MS-4736	100x100x40mm
JM/MS-4670	120x155x50mm
JM/MS-4676	178x127x25mm
JM/MS-4677	178x229x25mm

ACCESSORIES

CODE	SIZE
JM/CR-4723	Holding Pins



Perforated Trinket Boxes

V1:2016

01527 559010

