


# C-VAC TORNADO<sup>™</sup> Smoke Evacuation Electrosurgical Pencil




CIMPAX<sup>®</sup> C-VAC TORNADO<sup>™</sup>  
Electrosurgical Pencil with Smoke Evacuation  
LATEX, LEAD, PVC and PHTHALATE FREE

# CIMPAX<sup>®</sup>

CE<sub>0543</sub>

CIMPAX<sup>®</sup> C-VAC TORNADO™ the ultimate combination of electrosurgery and removal of hazardous smoke in one instrument. The pencil has a unique triangular design which combined with the silicone cover offers the Surgeon perfect ergonomic and safety grip. The ultra soft micro switch push buttons gives further comfort for the user due the exceptional tactile feedback. Since the C-VAC TORNADO™ is Latex, lead, PVC and phthalate free you can perform surgery without concerns for your health, your colleagues and patient's health or the environment.

- Unique triangular design from tip to end.
- Ultra soft push buttons for enhanced comfort.
- Maximum safety grip silicone cover.
- Transparent tip for clear view of the electrode.
- LATEX, LEAD, PVC and PHTHALATE FREE.
- Flexible tubing with integrated ESU cable.
- Ready to use, no assembling of suction tip or electrode required.



#### Ordering information :-

**C-VAC TORNADO™** Integrated smoke evacuation pencil  
4 meter lead, PVC and phthalate free ESU cable with  
US-3-pin connector and 3 meter PVC and phthalate free  
tubing with combined 22/15/9mm connector. Single use.

With 70 mm standard blade electrode  
With 70 mm coated blade electrode  
Packaging: Box of 10 pcs., sterile.  
Extender, box of 25 pcs., sterile.

Ref. C31-19-001  
Ref. C31-19-101  
Ref. C31-15-145

**Judd**  
medical

Distributed in the UK by:

Judd Medical Limited  
Unity House, Buntsford Park Road, Bromsgrove,  
Worcestershire, B60 3DX  
Tel: 01527 559010 Fax: 01527 559080  
Email: sales@judd-medical.co.uk Web: www.judd-medical.co.uk