

FOOTSWITCHES

This latest range of Prima "Easyclean" electrosurgery footswitches is produced in a special latex free, synthetic rubber - easy to clean!

Non-slip and water tight, the 'Easyclean' range of footswitches are designed to resist the most caustic cleaning agents used in hospitals today.

Safety features

Unique serial number for tracing



Safety features

Individually labelled to identify product

Single & Twin Footswitches

- Easy to clean
- Non slip
- Latex free
- Resists caustic cleaning agents



All cables 5m long

SINGLE Footswitches

Yellow (cut)	Blue (coag)	Black	White

Generator connectors

3pin Valleylab, Megadyne Bipolar					
Erbe Bipolar					
Eschmann 300 & 400 Series		FS109	FS117		FS719
Eschmann 800 Series		FS209	FS217		
Unterminated footswitches without connectors					FS176

TWIN Footswitches



All cables 5m long

Generator connectors

4pin Conmed /Aspen					FS250
4pin Valleylab /Megadyne					FS206
Erbe					FS720
Eschmann 300 & 400 Series					FS207
Eschmann 800 Series					FS307