

UNIVERSAL BIPOLAR CABLE

Reusable

Judd
medical

- Fits ALL US & European Bipolar Instruments
- Silicone sheathed
- Takes twin pin or block connector
- 100 uses

Grey

V-30

Bowa, EMC, Erbe Int.
Ethicon, Lamidey,
Valleylab Coa comp,
Olympus,
Soering

V-30s
4m cable



Blue

S-30/B

Orange

S-30/O

Fuchsia

S-30/F

Grey

S-30

Twin 4mm pin,
Erbe Int.
Eschmann,
Alsa MCD series

S-30s
4m cable



UNIVERSAL MONOPOLAR CABLE

Reusable

- Double insulated cable
- Built-in cable restraints
- Highly flexible cable/connector
- 100 uses

Orange

B-50/O

C-50/O

Fuchsia

B-50/F

C-50/F

Purple

B-50/P

C-50/P

Blue


B-50/B

C-50/B

B
4mm



C
8mm



Generator connectors



UNIVERSAL 2:1 MONOPOLAR CABLE

Reusable

- Always have the correct cable
- Reduce stock level requirements
- Economic, efficient, environmentally effective



4mm pin



8mm pin



Upto 300 autoclave cycles

**A truly 'Universal'
Monopolar Cable fits all
instrument connector pins**

Instrument connector - Fits all male
pins 4-5mm



Product Code

Description

Pk
Size

PMS8700

Universal Monopolar 2:1, 4m Cable

10