


BIPOLAR FORCEPS

Single Use

Judd
medical

Features & Benefits

Highly polished tip to minimise sticking



Integral Cable

Pre-attached
3m Cable

Standard Bipolar Forceps c/w Integral Cable

Product Code	Description	Tip size	Length	Single Use	Pk Size
S310-110-050/AP	McPherson disposable forceps	0.5mm	102mm	•	10
S310-100-050/AP	Jeweller #5 disposable forceps	0.5mm	115mm	•	10
S310-190-100/AP	Adson disposable forceps	1.0mm	120mm	•	10
S310-120-050/AP	Straight disposable forceps	0.5mm	150mm	•	10
S310-120-100/AP	Straight disposable forceps	1.0mm	150mm	•	10
S310-160-100/AP	Straight disposable forceps	1.0mm	180mm	•	10
S310-140-100/AP	Straight disposable forceps	1.0mm	200mm	•	10
S310-170-100/AP	Straight angled tip disposable forceps	1.0mm	200mm	•	10
S310-210-100/AP	Bayonet disposable forceps	1.0mm	165mm	•	10
S310-300-100/AP	Bayonet disposable forceps	1.0mm	200mm	•	10
S310-310-100/AP	Bayonet disposable forceps	1.0mm	220mm	•	10
S310-200-100/AP	Bayonet disposable forceps	1.0mm	240mm	•	10

BIPOLAR FORCEPS Non-Stick

Single Use

Features & Benefits

Safe precise grasping and manipulation of tissue

Flexible bridge connector to allow for greater alignment, with 2 x 4mm plugs.
Suitable for the most popular diathermy machines

CE marked

Gold Coloured Tips

Titanium Nitride Coating,
Non-Stick finish



Integral Cable

Pre-attached
3m Cable

Product Code	Description	Tip size	Length	Single Use	Pk Size
JNS-5110M	Non-Stick McPherson Bipolar Forcep	0.5mm	110mm	•	10
JNS-5115J	Non-Stick Jeweller #5 Bipolar Forcep	0.5mm	115mm	•	10
JNS-1120A	Non-Stick Adson Bipolar Forcep	1.0mm	120mm	•	10
JNS-5150S	Non-Stick Straight Bipolar Forcep	0.5mm	150mm	•	10
JNS-1150S	Non-Stick Straight Bipolar Forcep	1.0mm	150mm	•	10
JNS-1160S	Non-Stick Straight Bipolar Forcep	1.0mm	160mm	•	10
JNS-1200S	Non-Stick Straight Bipolar Forcep	1.0mm	200mm	•	10
JNS-1200SA	Non-Stick Straight, Angled Bipolar Forcep	1.0mm	200mm	•	10
JNS-1165B	Non-Stick Bayonet Bipolar Forcep	1.0mm	165mm	•	10
JNS-1200B	Non-Stick Bayonet Bipolar Forcep	1.0mm	200mm	•	10
JNS-1220B	Non-Stick Bayonet Bipolar Forcep	1.0mm	220mm	•	10
JNS-1240B	Non-Stick Bayonet Bipolar Forcep	1.0mm	240mm	•	10