

Faster, cleaner and safer procedures during laparoscopic surgery.

- Lower power settings achievable
- Reintroduction during surgery is minimised

Insulation tested

Meets requirements of EN60601 European safety standards



Clearer vision

Advanced PTFE coating for reduced eschar build-up, less smoke & clearer vision



4mm banana plug fitting offers secure connection to most monopolar cables



nonstick Laparoscopic Electrodes

Product Code	Description	Length	Single Use	Pk Size
LES104/G	nonstick Laparoscopic Hook electrode	330mm	•	10
LEL104/G	nonstick Laparoscopic Hook electrode	450mm	•	10
LES184/G	nonstick Laparoscopic Spatula electrode	330mm	•	10
LEL184/G	nonstick Laparoscopic Spatula electrode	450mm	•	10
LES204/G	nonstick Laparoscopic L-Hook electrode	330mm	•	10
LEL204/G	nonstick Laparoscopic L-Hook electrode	450mm	•	10
LES104HC/G	nonstick Laparoscopic L-wire electrode, with hand control	330mm	•	10
LES184HC/G	nonstick Laparoscopic Spatula electrode, with hand control	330mm	•	10
LES204HC/G	nonstick Laparoscopic L-Hook electrode, with hand control	330mm	•	10