

Cotswold microdissection needles offer improved cutting and coagulation performance, whilst minimising thermal spread, offering the surgeon the precision required for exceptionally fine dissection in a wide variety of procedures.

- Ultra sharp tip for clean, precise soft tissue dissection.
- Heat resistant tungsten alloy maintains sharpness
- Effective at low power settings
- Highly polished tip to facilitate cleaning

**Greater accuracy**


The fine tungsten tip allows for greater accuracy to be achieved at very low power settings



**Precise electrosurgery**


Needle tip remains sharp throughout the procedure. With a melting point greater than 3420°C, the needle will not become blunt, especially at the lower power settings

Product Code	Description	Shaft	Tip	Single Use	Pk Size
<b>PMS 1650</b>	Cotswold needle, 10mm straight needle	30mm	3mm	•	10
<b>PMS 1651</b>	Cotswold needle, 20mm straight needle	30mm	3mm	•	10
<b>PMS 1614</b>	Cotswold needle, 20mm straight needle	40mm	3mm	•	10
<b>PMS 1654</b>	Cotswold needle, 10mm straight needle	55mm	3mm	•	10
<b>PMS 1602</b>	Cotswold needle, 10mm straight needle	65mm	3mm	•	10
<b>PMS 1656</b>	Cotswold needle, 10mm straight needle	75mm	3mm	•	10
<b>PMS 1613</b>	Cotswold needle, 10mm straight needle	100mm	3mm	•	10
<b>PMS 1612</b>	Cotswold needle, 10mm straight needle	130mm	3mm	•	5
<b>PMS 1611</b>	Cotswold needle, 10mm straight needle	150mm	3mm	•	5
<b>PMS 1610</b>	Cotswold needle, 10mm straight needle	180mm	3mm	•	5




**5mm 45° Angled**

Product Code	Description	Shaft	Tip	Angled	Single Use	Pk Size
<b>PMS 1692</b>	Cotswold needle, 5mm 45° angled needle	30mm	3mm	45°	•	10
<b>PMS 1693</b>	Cotswold needle, 5mm 45° angled needle	40mm	3mm	45°	•	10
<b>PMS 1694</b>	Cotswold needle, 5mm 45° angled needle	55mm	3mm	45°	•	10
<b>PMS 1695</b>	Cotswold needle, 5mm 45° angled needle	75mm	3mm	45°	•	10
<b>PMS 1696</b>	Cotswold needle, 5mm 45° angled needle	100mm	3mm	45°	•	10
<b>PMS 1697</b>	Cotswold needle, 5mm 45° angled needle	130mm	3mm	45°	•	10




**10mm 45° Angled**

Product Code	Description	Shaft	Tip	Angled	Single Use	Pk Size
<b>PMS 1653</b>	Cotswold needle, 10mm 45° angled needle	30mm	3mm	45°	•	10
<b>PMS 1631</b>	Cotswold needle, 10mm 45° angled needle	40mm	3mm	45°	•	10
<b>PMS 1655</b>	Cotswold needle, 10mm 45° angled needle	55mm	3mm	45°	•	10
<b>PMS 1603</b>	Cotswold needle, 10mm 45° angled needle	65mm	3mm	45°	•	10
<b>PMS 1605</b>	Cotswold needle, 10mm 45° angled needle	75mm	3mm	45°	•	10
<b>PMS 1623</b>	Cotswold needle, 10mm 45° angled needle	100mm	3mm	45°	•	10



**10mm 90° Angled**

Product Code	Description	Shaft	Tip	Angled	Single Use	Pk Size
<b>PMS 1660</b>	Cotswold needle, 10mm 90° angled needle	30mm	3mm	90°	•	10
<b>PMS 1661</b>	Cotswold needle, 10mm 90° angled needle	40mm	3mm	90°	•	10
<b>PMS 1662</b>	Cotswold needle, 10mm 90° angled needle	55mm	3mm	90°	•	10



Please note: Some none stock items may have minimum order quantity requirements.