

BIPOLAR CABLE

Single Use

Product Code	Description	Pk Size
PMS3701/AP	Bipolar Cable, Single Use, twin 4mm Pin to USA, 3m Cable, fly leads	10
PMS3703/AP	Bipolar Cable, Single Use, twin 4mm pin to Euro, 3m Cable, fly leads	10
PMS3708/AP	Bipolar Cable, Single Use, twin 4mm pin to Euro, 3m Cable, Valleylab connector	10
PMS3709/AP	Bipolar Cable, Single Use, twin 4mm pin to US twin pin connector, 3m cable, Valleylab connector	10



MONOPOLAR CABLE

Single Use

Product Code	Description	Pk Size
PMS3601/AP	Monopolar Cable, Single Use, 4mm pin to 4.8mm U.K forcep socket, 3m cable	10
PMS3602/AP	Monopolar Cable, Single Use, 8mm pin to 4.8mm U.K forcep socket, 3m cable	10
PMS3603/AP	Monopolar Cable, Single Use, 4mm pin to 4mm endoscopic instrument socket, 3m cable	10
PMS3604/AP	Monopolar Cable, Single Use, 8mm pin to 4mm endoscopic instrument socket, 3m cable	10

