


This updated range of **SINGLE USE Bipolar Forceps** designed to be used in soft tissue surgical procedures, has **GOLD COLOURED TIPS, Titanium Nitride Coating** which allows the heat to be efficiently diffused around the tips entire surface eliminating hot spots and promoting a **NON-STICK** effect.



# NEW Single Use Non-Stick Bipolar Forceps

Available in 0.5mm & 1mm Tips



## Features & Benefits:

- ✓ Titanium Nitride Coating, Non-stick finish
- ✓ Safe precise grasping and manipulation of tissue
- ✓ Proven functionality as with reusable forceps, replicating balance & design
- ✓ Flexible bridge connector to allow for greater alignment, with 2 x 4mm plugs & pre-attached 3m cable. Suitable for the most popular diathermy machines.
- ✓ Variety of patterns, tip sizes and lengths to suit a range of surgical applications
- ✓ All products are CE marked

### SINGLE USE NON-STICK BIPOLAR FORCEPS Sold in Pks of 10

Product Code	Description	Tip size	Length
<b>JNS-5110M</b>	<b>McPherson</b> Bipolar Forcep	0.5mm	110mm
<b>JNS-5115J</b>	<b>Jeweller #5</b> Bipolar Forcep	0.5mm	115mm
<b>JNS-1120A</b>	<b>Adson</b> Bipolar Forcep	1.0mm	120mm
<b>JNS-5150S</b>	Straight Bipolar Forcep	0.5mm	150mm
<b>JNS-1150S</b>	Straight Bipolar Forcep	1.0mm	150mm
<b>JNS-1160S</b>	Straight Bipolar Forcep	1.0mm	160mm
<b>JNS-1200S</b>	Straight Bipolar Forcep	1.0mm	200mm
<b>JNS-1200SA</b>	Straight, Angled Bipolar Forcep	1.0mm	200mm
<b>JNS-1165B</b>	Bayonet Bipolar Forcep	1.0mm	165mm
<b>JNS-1200B</b>	Bayonet Bipolar Forcep	1.0mm	200mm
<b>JNS-1220B</b>	Bayonet Bipolar Forcep	1.0mm	220mm
<b>JNS-1240B</b>	Bayonet Bipolar Forcep	1.0mm	240mm