

# Single Use Spackmans Intra Uterine Cannula

We've **ADDED** a **NEW**  
**Single Use Instrument**  
to our Product Range

This **SINGLE USE Spackmans Cannula** allows for the intra-operative manipulation of both the anteverted and retroverted uterus as well as facilitating quick, simple, effective hysterosalpingogram for tubal patency studies.

- **Cost effective single use version of the well established, but difficult to clean, reusable Spackman Cannula.**
- **Eliminates the significant cleaning and re-processing costs of a frequently used, fine lumen instrument. Ensuring a clean, sterile instrument every time.**
- **Anatomical cervical cone ensures effective, leak free seal.**
- **Atraumatic tip provides reduced uterine perforation risk.**
- **Convenient, adjustable tenaculum attachment.**
- **Weighted handle provides effective 'hands free' anteversion.**
- **CE Marked**



**Product Codes:**

**JMSCA32/12** 32cm x 12fg (Sold in Pks of 10)

**JMSCA37/12** 37cm x 12fg (Sold in Pks of 10)

