


# AUTOClabels

## Autoclave tracking system



Ideal for tracking all:

- Reusable cables
- Larger instruments
- Electrosurgery devices
- Each label has its own unique consecutive number
- Ideal for quick identification



The AUTOClabels punch & label are easy to attach to the product, no tools needed, just clip on.

No tie needed - just push and lock the punch on the cable and the label just feeds back into itself - takes seconds!


Pop a hole in the AUTOClabels every time the product is autoclaved with the permanently attached AUTOClabels punch

OR

Track your product's autoclave history by using the unique label number printed at the top of each AUTOClabels tag

OR

Use the blank AUTOClabels to apply your own barcode



### AUTOClabels

Product Code	Description	Autoclave cycles	Reusable	Pk Size
<b>ACL30</b>	AUTOClabel, Punch and label set	30	•	100
<b>ACL40</b>	AUTOClabel, Punch and label set	40	•	100
<b>ACL60</b>	AUTOClabel, Punch and label set	60	•	100
<b>ACL100</b>	AUTOClabel, Punch and label set	100	•	100
<b>ACL000</b>	AUTOClabel, Blank for barcode	100	•	100

