







DISPERSIVE CABLES



| Dispersive plate connector |
|---|
| Fits standard 25mm wide dispersive plate |
|  |

| Generator connectors | | |
|--|---|---------------------------|
| Valleylab (old generation) |  | DE1553 3m cable |
| | | DE1555 5m cable |
| Erbe ACC/ACC Erbe T Series Erbe VIO |  | DE1523 3m cable |
| | | DE1525 5m cable |
| Alsa / Bard / Berchtold Bovie-Aaron 950, 1200 Eschmann TD400, TD800 Martin / Meditom |  | DE1503 3m cable |
| | | DE1505 5m cable |
| GU Eschmann 300 series |  | DE1513 3m cable |
| | | DE1515 5m cable |
| Valleylab SSE2 K, L and SSE3 Conmed |  | DE1533 3m cable |
| | | DE1535 5m cable |
| Aspen / Bowa / Valleylab Erbe to 2.5 split-pad Conmed / Petas 500, 600 Bovie-Aaron 1250,2250,3250 |  | DE1543 3m cable |
| | | DE1545 5m cable |

