


# XODUS MEDICAL

## Gel Positioning Products

### Effective Solutions for Positioning



Pittsburgh, PA U.S.A • Tel: 724-337-5500  
Fax: 724-337-0555

[www.xodusmedical.com](http://www.xodusmedical.com)  
[info@xodusmedical.com](mailto:info@xodusmedical.com)



ISO 13485 CERTIFIED  
FDA REGISTERED

MADE IN THE  
**USA**

**41806**  
Adult Horseshoe



**41805**  
Adult Head Donut



**41808**  
Pedi Horseshoe



**41807**  
Pedi Head Donut



**41831**  
Candy Cane



**41834**  
Prone Head Rest



**41835**  
Supine Head Rest



**41822**  
Heel Protector



**41830**  
Chest Roll

**41825**  
Chest Roll

*Chest Rolls Available in Six Sizes*

**xodus**  
medical, inc.  
*Making Surgery Safer™*

Tel: 1-724-337-5500  
Fax: 1-724-337-0555

[www.xodusmedical.com](http://www.xodusmedical.com)  
[email: info@xodusmedical.com](mailto:info@xodusmedical.com)

CE 0499



**41819**  
Oval Ulnar Nerve Pad



**41820**  
Square Ulnar Nerve Pad

**41809**  
Hip Pad



**41823**  
Sand Bag



*Sand Bags  
Available in Two Sizes*

**41813**  
Arm Board Pads

*Available in Three Sizes*



**41811**  
Medium Knee Crutch



*Knee Crutches Available in  
Three Sizes*


**41801**  
3 Piece Table Set  
*Contains 41802, 41803 and 41804*

**41816**  
Full Table Pad  
*3/4 table pads  
available with or  
without cutout*

**41802**  
Head Section

**41803**  
Torso Section

**41804**  
Foot Section



- **Soft, Durable and Comfortable**
- **Protects Sensitive Skin Tissue from Pressure Sores**
- **Reusable – Can Be Cleaned by Standard Operating Room Disinfectant**
- **Bar Coded to Facilitate Inventory Management**
- **Wide Range of Shapes and Sizes**
- **Effective for Maintaining Patient Position**



## Gel Pad Positioning Products

<i>Reorder Number</i>	<i>Product</i>	<i>Size</i>	<i>Case Quantity</i>
41801	3 Piece Table Set	Includes 41802, 41803 and 41804	Set
41802	Head Section	1/2" X 9" X 20"	Each
41803	Torso Section	1/2" X 39 1/2" X 20"	Each
41804	Foot Section	1/2" X 20" X 20" w/cutout	Each
41834	Prone Head Rest	11" L X 9" W X 6" H	Each
41835	Supine Head Rest	8" L X 9" W X 4" H	Each
41805	Adult Head Donut	2" X 3" I.D. X 8" O.D.	Each
41806	Adult Horseshoe	1 3/4" X 3" I.D. X 8" O.D.	Each
41807	Pedi Head Donut	1 1/4" X 2 1/4" I.D. X 5 1/2" O.D.	Each
41808	Pedi Horseshoe	1 1/4" X 2 1/4" I.D. X 5 1/2" O.D.	Each
41809	Hip Pad	1/2" X 20" X 20"	Each
41810	Small Knee Crutch	3/8" X 9" X 11"	Each
41811	Medium Knee Crutch	1/4" X 11" X 17"	Each
41812	Large Knee Crutch	1/2" X 20" X 20"	Each
41813	Small Armboard	1/4" X 4" X 15"	Each
41814	Large Armboard	1/4" X 5" X 25"	Each
41815	Duraboard Pad	3/8" X 6" X 26"	Each
41816	Full Table Pad	1/2" X 20" X 72"	Each
41817	3/4 Table Pad w/	1/2" X 20" X 46"	Each
41818	3/4 Table Pad w/o	1/2" X 20" X 46"	Each
41819	Oval Ulnar Nerve Pad	1/4" X 6" X 18"	Each
41820	Square Ulnar Nerve Pad	1/4" X 12" X 18"	Each
41821	Foot/Heel Protector	1/4" X 4" X 15"	Each
41822	Heel Protector	Clam Heel Protector	Each
41823	Small Sandbag	2" X 4 1/2" X 10"	Each
41824	Large Sandbag	2 1/4" X 5" X 14"	Each
41825	Chest Roll	12" L X 4" W X 4 H"	Each
41826	Chest Roll	12" L X 6" W X 6 H"	Each
41827	Chest Roll	16" L X 4" W X 4 H"	Each
41828	Chest Roll	16" L X 6" W X 6 H"	Each
41829	Chest Roll	20" L X 4" W X 4 H"	Each
41830	Chest Roll	20" L X 6" W X 6 H"	Each
41831	Candy Cane	2" O.D. X 24"	Each



Pittsburgh, PA U.S.A • Tel: 724-337-5500  
Fax: 724-337-0555

[www.xodusmedical.com](http://www.xodusmedical.com)  
[info@xodusmedical.com](mailto:info@xodusmedical.com)



ISO 13485 CERTIFIED  
FDA REGISTERED

