

Judd
medical

Ear, Nose & Throat Instrumentation Catalogue

www.judd-medical.co.uk

REUSABLE

We are an Approved
Supplier to the
NHS Supply Chain

NHS


Supply Chain

Index




	Pages		
Artery Forceps	110-111		
Aural Dressing Forceps	111-112		
Aural Forceps	113-118		
Aural Instruments	118-127		
Aural Snares	127		
Aural Specula	128-131		
Biopsy Forceps	131-132		
Bipods	132	Scoops	161
Bone Instruments	133	Suction Tubes - Aspiration	161
Bone Nibbling Forceps	133	Suction Tubes	162-165
Bronchoscopes	134	Syringes	165
Cannulae	135	Tongue Depressors	166
Cerumen Hooks	135	Tonsil Forceps	166
Cocanising Wires	136	Tonsil Holding Forceps	166-167
Conchotomes	136	Tonsil Snares	167
Curettes	136-138	Tracheal Dilators	167
Dissectors	138-139	Trocars	167-168
Elevators	139	Turbinate Forceps	168-169
Gouges	139	Wool Carriers	169
Knives	149	Misc	169
Laryngeal Forceps	140-141		
Laryngoscopes	141-142		
Mirrors	142		
Mouth Gags	143-145		
Mouth Jacks	145		
Myringotomes	146		
Nasal Forceps	146-147		
Nasal Instruments	147-148		
Nasal Rasps	150		
Nasal Snares	150		
Nasal Specula	151-153		
Nodule Forceps	153		
Perforators	154		
Punch Forceps	154-157		
Retractors	157-158		
Scissors - Nasal/Aural	158-160		

Artery Forceps




G052-014 Birkett Artery Forcep, curved, 18cm

G052-016 Birkett Artery Forcep, curved to side, 18cm




G052-012 Birkett Artery Forcep, straight, 18cm



G052-237 Negus Artery Forcep, large curve, 18cm


G052-233 Negus Artery Forcep, small curve, 18cm



G052-340 Sawtell Tonsil Artery Forcep, curved, 18cm

Artery Forceps

G052-350 Schnidt Tonsil Artery Forcep, curved, 19cm



G052-420 Wilson Tonsil Artery Forcep, curved, 19cm



Aural Dressing Forceps


E038-305 Pritchard Aural Forcep, slender jaw, square serrated ends, effective length 5cm, 14.5cm



E040-305 Hunter Todd Aural Forcep, bulbous serrated jaw, 12.5cm



Aural Dressing Forceps




E038-105 Hartman Aural Forcep, serrated jaw, heavy pattern, screw joint, 14.5cm




E040-205 Lucae(Fagge) Aural Forcep, bayonet, serrated jaw, 14cm

E040-210 Lucae Aural Forcep, bayonet, 1:2 teeth, 14cm



E038-405 Tilley Aural Forcep, serrated jaw, 6.7cm working length, 14.5cm overall length



E040-515 Wilde (Troeltsch) Dressing Forcep, 1:2 teeth, 12cm




E040-505 Wilde (Troeltsch) Dressing Forcep, serrated ends, 12cm

E040-520 Wilde (Troeltsch) Dressing Forcep, serrated ends, 15cm


E040-525 Wilde (Troeltsch) Dressing Forcep, serrated ends, 18cm

- E036-105** Beales Aural Forcep, crocodile action, triangular serrated jaw, straight, 4.0mm x 2.5mm




- E036-205** Cawthorne Aural Forcep, crocodile action, smooth straight jaw, 4.0mm x 0.8mm

- E046-050** Cawthorne Micro Aural Forcep, smooth jaw, straight, 3.0mm x 0.5mm



- E036-305** Fuller Aural Forcep, crocodile action, straight jaw, piston holding, 4.0mm x 0.8mm




- E036-420** Hartman Aural Forcep, crocodile action, extra fine serrated jaw, 4.0mm x 0.8mm


- E036-430** Hartman Aural Forcep, crocodile action, extra fine serrated jaw, curved upwards, 4.0mm x 0.8mm

- E036-432** Hartman Aural Forcep, crocodile action, extra fine serrated jaw, curved to left, 4.0mm x 0.8mm


- E036-434** Hartman Aural Forcep, crocodile action, extra fine serrated jaw, curved to right, 4.0mm x 0.8mm




Aural Forceps




E036-425 Hartman Aural Forcep, crocodile action, tenaculum jaw, 7.0mm x 1.25mm



E036-405 Hartman Aural Forcep, crocodile action, grunwald punch, 7.0mm x 1.25mm




E036-410 Hartman Aural Forcep, crocodile action, round cup jaw, 7.0mm x 2.25mm




E036-415 Hartman Aural Forcep, crocodile action, standard serrated jaw, 7.0mm x 1.25mm, 7.5cm working length


- E038-205** Henckel (Strumpel) Aural Forcep, crocodile action, fine fenestrated blade, effective length 7.5cm




- E046-310** Hough Micro Aural Forcep, serrated jaw, straight, 4.0mm x 0.5mm




- E046-315** Hough Micro Aural Forcep, serrated, curved to left, 4.0mm x 0.5mm



- E046-320** Hough Micro Aural Forcep, serrated, curved to right, 4.0mm x 0.5mm



Aural Forceps




E046-305 Hough Micro Aural Forcep, oval cup jaw, straight, 4.0mm x 0.5mm



E042-105 House-Dieter Aural Nipper Forcep, crocodile action, straight up cutting, 7.5cm


E042-110 House-Dieter Aural Nipper Forcep, crocodile action, right side cutting, 7.5cm

E042-115 House-Dieter Aural Nipper Forcep, crocodile action, left side cutting, 7.5cm



E036-505 McGee Aural Forcep, crocodile action, wire closing, straight jaw, 4.0mm x 0.8mm


E036-506 McGee Aural Forcep, crocodile action, wire closing, straight jaw, 8.0mm x 0.8mm




E036-615 Ormerod Aural Forcep, crocodile action, serrated jaw, angled upwards

E036-610 Ormerod Aural Forcep, crocodile action, circular cup, straight jaw, 1.0mm dia. cup

E046-105 Ormerod Micro Aural Forcep, crocodile action, circular cup, straight jaw, 0.5mm dia. cup




E036-605 Ormerod Aural Forcep, crocodile action, elongated cup, straight jaw, 5.0mm x 1.0mm




E036-710 Shea Aural Forcep, crocodile action, oval cup, curved to left, 4.0mm x 0.8mm

E046-410 Shea Micro Aural Forcep, oval cup, curved to left, 3.0mm x 0.5mm




E036-715 Shea Aural Forcep, crocodile action, oval cup, curved to right, 4.0mm x 0.8mm

E046-415 Shea Micro Aural Forcep, oval cup, curved to right, 3.0mm x 0.5mm




Aural Forceps



E036-721 Shea Aural Forcep, crocodile action, oval cup, upward curved jaw, 4.0mm x 0.8mm

E046-420 Shea Micro Aural Forcep, oval cup, curved upward, 3.0mm x 0.5mm




E036-705 Shea Aural Forcep, crocodile action, oval cup, straight jaw, 4.0mm x 0.8mm

E046-405 Shea Micro Aural Forcep, oval cup, straight, 3.0mm x 0.5mm


Aural Instruments




E108-130 Beales Platinectomy & Vein Graft Instrument, knife, for stapedius tendon, no.1, 17cm



E108-145 Beales Platinectomy & Vein Graft Instrument, vein graft applicator, left ear, no.10, 17cm



E108-150 Beales Platinectomy & Vein Graft Instrument, vein graft applicator, right ear, no.9, 17cm




E108-105 Beales Platinectomy & Vein Graft Instrument, curette, for removal of superior bony rim, no.11, 17cm



E108-120 Beales Platinectomy & Vein Graft Instrument, hook, angled left, simpson hall, no.7, 17cm


E108-125 Beales Platinectomy & Vein Graft Instrument, hook, angled right, simpson hall, no.8, 17cm



E108-110 Beales Platinectomy & Vein Graft Instrument, hook, downward angle, for removal of foot plate, no.5, 17cm




E108-115 Beales Platinectomy & Vein Graft Instrument, hook, upward angle, for removal of foot plate, no.6, 17cm



E108-135 Beales Platinectomy & Vein Graft Instrument, osteotome, for elevation of fixed foot plate, no.3, 17cm




E108-155 Beales Platinectomy & Vein Graft Instrument, perforator, for foot plate, no.4, 17cm



E108-140 Beales Platinectomy & Vein Graft Instrument, raspatory, for removal of mucous membrane, no.2, 17cm



E112-110 Beales Tympanoplasty Raspatory, slight curve, 16.5cm, no.1




E112-115 Beales Tympanoplasty Raspatory, small curve, 16.5cm, no.2



E112-120 Beales Tympanoplasty Rugine, double ended, sharp edged, 16.5cm, no.3



E112-125 Beales Tympanoplasty Rugine, double ended, round sharp end, 16.5cm, no.4



Ear, Nose and Throat


Aural Instruments




E112-130 Beales Tympanoplasty Spatula, graft smoothing, 16.5cm, no.5




E112-105 Beales Tympanoplasty Probe, fine blunt end, 16.5cm, no.6



E108-205 Cawthorne (Wullstein) Middle Ear Hook, no.1, slight curve, 16.5cm



E108-210 Cawthorne (Wullstein) Middle Ear Hook, no.2, medium curve, 16.5cm



E108-215 Cawthorne (Wullstein) Middle Ear Hook, no.3, full curve, 16.5cm



E108-220 Cawthorne Middle Ear Hook, no.4, double curve, 16.5cm



E108-225 Cawthorne Middle Ear Hook, no.5, angled 90o, 0.5mm, 16.5cm



E108-230 Cawthorne Middle Ear Hook, no.6, angled 90o, 1.0mm, 16.5cm




E108-235 Cawthorne Middle Ear Hook, no.7, retro curve, 16.5cm



E108-240 Cawthorne Middle Ear Hook, no.8, angled 45o, 0.5mm, 16.5cm

E108-245 Cawthorne (Wullstein) Middle Ear Hook,
no.9, straight, 16.5cm



E120-105 Cawthorne First Incision Knife And Elevator,
16.5cm



E002-010 Cawthorne (Plester) Strut Applicator, 16.5cm



E108-305 Hough Middle Ear Excavator, no.1, whirlybird,
large, right, 15cm



E108-310 Hough Middle Ear Excavator, no.2, whirlybird,
large, left, 15cm




E108-315 Hough Middle Ear Excavator, no.3, whirlybird,
small, right, 15cm




E108-320 Hough Middle Ear Excavator, no.4, whirlybird,
small, left, 15cm




E108-325 Hough Middle Ear Pick, no.5, Stapedectomy, 90o
angle, 0.33mm, 16.5cm



E108-330 Hough Middle Ear Pick, no.6, Stapedectomy, 45o
angle, 0.5mm, 16.5cm




E108-335 Hough Middle Ear Pick, no.7, Stapedectomy, 15o
angle, 0.5mm, 16.5cm




Ear, Nose and Throat

Aural Instruments




E108-340 Hough Middle Ear Excavator, Hook, no.8,
90o angle, 16.5cm



E108-345 Hough Middle Ear Excavator, Hook, no.9,
45o angle, 16.5cm




E108-350 Hough Middle Ear Pick, no.10, curved, 16.5cm



E108-405 Ironside Middle Ear Knife, no.1, 16.5cm



E108-410 Ironside Middle Ear Knife, no.2, 16.5cm




E108-415 Ironside Middle Ear Elevator, no.3, 16.5cm




E108-420 Ironside Middle Ear Needle, sharp point, no.4,
16.5cm


E108-425 Ironside Middle Ear Hook, angled 45o up, no.5, 16.5cm




E108-430 Ironside Middle Ear Hook, angled 45o down, no.6, 16.5cm




E108-435 Ironside Middle Ear Hook, angled 90o, right, no.7, 16.5cm



E108-440 Ironside Middle Ear Hook, angled 90o, left, no.8, 16.5cm




E108-445 Ironside Middle Ear Measuring Probe, no.9, 16.5cm




E108-705 Rosen Middle Ear Knife, meatal incision, no.1, 16.5cm



E108-710 Rosen Middle Ear Elevator, broad, bayonet shape, no.2, 19cm



Aural Instruments




E108-715 Rosen Middle Ear Elevator, narrow, bayonet shape, no.3, 19cm



E108-720 Rosen Middle Ear Elevator, drum lifting, no.4, 16.5cm




E108-725 Rosen Middle Ear Curette, tympanic ring small, no.5, 16.5cm



E108-730 Rosen Middle Ear Curette, tympanic ring, large, no.6, 16.5cm



E108-735 Rosen Middle Ear Probe, fine, no.7, 16.5cm




E108-740 Rosen Middle Ear Mobiliser, slight curve, no.8, 16.5cm



E108-745 Rosen Middle Ear Mobiliser, full curve, no.9, 16.5cm




E108-780 Rosen Middle Ear, Fenestrator, oval window, no.12, 16.5cm




E108-760 Rosen Middle Ear Scoop, oval, size 000, no.11, 16.5cm

E108-765 Rosen Middle Ear Scoop, oval, size 00, no.11,
16.5cm




E108-770 Rosen Middle Ear Scoop, oval, size 0, no.11,
16.5cm




E108-775 Rosen Middle Ear Scoop, oval, size 1, no.11,
16.5cm




E108-815 Schuknecht Middle Ear Hook, upward, 0.4mm,
16cm




E108-810 Schuknecht Middle Ear Hook, downward, 0.4mm,
16cm




E108-825 Schuknecht Middle Ear Hook, upward, 0.7mm,
16cm




E108-820 Schuknecht Middle Ear Hook, downward, 0.7mm,
16cm




E108-805 Schuknecht Middle Ear Elevator, oval window,
16cm




Aural Instruments




E108-830 Schuknecht Middle Ear Hook, angled 90o, upward, 16cm



E108-835 Schuknecht Middle Ear Pick, straight, 16cm




E108-840 Schuknecht Middle Ear Spatula, 16cm



E108-905 Shea Middle Ear Knife, incising, angled, 16cm.(A)



E108-910 Shea Middle Ear Curette, straight, 16cm.(B)



E108-915 Shea Middle Ear Elevator, straight, left, 16cm.(C)



E108-920 Shea Middle Ear Elevator, straight, right, 16cm.(D)



E108-925 Shea Middle Ear Pick, sharp, straight, 16cm.(E)



E108-930 Shea Middle Ear Hook-Fenestra, straight, angled 25o, 16cm.(F)

Aural Instruments

E108-935 Shea Middle Ear Hook-Fenestra, straight, angled 45o, 16cm.(G)




E108-940 Shea Middle Ear Hook-Fenestra, straight, angled 90o, 16cm.(H)




E108-945 Shea Middle Ear Hook-Fenestra, straight, angled 90o, short point, 16cm.(I)



E108-950 Shea Middle Ear Pick, straight, angled 90o, semi blunt, 16cm.(J)



E108-955 Shea Middle Ear Knife, anterior curotomy, straight, 16cm.(K)



E088-205 Shea Piston Depth Gauge, combined depth of 3.5mm, 4.0mm & 4.5mm, 15.5cm



E088-135 House Depth Gauge, 3.5mm depth, 16.5cm

E088-140 House Depth Gauge, 4mm depth, 16.5cm

E088-145 House Depth Gauge, 4.5mm depth, 16.5cm

E088-150 House Depth Gauge, 5mm depth, 16.5cm




E136-105 Hughes Stapes Mobiliser, 16.5cm




Aural Snares

E178-210 Glegg Aural Snare, self forming, 16.5cm



E178-205 Glegg (Negus) Aural Snare, open ended, 16.5cm




Aural Speculum



- E034-612** Hartman Speculum, round end, stainless steel, no.1, 2mm bore
- E034-613** Hartman Speculum, round end, stainless steel, no.2, 3mm bore
- E034-614** Hartman Speculum, round end, stainless steel, no.3, 4mm bore
- E034-615** Hartman Speculum, round end, stainless steel, no.4, 5.5mm bore
- E034-616** Hartman Speculum, round end, stainless steel, no.5, 6.5mm bore
- E034-617** Hartman Speculum, round end, stainless steel, no.6, 7.5mm bore
- E034-618** Hartman Speculum, round end, stainless steel, no.7, 9mm bore





- E034-023** Gruber Speculum, oval end, black finish, no.1, 2.5mm x 3.5mm bore
- E034-024** Gruber Speculum, oval end, black finish, no.2, 3.5mm x 4.5mm bore
- E034-025** Gruber Speculum, oval end, black finish, no.3, 4.0mm x 6.0mm bore
- E034-026** Gruber Speculum, oval end, black finish, no.4, 6.0mm x 8.0mm bore
- E034-021** Gruber Speculum, oval end, black finish, no's 2-4, set /3
- E034-020** Gruber Speculum, oval end, black finish, no's 1-4, set /4
- E034-013** Gruber Speculum, round end, black finish, no.1, 3.0mm bore
- E034-014** Gruber Speculum, round end, black finish, no.2, 4.0mm bore
- E034-015** Gruber Speculum, round end, black finish, no.3, 5.0mm bore
- E034-017** Gruber Speculum, round end, black finish, no.4, 7.0mm bore
- E034-011** Gruber Speculum, round end, black finish, no's 2-4, set /3 sizes, 4, 5 & 7mm bore
- E034-010** Gruber Speculum, round end, black finish, no's 1-4, set /4 sizes, 3, 4, 5 & 7mm bore





- E034-034** Heath Speculum, oval end, black finish, no.1, 4mm x 6.0mm bore
- E034-035** Heath Speculum, oval end, black finish, no.2, 5mm x 7.0mm bore
- E034-036** Heath Speculum, oval end, black finish, no.3, 6mm x 8.0mm bore
- E034-037** Heath Speculum, oval end, black finish, no.4, 7mm x 9.0mm bore
- E034-030** Heath Speculum, oval end, black finish, no's 1-4, set /4 sizes, 4, 5, 6 & 7mm

- E034-155** Rosen Speculum, half slotted, oval end, black finish, no.1, 2.5mm x 3.5mm bore
- E034-175** Rosen Speculum, half slotted, oval end, black finish, no.2, 3.5mm x 4.5mm bore
- E034-185** Rosen Speculum, half slotted, oval end, black finish, no.3, 4mm x 6mm bore
- E034-195** Rosen Speculum, half slotted, oval end, black finish, no.4, 5mm x 7mm bore
- E034-196** Rosen Speculum, half slotted, oval end, black finish, no.5, 6mm x 8mm bore
- E034-100** Rosen Speculum, half slotted, oval end, black finish, no's 2-5, set /4
- E034-105** Rosen Speculum, half slotted, oval end, black finish, no's 1-5, set /5
- E034-083** Rosen Speculum, half slotted, round end, black finish, 3mm bore, no.1
- E034-084** Rosen Speculum, half slotted, round end, black finish, 4mm bore, no.2
- E034-085** Rosen Speculum, half slotted, round end, black finish, 5mm bore, no.3
- E034-086** Rosen Speculum, half slotted, round end, black finish, 6mm bore, no.4
- E034-087** Rosen Speculum, half slotted, round end, black finish, 7mm bore, no.5
- E034-080** Rosen Speculum, half slotted, round end, black finish, no's 2-5, Set/4 sizes, 4, 5, 6 & 7mm bore
- E034-081** Rosen Speculum, half slotted, round end, black finish, no's 1-5, set /5 sizes, 3, 4, 5, 6 & 7mm bore



- E034-092** Rosen Speculum, unslotted, oval end, black finish, no.1, 2.5mm x 3.5mm bore
- E034-094** Rosen Speculum, unslotted, oval end, black finish, no.2, 3.5mm x 4.5mm bore
- E034-095** Rosen Speculum, unslotted, oval end, black finish, no.3, 4mm x 6mm bore
- E034-096** Rosen Speculum, unslotted, oval end, black finish, no.4, 5mm x 7mm bore
- E034-097** Rosen Speculum, unslotted, oval end, black finish, no.5, 6mm x 8mm bore
- E034-090** Rosen Speculum, unslotted, oval end, black finish, no's 2-5 set /4
- E034-089** Rosen Speculum, unslotted, oval end, black finish, no's 1-5 set /5
- E034-073** Rosen Speculum, unslotted, round end, black finish, 3mm bore, no.1
- E034-074** Rosen Speculum, unslotted, round end, black finish, 4mm bore, no.2
- E034-075** Rosen Speculum, unslotted, round end, black finish, 5mm bore, no.3
- E034-076** Rosen Speculum, unslotted, round end, black finish, 6mm bore, no.4
- E034-077** Rosen Speculum, unslotted, round end, black finish, 7mm bore, no.5
- E034-071** Rosen Speculum, unslotted, round end, black finish, no's 2-5, set /4, 4, 5, 6 & 7mm bore
- E034-070** Rosen Speculum, unslotted, round end, black finish, no's 1-5, set /5 sizes, 3, 4, 5, 6 & 7mm bore




Aural Speculum



- E034-245** Shea Speculum, oblique end, black finish, 3.3mm x 4.1mm bore, no.1
- E034-255** Shea Speculum, oblique end, black finish, 4mm x 5.5mm bore, no.2
- E034-275** Shea Speculum, oblique end, black finish, 5mm x 7.5mm bore, no.3
- E034-285** Shea Speculum, oblique end, black finish, 6mm x 8.5mm bore, no.4
- E034-295** Shea Speculum, oblique end, black finish, 7mm x 9.5mm bore, no.5
- E034-205** Shea Speculum, oblique end, black finish, no.s 1-5, set /5





- E034-057** Peters (Siegle) Speculum set, main body only
- E034-059** Peters (Siegle) Speculum set, magnified lens only
- E034-058** Peters (Siegle) Speculum set, rubber tube and bulb only
- E034-049** Peters (Siegle) Pneumatic Speculum Set Complete With Magnified Glass, 5 speculum, rubber tube and bulb, black finish
- E034-056** Peters (Siegle) Speculum Set, plain lens only
- E034-050** Peters (Siegle) Pneumatic Speculum Set Complete With Plain Glass, 5 speculum, rubber tube and bulb, black finish
- E034-051** Peters (Siegle) Speculum Only, black finish, no.0, oval
- E034-052** Peters (Siegle) Speculum Only, black finish, no.1, oval
- E034-053** Peters (Siegle) Speculum Only, black finish, no.2, oval
- E034-054** Peters (Siegle) Speculum Only, black finish, no.3, oval
- E034-055** Peters (Siegle) Speculum Only, black finish, no.4, oval



- E034-063** Toynbee Speculum, round end, chrome finish, no.0, 2mm bore
- E034-065** Toynbee Speculum, round end, chrome finish, no.1, 3mm bore
- E034-066** Toynbee Speculum, round end, chrome finish, no.2, 4mm bore
- E034-067** Toynbee Speculum, round end, chrome finish, no.3, 5mm bore
- E034-068** Toynbee Speculum, round end, chrome finish, no.4, 6mm bore


Aural Speculum

- E034-345** Tumarkin Transmeatal Speculum, black finish, round, slotted, 3.5mm bore, no.0
- E034-355** Tumarkin Transmeatal Speculum, black finish, round, slotted, 4.5mm bore, no.0.5
- E034-365** Tumarkin Transmeatal Speculum, black finish, round, slotted, 5.5mm bore dia, no.1
- E034-370** Tumarkin Transmeatal Speculum, black finish, round, slotted, 6.0mm bore, no.2
- E034-375** Tumarkin Transmeatal Speculum, black finish, round, slotted, 6.5mm bore, no.3
- E034-380** Tumarkin Transmeatal Speculum, black finish, round, slotted, 7.0mm bore, no.4
- E034-385** Tumarkin Transmeatal Speculum, black finish, round, slotted, 7.5mm bore, no.5
- E034-390** Tumarkin Transmeatal Speculum, black finish, round, slotted, 8.0mm bore, no.6
- E034-391** Tumarkin Transmeatal Speculum, black finish, round, slotted, 10.0mm bore, no.7
- E034-300** Tumarkin Transmeatal Speculum, black finish, round, slotted, set /8




Biopsy Forceps


- E050-305** Chevalier Jackson Laryngeal Biopsy Forcep, round basket jaw, 30cm
- E050-307** Chevalier Jackson Laryngeal Biopsy Forcep, round basket jaw, 40cm
- E050-310** Chevalier Jackson Laryngeal Biopsy Forcep, round basket jaw, 50cm
- E050-312** Chevalier Jackson Laryngeal Biopsy Forcep, round basket jaw, 60cm



- E050-105** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 30cm
- E050-110** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 40cm
- E050-115** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 45cm
- E050-120** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 50cm
- E050-125** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 55cm
- E050-130** Chevalier Jackson Laryngeal Biopsy Forcep, square basket jaw, 60cm



- E050-520** Paterson Laryngeal Biopsy Forcep, 20cm
- E050-525** Paterson Laryngeal Biopsy Forcep, 25cm
- E050-530** Paterson Laryngeal Biopsy Forcep, 30cm
- E050-535** Paterson Laryngeal Biopsy Forcep, 35cm
- E050-540** Paterson Laryngeal Biopsy Forcep, 40cm
- E050-545** Paterson Laryngeal Biopsy Forcep, 45cm
- E050-550** Paterson Laryngeal Biopsy Forcep, 50cm
- E050-555** Paterson Laryngeal Biopsy Forcep, 55cm
- E050-560** Paterson Laryngeal Biopsy Forcep, 60cm



Ear, Nose and Throat

Biopsy Forceps



E050-630 Paterson Laryngeal Biopsy Forcep, extra fine, 30cm

E050-635 Paterson Laryngeal Biopsy Forcep, extra fine, 35cm


E050-640 Paterson Laryngeal Biopsy Forcep, extra fine, 40cm

E050-645 Paterson Laryngeal Biopsy Forcep, extra fine, 45cm




E050-802 Bronchoscopic and Oesophagoscopic Forceps, round cup biopsy, 35cm

Bipods




E166-105 Draffin Suspension Rods, child, three ring, 42.5cm, pair/2

E166-110 Draffin Suspension Rods, adult, four Ring, 48cm, pair/2




E166-115 Magauran Plate, for use with Draffin suspension rods 48cm

E016-050 Bone Crusher, 7cm



Bone Nibbling Forceps


E068-080 Cawthorne Bone Nibbling Forcep, 3mm jaw, 16cm



E068-105 Freidman Nibbling Forcep, slight curve, 14cm




E068-205 Jansen Middleton Nibbling Forcep, compound action, angled, 19cm



E068-300 Lempert Nibbling Forcep, straight, 2.5mm bite, 16cm

E068-305 Lempert Nibbling Forcep, curved, 2.5mm bite, 16cm




Bronchoscopes



- E006-105** Bronchoscope Fibre Optic, polished brass, suckling, short, 3.7mm I.D, 18cm, complete with light carrier
- E006-110** Bronchoscope Fibre Optic, polished brass, suckling, small, 3.7mm I.D, 27.5cm, complete with light carrier
- E006-115** Bronchoscope Fibre Optic, polished brass, suckling, 4.1mm I.D, 27.5cm, complete with light carrier
- E006-120** Bronchoscope Fibre Optic, polished brass, infant, 5.4mm I.D, 27.5cm, complete with light carrier
- E006-125** Bronchoscope Fibre Optic, polished brass, child, 7mm I.D, 30cm, complete with light carrier
- E006-130** Bronchoscope Fibre Optic, polished brass, adolescent small, 7mm I.D, 40cm, complete with light carrier
- E006-135** Bronchoscope Fibre Optic, polished brass, adolescent, 8mm I.D, 40cm, complete with light carrier
- E006-140** Bronchoscope Fibre Optic, polished brass, adult small, 9mm I.D, 40cm, complete with light carrier
- E006-145** Bronchoscope Fibre Optic, polished brass, adult, 10mm I.D, 40cm, complete with light carrier
- E006-150** Bronchoscope Fibre Optic, polished brass, lower bronchus, 6.4mm I.D, 45cm, complete with light carrier
- E006-155** Bronchoscope Fibre Optic, aspirating, polished brass, child, 6.5mm I.D, 30cm, complete with light carrier
- E006-160** Bronchoscope Fibre Optic, aspirating, polished brass, adult, 7.5mm I.D, 40cm, complete with light carrier
- E006-180** Bronchoscope Fibre Optic, stainless steel, suckling, 4.1mm I.D, 27.5cm, complete with light carrier
- E006-182** Bronchoscope Fibre Optic, stainless steel, infant, 5.4mm I.D, 27.5cm, complete with light carrier
- E006-184** Bronchoscope Fibre Optic, stainless steel, child, 7mm I.D, 30cm, complete with light carrier
- E006-186** Bronchoscope Fibre Optic, stainless steel, adolescent small, 7mm I.D, 40cm, complete with light carrier
- E006-188** Bronchoscope Fibre Optic, stainless steel, adolescent, 8mm I.D, 40cm, complete with light carrier
- E006-190** Bronchoscope Fibre Optic, stainless steel, adult small, 9mm I.D, 40cm, complete with light carrier
- E006-192** Bronchoscope Fibre Optic, stainless steel, adult, 10mm I.D, 40cm, complete with light carrier
- E006-194** Bronchoscope Fibre Optic, stainless steel, lower bronchus, 6.4mm I.D, 45cm, complete with light carrier
- E006-196** Bronchoscope Fibre Optic, aspirating, stainless steel, child, 6.5mm I.D, 30cm, complete with light carrier
- E006-198** Bronchoscope Fibre Optic, aspirating, stainless steel, adult, 7.5mm I.D, 40cm, complete with light carrier
- E006-200** Fibre Light Carrier Only, for 40cm Bronchoscope
- E006-202** Fibre Light Carrier Only, for 30cm Bronchoscope
- E006-204** Fibre Light Carrier Only, for 27.5cm Bronchoscope infant
- E006-206** Fibre Light Carrier Only, for 27.5cm Bronchoscope Suckling
- E006-208** Fibre Light Carrier Only, for 18cm Bronchoscope
- E006-210** Fibre Light Carrier Only, for 45cm Bronchoscope

Cannula

- E008-120** Eicken Sinus Suction / Washout Cannula, 3mm dia, shallow curve, fig.A, 14cm
- E008-122** Eicken Sinus Suction / Washout Cannula, 3mm dia, medium curve, fig.B, 14cm
- E008-124** Eicken Sinus Suction / Washout Cannula, 3mm dia, full curve, fig.C, 14cm
- E008-126** Eicken Sinus Suction / Washout Cannula, 4mm dia, shallow curve, fig.A, 14cm
- E008-128** Eicken Sinus Suction / Washout Cannula, 4mm dia, medium curve, fig.B, 14cm
- E008-129** Eicken Sinus Suction / Washout Cannula, 4mm dia, full curve, fig.C, 14cm
- E008-130** Eicken Sinus Suction / Washout Cannula, 2.5mm dia, shallow curve, fig.A, 14cm
- E008-132** Eicken Sinus Suction / Washout Cannula, 2.5mm dia, medium curve, fig.B, 14cm
- E008-134** Eicken Sinus Suction / Washout Cannula, 2.5mm dia, full curve, fig.C, 14cm



Cerumen Hooks

- E096-105** Cawthorne Cerumen Hook, 14.5cm
- E096-205** Formby Scoop and Hook, double ended, 18cm
- E096-405** Hartman Hook, sharp, 15cm
- E096-624** Lucae Button Hook, probe tip, large, 4mm, 14cm
- E096-622** Lucae Button Hook, probe tip, medium, 3mm, 14cm
- E096-620** Lucae Button Hook, probe tip, small, 2mm, 14cm
- E096-505** St. Barts Hook, blunt, 17cm




Cocanising Wires




- E150-405** Kings College Hospital Cocaine Wire Probe, 9cm
- E150-410** Kings College Hospital Cocaine Wire Probe, 13.5cm

Conchotomes



- E014-105** Hartman Conchotome, crocodile action, 5mm dia cut, 12cm effective length
- E014-107** Hartman Conchotome, crocodile action, 7mm dia cut, 12cm effective length
- E014-109** Hartman Conchotome, crocodile action, 9mm dia cut, 12cm effective length
- E014-111** Hartman Conchotome, crocodile action, 11mm dia cut, 12cm effective length


Curettes



- E018-110** Barnhill Adenoid Curette, length 21.5cm, no.0
- E018-112** Barnhill Adenoid Curette, length 21.5cm, no.1
- E018-114** Barnhill Adenoid Curette, length 21.5cm, no.2
- E014-116** Barnhill Adenoid Curette, length 21.5cm, no.3
- E014-118** Barnhill Adenoid Curette, length 21.5cm, no.4




- E018-208** Beckmann Adenoid Curette, no.1, 10mm, 21.5cm
- E018-210** Beckmann Adenoid Curette, no.2, 11mm, 21.5cm
- E018-212** Beckmann Adenoid Curette, no.3, 12mm, 21.5cm
- E018-214** Beckmann Adenoid Curette, no.4, 13mm, 21.5cm
- E018-216** Beckmann Adenoid Curette, no.5, 15mm, 21.5cm
- E018-218** Beckmann Adenoid Curette, no.6, 17mm, 21.5cm



- E020-220** Buck Curette, angled, blunt, size 00, 2mm
- E020-225** Buck Curette, angled, blunt, size 0, 2.5mm
- E020-230** Buck Curette, angled, blunt, size 1, 3mm
- E020-235** Buck Curette, angled, blunt, size 2, 3.5mm
- E020-240** Buck Curette, angled, blunt, size 3, 4mm
- E020-320** Buck Curette, angled, sharp, size 00, 2mm
- E020-325** Buck Curette, angled, sharp, size 0, 2.5mm
- E020-330** Buck Curette, angled, sharp, size 1, 3mm
- E020-335** Buck Curette, angled, sharp, size 2, 3.5mm
- E020-340** Buck Curette, angled, sharp, size 3, 4mm

- E020-260** Buck Curette, straight, blunt, size 00, 2mm
E020-265 Buck Curette, straight, blunt, size 0, 2.5mm
E020-270 Buck Curette, straight, blunt, size 1, 3mm
E020-275 Buck Curette, straight, blunt, size 2, 3.5mm
E020-280 Buck Curette, straight, blunt, size 3, 4mm
E020-360 Buck Curette, straight, sharp, size 00, 2mm
E020-365 Buck Curette, straight, sharp, size 0, 2.5mm
E020-370 Buck Curette, straight, sharp, size 1, 3mm
E020-375 Buck Curette, straight, sharp, size 2, 3.5mm
E020-380 Buck Curette, straight, sharp, size 3, 4mm



- E022-350** Williger Curette, double ended, small, 13.5cm
E022-352 Williger Curette, double ended, medium, 13.5cm
E022-354 Williger Curette, double ended, large, 13.5cm



- E023-822** Antrum Curette, narrow oval cup, 19cm
E023-824 Antrum Curette, narrow oval cup, curved 55o, 19cm
E023-826 Antrum Curette, narrow oval cup, curved 90o, 19cm




- E022-705** Volkmann Curette, double ended, small, 3mm / 5mm wide, 22cm
E022-710 Volkmann Curette, double ended, medium, 4mm / 6mm wide, 22cm
E022-715 Volkmann Curette, double ended, large, 5mm / 7mm wide, 22cm



- E022-301** House Curette, 1.5mm / 1.8mm cup, slightly curved, double ended, 15cm
E022-304 House Curette, 1.5mm / 1.8mm cup, fully curved, double ended, 15cm




- E022-305** House Curette, 2mm / 2.5mm cup, slightly curved, double ended, 18cm
E022-310 House Curette, 2mm / 2.5mm cup, strong curved, double ended, 18cm




Ear, Nose and Throat

Curettes




- E022-406** Mackie Ethmoid Curette, 6mm diameter, 19.5cm
- E022-408** Mackie Ethmoid Curette, 8mm diameter, 19.5cm
- E022-409** Mackie Ethmoid Curette, 9mm diameter, 19.5cm



- E018-408** St. Clair Thomson Adenoid Curette, 8mm, 23.5cm
- E018-410** St. Clair Thomson Adenoid Curette, 10mm, 23.5cm
- E018-412** St. Clair Thomson Adenoid Curette, 12mm, 23.5cm
- E018-414** St. Clair Thomson Adenoid Curette, 14mm, 23.5cm
- E018-416** St. Clair Thomson Adenoid Curette, 16mm, 23.5cm
- E018-418** St. Clair Thomson Adenoid Curette, 18mm, 23.5cm


Dissectors




- E030-105** Beavis Tonsil Dissector & Pillar Retractor, double ended, 20cm



- E028-306** Gwynne Evans Tonsil Dissector, double ended, 6mm wide
- E028-308** Gwynne Evans Tonsil Dissector, double ended, 8mm wide



- E030-110** Hurd Tonsil Dissector & Pillar Retractor, double ended, 20cm



- E030-015** Mollison Tonsil Dissector & Pillar Retractor, double ended, 20cm



- E028-405** Negus Aspiring Dissector, 21cm



- E032-225** Paulsen Dissector, curved, 2.5mm
- E032-227** Paulsen Dissector, curved, 3.5mm
- E032-229** Paulsen Dissector, curved, 4.5mm



- E028-510** St. Clair Thomson Dissector, blunt, 16cm

Dissectors

E028-505 St. Clair Thomson Dissector, sharp, 16cm



E028-605 Syme Periosteal Elevator / Dissector, double ended, 16.5cm




E196-950 Weir Suction / Dissecting Tube, 19cm




Elevators

E032-210 Cottle Elevator, single ended, 16cm




E032-205 Cottle Elevator, double ended, graduated, 21cm




E032-215 Cottle Septum Elevator / Dissector, double ended, curved, 23.5cm




E032-315 Freer Elevator, double ended, sharp / blunt, 18cm



E032-320 Freer Elevator, double ended, with knife curved to side, 19cm



E158-105 Hill Elevator / Raspatory, double ended, 21.5cm



Gouges

E092-105 Guys Nasal Gouge, bayonet shape, 12.5cm



Ear, Nose and Throat

Gouges



- E090-122** Jenkins Gouge, 2mm wide 14.5cm
- E090-123** Jenkins Gouge, 3mm wide 14.5cm
- E090-124** Jenkins Gouge, 4mm wide 14.5cm
- E090-125** Jenkins Gouge, 5mm wide 14.5cm
- E090-126** Jenkins Gouge, 6mm wide 14.5cm
- E090-127** Jenkins Gouge, 7mm wide 14.5cm
- E090-128** Jenkins Gouge, 8mm wide 14.5cm
- E090-129** Jenkins Gouge, 9mm wide 14.5cm
- E090-130** Jenkins Gouge, 10mm wide 14.5cm
- E090-132** Jenkins Gouge, 12mm wide 14.5cm
- E090-134** Jenkins Gouge, 14mm wide 14.5cm




- E092-205** Killian Nasal Gouge, bayonet, 16cm




- E092-305** Tilley Nasal Gouge, Bayonet Shape, 16cm

Laryngeal Forceps - Micro 25cm length



Standard Micro Laryngeal Forcep

- E061-105** 2mm shaft diameter, oval cup straight, 2.5mm x 4.0mm cup, 25cm length
- E061-205** 2mm shaft diameter, round cup straight, 2.5mm cup, 25cm length
- E061-210** 2mm shaft diameter, round cup curved left, 2.5mm cup, 25cm length
- E061-215** 2mm shaft diameter, round cup curved right, 2.5mm cup, 25cm length
- E061-220** 2mm shaft diameter, round cup curved upwards, 2.5mm cup, 25cm length
- E061-225** 2mm shaft diameter, round cup curved downwards, 2.5mm cup, 25cm length



Standard Micro Laryngeal Scissor

- E061-305** 2mm shaft diameter, straight, 25cm length
- E061-310** 2mm shaft diameter, curved left, 25cm length
- E061-315** 2mm shaft diameter, curved right, 25cm length
- E061-329** 2mm shaft diameter, straight horizontal, 25cm length



Standard Micro Laryngeal Grasping Forcep


- E061-405** 2mm shaft diameter, serrated jaw straight, 25cm length
- E061-410** 2mm shaft diameter, serrated jaw curved left, 25cm length
- E061-415** 2mm shaft diameter, serrated jaw curved right, 25cm length
- E061-420** 2mm shaft diameter, serrated jaw straight horizontal, 25cm length

Laryngeal Forceps


- E060-230** Chevalier Jackson Laryngeal Grasping Forcep, 30cm
- E060-235** Chevalier Jackson Laryngeal Grasping Forcep, 35cm
- E060-240** Chevalier Jackson Laryngeal Grasping Forcep, 40cm
- E060-245** Chevalier Jackson Laryngeal Grasping Forcep, 45cm
- E060-255** Chevalier Jackson Laryngeal Grasping Forcep, 55cm



- E060-130** Chevalier Jackson Laryngeal Grasping Forcep, flat round serrated jaw, 30cm




- E060-430** Chevalier Jackson Laryngeal Grasping Forcep, rotation, 30cm
- E060-435** Chevalier Jackson Laryngeal Grasping Forcep, rotation, 35cm
- E060-440** Chevalier Jackson Laryngeal Grasping Forcep, rotation, 40cm
- E060-455** Chevalier Jackson Laryngeal Grasping Forcep, rotation, 55cm




Laryngoscopes

- E128-520** Negus Laryngoscope, polished brass, complete with twin fibre light carrier, infant, effective length 82mm
- E128-522** Negus Laryngoscope, stainless steel, complete with twin fibre light carrier, infant, effective length 82mm
- E128-525** Negus Laryngoscope, polished brass, complete with twin fibre light carrier, child, effective length 99mm
- E128-527** Negus Laryngoscope, stainless steel, complete with twin fibre light carrier, child, effective length 99mm
- E128-530** Negus Laryngoscope, polished brass, complete with twin fibre light carrier, adult, effective length 135mm
- E128-532** Negus Laryngoscope, stainless steel, complete with twin fibre light carrier, adult, effective length 135mm
- E128-620** Negus Fibre Light Carrier only, infant
- E128-625** Negus Fibre Light Carrier only, child
- E128-630** Negus Fibre Light Carrier only, adult



Ear, Nose and Throat

Laryngoscopes



- E128-305** RNTNE Laryngoscope, black finish, small, 20mm x 15mm proximal, 15mm x 11mm distal, effective length 120mm (fibre light carrier not included)
- E128-310** RNTNE Laryngoscope, black finish, medium, 24mm x 15.5mm proximal, 15mm x 11mm distal, effective length 150mm (fibre light carrier not included)
- E128-315** RNTNE Laryngoscope, black finish, large, 28mm x 19mm proximal, 18mm x 13mm distal, effective length 152mm (fibre light carrier not included)
- E128-405** RNTNE Anterior Commissure Laryngoscope, black finish, small, 21mm x 19mm proximal, 15mm x 14mm distal, effective length 123mm, (fibre light carrier not included)
- E128-410** RNTNE Anterior Commissure Laryngoscope, black finish, medium, 22mm x 20mm proximal, 15mm x 14mm distal, effective length 149mm, (fibre light carrier not included)
- E128-415** RNTNE Anterior Commissure Laryngoscope, black finish, large, 28mm x 26mm proximal, 18mm x 15mm distal, effective length 149mm, (fibre light carrier not included)
- E128-500** Fibre Light Carrier Only, for use with RNTNE laryngoscopes

Mirrors



- E134-106** Laryngeal Mirror, fixed handle, 6mm, no.000
- E134-108** Laryngeal Mirror, fixed handle, 8mm, no.00
- E134-110** Laryngeal Mirror, fixed handle, 10mm, no.0
- E134-112** Laryngeal Mirror, fixed handle, 12mm, no.1
- E134-114** Laryngeal Mirror, fixed handle, 14mm, no.2
- E134-116** Laryngeal Mirror, fixed handle, 16mm, no.3
- E134-118** Laryngeal Mirror, fixed handle, 18mm, no.4
- E134-120** Laryngeal Mirror, fixed handle, 20mm, no.5
- E134-122** Laryngeal Mirror, fixed handle, 22mm, no.6
- E134-124** Laryngeal Mirror, fixed handle, 24mm, no.7
- E134-126** Laryngeal Mirror, fixed handle, 26mm, no.8
- E134-128** Laryngeal Mirror, fixed handle, 28mm, no.9
- E134-130** Laryngeal Mirror, fixed handle, 30mm, no.10




- E134-306** Laryngeal Mirror, without handle, 6mm, no.000
- E134-308** Laryngeal Mirror, without handle, 8mm, no.00
- E134-310** Laryngeal Mirror, without handle, 10mm, no.0
- E134-312** Laryngeal Mirror, without handle, 12mm, no.1
- E134-314** Laryngeal Mirror, without handle, 14mm, no.2
- E134-316** Laryngeal Mirror, without handle, 16mm, no.3
- E134-318** Laryngeal Mirror, without handle, 18mm, no.4
- E134-320** Laryngeal Mirror, without handle, 20mm, no.5
- E134-322** Laryngeal Mirror, without handle, 22mm, no.6
- E134-324** Laryngeal Mirror, without handle, 24mm, no.7
- E134-326** Laryngeal Mirror, without handle, 26mm, no.8
- E134-328** Laryngeal Mirror, without handle, 28mm, no.9
- E134-330** Laryngeal Mirror, without handle, 30mm, no.10

- E084-105** Boyle Davis Mouth Gag, child, 13.5cm x 7cm
(15mm wide fitting)
- E084-110** Boyle Davis Mouth Gag, adult, 15cm x 9cm
(15mm wide fitting)
- E084-150** Boyle Davis replaceable tooth guard, silicone,
pair, for use with E084-105 and E084-110



- E084-205** Boyle Davis Tongue Plate, plain, 6.5cm,
(blade interchangeable with all 15mm mouth gags)
- E084-210** Boyle Davis Tongue Plate, plain, 7.5cm,
(blade interchangeable with all 15mm mouth gags)
- E084-215** Boyle Davis Tongue Plate, plain, 9.0cm,
(blade interchangeable with all 15mm mouth gags)
- E084-220** Boyle Davis Tongue Plate, plain, 10.0cm,
(blade interchangeable with all 15mm mouth gags)
- E084-225** Boyle Davis Tongue Plate, plain, 11.5cm,
(blade interchangeable with all 15mm mouth gags)



- E084-405** Doughty Tongue Plate, slotted, 6.5cm,
(blade interchangeable with all 15mm mouth gags)
- E084-410** Doughty Tongue Plate, slotted, 7.5cm,
(blade interchangeable with all 15mm mouth gags)
- E084-415** Doughty Tongue Plate, slotted, 9.0cm,
(blade interchangeable with all 15mm mouth gags)
- E084-420** Doughty Tongue Plate, slotted, 10.0cm,
(blade interchangeable with all 15mm mouth gags)
- E084-425** Doughty Tongue Plate, slotted, 11.5cm,
(blade interchangeable with all 15mm mouth gags)




Mouth Gags



- E084-520** Dingman Tongue Plate, large,
(blade interchangeable with all 15mm mouth gags)
- E084-515** Dingman Tongue Plate, medium,
(blade interchangeable with all 15mm mouth gags)
- E084-510** Dingman Tongue Plate, small,
(blade interchangeable with all 15mm mouth gags)
- E084-525** Dingman Mouth Gag, frame only, (15mm wide fitting)
- E084-505** Dingman Mouth Gag, complete, 1 frame and 3 tongue plates (15mm wide fitting)




- E086-205** Fergusson Mouth Gag, ackland jaw, child, 17cm
- E086-210** Fergusson Mouth Gag, ackland jaw, adult, 20cm




- E086-110** Doyen Mouth Gag, with cam rack, adult, 12.5cm
- E086-105** Doyen Mouth Gag, with cam rack, child, 11.5cm

Mouth Gags

- E084-620** Kilner Dott (East Grinstead) Tongue Plate, large, (blade interchangeable with all 10mm mouth gags)
- E084-615** Kilner Dott (East Grinstead) Tongue Plate, medium, (blade interchangeable with all 10mm mouth gags)
- E084-610** Kilner Dott (East Grinstead) Tongue Plate, small, (blade interchangeable with all 10mm mouth gags)
- E084-625** Kilner Dott Mouth Gag, frame only (10mm wide fitting)
- E084-605** Kilner Dott Mouth Gag, complete set, 1 frame and 3 tongue plates (10mm wide fitting)




- E084-720** Kilner Dott (East Grinstead) Tongue Plate, insulated, large, (blade interchangeable with all 10mm mouth gags)
- E084-715** Kilner Dott (East Grinstead) Tongue Plate, insulated, medium, (blade interchangeable with all 10mm mouth gags)
- E084-710** Kilner Dott (East Grinstead) Tongue Plate, insulated, small, (blade interchangeable with all 10mm mouth gags)
- E084-725** Kilner Dott Mouth Gag, frame only, insulated (10mm wide fitting)
- E084-705** Kilner Dott Mouth Gag, complete set, insulated, 1 frame and 3 tongue plates (10mm wide fitting)



Mouth Jacks

- E116-105** Negus Mouth Jack, with rack action to raise and lower



Myringotomes



E138-105 Agnew Politzer Myringotome, angled, round handle, 19cm




E138-220 Lucae Myringotome, bayonet shape, large lancet for adults, 18cm


E138-221 Lucae Myringotome, bayonet shape, small lancet for children, 18cm



E138-315 Sexton Knife, bayonet shape, 17cm



E138-410 Trautmann Myringotome, downward cutting, 18cm



E138-405 Trautmann Myringotome, upward cutting, 18cm

Nasal Forceps



E040-275 Grunwald (Jansen) Forcep, bayonet, serrated, 16cm

E040-277 Grunwald (Jansen) Forcep, bayonet, serrated, 18cm

E040-280 Grunwald (Jansen) Forcep, bayonet, serrated, 21.5cm


E040-285 Grunwald (Jansen) Forcep, bayonet, 1:2 teeth, 16cm

E040-295 Grunwald (Jansen) Forcep, bayonet, 1:2 teeth, 21.5cm


Nasal Forceps

E054-115 Hartman Nasal / Dressing Forcep, serrated, angled, 16cm

E054-110 Hartman Nasal / Dressing Forcep, serrated, angled, 20cm



E054-105 Tilley Nasal / Dressing Forcep, serrated jaw, 8.2cm working length, 16.5cm overall length



Nasal Instruments

E160-105 Cottle Alar Protector, 15cm



E160-605 Cottle Retractor, upper lateral exposing, 16cm



E162-175 Cottle Tenaculum Hook, retrograde small, 15cm



E162-177 Cottle Tenaculum Hook, retrograde large, 15cm




E158-305 McKenty Elevator / Raspatory, double ended, 22cm




Ear, Nose and Throat


Nasal Instruments




E012-175 Freer Nasal Chisel, straight, 4mm wide, 16cm



E032-310 Freer Nasal Elevator, single ended, blunt, 19cm




E032-305 Freer Nasal Elevator, single ended, sharp, 19cm




E032-320 Freer Nasal Elevator, double ended with knife, curved to side, 19cm


E032-312 Freer Elevator, single ended, curved to side, 19cm




E122-110 Freer Cartilage Dissecting Knife, flat round edge, 19cm



E122-115 Freer L-Flap Dissecting Knife, round edge, curved on flat, 19cm



E122-105 Freer Incision Knife, round edged, 19cm



E032-350 Freer Nasal Separator And Elevator, double ended, 19cm



E168-210 Lempert Rugine, 7mm broad curved blade, 18cm



E168-205 Lempert Rugine, 3mm narrow curved blade, 19cm

E118-120 Ballenger Nasal Cartilage Knife, cranked stem, revolving blade, 3mm wide

E118-125 Ballenger Nasal Cartilage Knife, cranked stem, revolving blade, 4mm wide

E118-130 Ballenger Nasal Cartilage Knife, cranked stem, revolving blade, 5mm wide




E118-105 Ballenger Nasal Cartilage Knife, straight stem, revolving blade, 3mm wide

E118-110 Ballenger Nasal Cartilage Knife, straight stem, revolving blade, 4mm wide

E118-115 Ballenger Nasal Cartilage Knife, straight stem, revolving blade, 5mm wide



E118-135 Cottle Knife, single ended, straight, 13.5cm



E122-410 Joseph Nasal Knife, straight blade, sharp tip, 16cm

E122-412 Joseph Nasal Knife, curved blade, sharp tip, 16cm

E122-414 Joseph Nasal Knife, straight blade, blunt tip, 16cm

E122-416 Joseph Nasal Knife, curved blade, blunt tip, 16cm

E122-417 Joseph Button-End Nasal Knife, angular blade, 17cm



E125-110 Sickle Knife, sharp, 19cm



E126-205 Thorburn (Plester) Sickle Knife, double edged, large blade, 16cm

E126-210 Thorburn (Plester) Sickle Knife, double edged, small blade, 16cm




Ear, Nose and Throat

Nasal Raspers



E156-105 Joseph Nasal Rasp, fine 1mm teeth, 16.5cm



E156-110 Joseph Nasal Rasp, coarse 1.5mm teeth, 16.5cm



E156-205 Maltz Nasal Rasp, straight, 18.5cm

E156-210 Maltz Nasal Rasp, curved, 18.5cm



E156-215 Maltz Nasal Rasp, forward cutting teeth, 17cm


E156-220 Maltz Nasal Rasp, reverse cutting teeth, 17cm



E156-115 Joseph Nasal Raspatory / Elevator, smooth end, straight, 16.5cm

E156-120 Joseph Nasal Raspatory / Elevator, smooth end, curved, 16.5cm

Nasal Snares




E178-220 Glegg (Negus) Nasal Snare, self-forming, 23cm




E178-215 Glegg (Negus) Nasal Snare, open ended, 23cm


- E182-250** Killian Nasal Speculum, 5.0cm insulated blade
- E182-265** Killian Nasal Speculum, 6.5cm insulated blade
- E182-275** Killian Nasal Speculum, 7.5cm insulated blade
- E182-290** Killian Nasal Speculum, 9.0cm insulated blade
- E182-200** Killian Nasal Speculum, insulated blade, set /4



- E182-150** Killian Nasal Speculum, 5.0cm blade
- E182-165** Killian Nasal Speculum, 6.5cm blade
- E182-175** Killian Nasal Speculum, 7.5cm blade
- E182-190** Killian Nasal Speculum, 9.0cm blade
- E182-100** Killian Nasal Speculum, set of 4 blade sizes, 5.0cm, 6.5cm, 7.5cm and 9.0cm



- E182-010** Cottle Nasal Speculum, satin finish, no.1, 20mm blade
- E182-015** Cottle Nasal Speculum, satin finish, no.2, 32mm blade
- E182-020** Cottle Nasal Speculum, satin finish, no.3, 41mm blade
- E182-025** Cottle Nasal Speculum, satin finish, no.4, 60mm blade
- E182-050** Cottle Nasal Speculum, set /4




- E182-305** Shah Nasal Speculum, left
- E182-310** Shah Nasal Speculum, right




- E182-315** Shah Nasal Speculum, insulated, left
- E182-320** Shah Nasal Speculum, insulated, right




Nasal Speculum



- E182-438** St. Clair Thomson Nasal Speculum, 38mm blade
- E182-451** St. Clair Thomson Nasal Speculum, 51mm blade
- E182-464** St. Clair Thomson Nasal Speculum, 64mm blade
- E182-476** St. Clair Thomson Nasal Speculum, 76mm blade
- E182-400** St. Clair Thomson Nasal Speculum, set /4




- E182-501** Thudicum Nasal Speculum, no.1, 12mm deep blade
- E182-502** Thudicum Nasal Speculum, no.2, 14mm deep blade
- E182-503** Thudicum Nasal Speculum, no.3, 16mm deep blade
- E182-504** Thudicum Nasal Speculum, no.4, 18mm deep blade
- E182-505** Thudicum Nasal Speculum, no.5, 20mm deep blade
- E182-506** Thudicum Nasal Speculum, no.6, 22mm deep blade
- E182-507** Thudicum Nasal Speculum, no.7, 24mm deep blade
- E182-508** Thudicum Nasal Speculum, no.8, 27mm deep blade
- E182-500** Thudicum Nasal Speculum, set /7, no.1 to 7




- E182-601** Thudicum Nasal Speculum, insulated, no.1, 12mm deep blade
- E182-602** Thudicum Nasal Speculum, insulated, no.2, 14mm deep blade
- E182-603** Thudicum Nasal Speculum, insulated, no.3, 16mm deep blade
- E182-604** Thudicum Nasal Speculum, insulated, no.4, 18mm deep blade
- E182-605** Thudicum Nasal Speculum, insulated, no.5, 20mm deep blade
- E182-606** Thudicum Nasal Speculum, insulated, no.6, 22mm deep blade
- E182-607** Thudicum Nasal Speculum, insulated, no.7, 24mm deep blade
- E182-608** Thudicum Nasal Speculum, insulated, no.8, 27mm deep blade
- E182-600** Thudicum Nasal Speculum, insulated, set /7, no.1 to 7

Nasal Speculum


- E182-055** Tieck Halle Nasal Speculum, small blade, fixed joint, 13.5cm
- E182-057** Tieck Halle Nasal Speculum, large blade, fixed joint, 13.5cm



- E182-060** Hartman Nasal Speculum, small blade, fixed joint, 15cm
- E182-062** Hartman Nasal Speculum, medium blade, fixed joint, 15cm
- E182-064** Hartman Nasal Speculum, large blade, fixed joint, 15cm




- E182-295** Beckmann (Lempert) Endaural Speculum, medium, 13.5cm




Nodule Forceps

- E059-030** Chevalier Jackson Foreign Body Forcep, fenestrated jaw, lightweight, 30cm
- E059-035** Chevalier Jackson Foreign Body Forcep, fenestrated jaw, lightweight, 35cm
- E059-040** Chevalier Jackson Foreign Body Forcep, fenestrated jaw, lightweight, 40cm
- E059-045** Chevalier Jackson Foreign Body Forcep, fenestrated jaw, lightweight, 45cm
- E059-050** Chevalier Jackson Foreign Body Forcep, fenestrated jaw, lightweight, 50cm




Perforators




- E148-105** Myles Naso Antrum Perforator, 17cm, 16fg
E148-110 Myles Naso Antrum Perforator, 17cm, 19fg

Punch Forceps



- E033-305** Ostrum Antrum Punch Forcep, classic pattern, reverse cut, 18cm




- E033-310** Ostrum Antrum Punch Forcep, tubular shaft, backward cut, straight 2.5mm x 8mm bite, adult
E033-312 Ostrum Antrum Punch Forcep, tubular shaft, backward cut, right cutting 2.5mm x 8mm bite, adult
E033-314 Ostrum Antrum Punch Forcep, tubular shaft, backward cut, left cutting 2.5mm x 8mm bite, adult
E033-316 Ostrum Antrum Punch Forcep, tubular shaft, backward cut, straight 1.5mm x 6mm bite, paediatric
E033-318 Ostrum Antrum Punch Forcep, tubular shaft, backward cut, right cutting 1.5mm x 6mm bite, paediatric
E033-320 Ostrum Antrum Punch Forcep, tubular shaft, backward cut, left cutting 1.5mm x 6mm bite, paediatric
E033-325 Ostrum Antrum Punch Forcep, rotating 360o, backward cut, 2.5mm x 8mm bite, adult
E033-335 Ostrum Antrum Punch Forcep, rotating 360o, backward cut, 1.5mm x 6mm bite, paediatric

Punch Forceps

- E033-102** Citelli Antrum Punch Forcep, 2mm bite, 7.5cm
E033-104 Citelli Antrum Punch Forcep, 4mm bite, 7.5cm




- E074-698** Ferris Smith Kerrison Punch Forcep, ring handle, 90o up cutting, 3mm, 15cm effective length
E074-699 Ferris Smith Kerrison Punch Forcep, ring handle, 90o up cutting, 5mm, 15cm effective length
E074-700 Ferris Smith Kerrison Punch Forcep, ring handle, 90o down cutting, 3mm, 15cm effective length
E074-701 Ferris Smith Kerrison Punch Forcep, ring handle, 90o down cutting, 5mm, 15cm effective length
E074-708 Ferris Smith Kerrison Punch Forcep, ring handle, 40o up cutting, 3mm, 15cm effective length
E074-709 Ferris Smith Kerrison Punch Forcep, ring handle, 40o up cutting, 5mm, 15cm effective length



- E074-901** Weil Blakesley Punch Forcep, straight, 3.5mm bite, 11cm effective length
E074-902 Weil Blakesley Punch Forcep, straight, 4.5mm bite, 11cm effective length
E074-903 Weil Blakesley Punch Forcep, straight, 5mm bite, 11cm effective length
E074-904 Weil Blakesley Punch Forcep, straight, 5.5mm bite, 11cm effective length
E074-905 Weil Blakesley Punch Forcep, straight, 6mm bite, 11cm effective length
E074-906 Weil Blakesley Punch Forcep, straight, 6.5mm bite, 11cm effective length
E074-911 Weil Blakesley Punch Forcep, 45o angle, 3.5mm bite, 11cm effective length
E074-912 Weil Blakesley Punch Forcep, 45o angle, 4.5mm bite, 11cm effective length
E074-913 Weil Blakesley Punch Forcep, 45o angle, 5mm bite, 11cm effective length
E074-914 Weil Blakesley Punch Forcep, 45o angle, 5.5mm bite, 11cm effective length
E074-915 Weil Blakesley Punch Forcep, 45o angle, 6mm bite, 11cm effective length
E074-916 Weil Blakesley Punch Forcep, 45o angle, 6.5mm bite, 11cm effective length
E074-921 Weil Blakesley Punch Forcep, 90o angle, 3.5mm bite, 11cm effective length
E074-922 Weil Blakesley Punch Forcep, 90o angle, 4.5mm bite, 11cm effective length
E074-923 Weil Blakesley Punch Forcep, 90o angle, 5mm bite, 11cm effective length
E074-924 Weil Blakesley Punch Forcep, 90o angle, 5.5mm bite, 11cm effective length
E074-925 Weil Blakesley Punch Forcep, 90o angle, 6mm bite, 11cm effective length
E074-926 Weil Blakesley Punch Forcep, 90o angle, 6.5mm bite, 11cm effective length



Punch Forceps



E033-205 Hajek Antrum Punch Forcep, 19.5cm




E074-230 Hajek Sphenoidal Punch Forcep, 4mm downward bite


E074-240 Hajek Sphenoidal Punch Forcep, 5mm downward bite

E074-225 Hajek Sphenoidal Punch Forcep, 4mm upward bite

E074-235 Hajek Sphenoidal Punch Forcep, 5mm upward bite



E033-405 Seymour Jones Antrum Punch Forcep, 16cm



E074-515 St. Barts Ethmoid Punch Forcep, downward cutting, 3mm bite, 11.5cm

E074-505 St. Barts Ethmoid Punch Forcep, straight, 3mm bite, 11.5cm

E074-510 St. Barts Ethmoid Punch Forcep, upward cutting, 3mm bite, 11.5cm




E075-025 Weil Blakesley Punch Forcep, cut through jaw, straight, 2.5mm bite, 12cm effective length

E075-035 Weil Blakesley Punch Forcep, cut through jaw, straight, 3.5mm bite, 12cm effective length

E075-048 Weil Blakesley Punch Forcep, cut through jaw, straight, 4.5mm bite, 12cm effective length

Punch Forceps

- E074-405** Tilley Henckel Ethmoid Punch Forcep, graduated upper shaft, 4.5mm bite, 11.5cm
- E074-410** Tilley Henckel Ethmoid Punch Forcep, graduated upper shaft, 5.5mm bite, 11.5cm
- E074-415** Tilley Henckel Ethmoid Punch Forcep, graduated upper shaft, 6.5mm bite, 11.5cm
- E074-420** Tilley Henckel Ethmoid Punch Forcep, graduated upper shaft, 8.0mm bite, 11.5cm



Retractors

- E160-205** Aufricht Nasal Retractor, double angled, 17cm



- E162-105** Tracheal Hook Retractor, single, blunt, 16cm
- E162-110** Tracheal Hook Retractor, single, sharp, 16cm
- E162-115** Tracheal Hook Retractor, double, blunt, 16cm
- E162-120** Tracheal Hook Retractor, double, sharp, 16cm



- E160-410** Cottle Four Prong Retractor, blunt, 14cm



- E160-405** Cottle (Foman) Two Prong Retractor, blunt, 14cm



- E160-420** Cottle Two Prong Retractor, sharp / blunt, left blunt, 14cm



- E160-415** Cottle Two Prong Retractor, sharp / blunt, right blunt, 14cm



- E160-505** Cottle Nievert Retractor, double ended, 19cm



Ear, Nose and Throat

Retractors



E160-605 Cottle Retractor, upper lateral exposing, 16cm


Retractors - Self Retaining



E164-105 Mollison Self Retaining Retractor, 2:2 sharp teeth, child, 13.5cm

E164-110 Mollison Self Retaining Retractor, 4:4 sharp teeth, adult, 16.5cm

Scissors - Nasal



E033-500 Nasal Scissor, straight, 4.5mm serrated blade

E033-502 Nasal Scissor, straight, 10mm serrated blade, adult


E033-504 Nasal Scissor, curved right, 10mm serrated blade, adult

E033-506 Nasal Scissor, curved left, 10mm serrated blade, adult

E033-512 Nasal Scissor, straight, 6mm serrated blade, paediatric

E033-514 Nasal Scissor, curved right, 6mm serrated blade, paediatric

E033-516 Nasal Scissor, curved left, 6mm serrated blade, paediatric




E033-520 Nasal Scissor, rotating 360o, straight, 10mm serrated blade, adult

E033-525 Nasal Scissor, rotating 360o, straight, 6mm serrated blade, paediatric

E044-105 Cawthorne (Bellucci) Aural Scissor, extra fine, straight, 4mm x 1.5mm

E044-110 Cawthorne (Bellucci) Aural Scissor, extra fine, curved to left, 4mm x 1.5mm

E044-115 Cawthorne (Bellucci) Aural Scissor, extra fine, curved to right, 4mm x 1.5mm




E044-205 Ormerod Aural Scissor, extra fine, curved upward, 4mm x 1.5mm



E044-300 Bellucci (Shea) Aural Scissor, straight, 8mm x 1.5mm

E044-305 Bellucci (Shea) Aural Scissor, curved to left, 8mm x 1.5mm


E044-310 Bellucci (Shea) Aural Scissor, curved to right, 8mm x 1.5mm



E048-105 Cawthorne Micro Aural Scissor, straight, 3mm x 0.8mm

E048-110 Cawthorne Micro Aural Scissor, curved to left, 3mm x 0.8mm


E048-115 Cawthorne Micro Aural Scissor, curved to right, 3mm x 0.8mm




Scissors




E172-405 Heymann Nasal Scissor, angled laterally, 18cm



E172-710 Scott Stevenson Tonsil Scissor, 21.5cm



E172-705 Wilson Tonsil Scissor, blunt / blunt, straight, 20cm




E172-905 Wood Tonsil Scissor, triangular blades, curved, 19cm



E172-305 Foman Scissor, blunt / blunt, angled, 14cm

E172-310 Foman Scissor, upper lateral, curved, 14cm


E172-315 Foman Scissor, tungsten carbide blade, blunt / blunt, angled, 14.5cm



E172-505 Knapp Scissor, blunt points, satin finish, curved, 11.5cm

E172-510 Knapp Scissor, blunt points, satin finish, straight, 11.5cm

- E174-205** Lempert Scoop, no.0000, 20cm
- E174-210** Lempert Scoop, no.000, 20cm
- E174-215** Lempert Scoop, no.00, 20cm
- E174-220** Lempert Scoop, no.0, 20cm
- E174-225** Lempert Scoop, no.1, 20cm
- E174-230** Lempert Scoop, no.2, 20cm
- E174-235** Lempert Scoop, set /6



Suction Tubes - Aspiration

- E194-105** Aspiration Tube, curved, 2mm dia, 20cm
- E194-110** Aspiration Tube, curved, 2mm dia, 25cm
- E194-115** Aspiration Tube, curved, 2mm dia, 30cm
- E194-120** Aspiration Tube, curved, 2mm dia, 35cm
- E194-125** Aspiration Tube, curved, 2mm dia, 40cm
- E194-130** Aspiration Tube, curved, 2mm dia, 45cm
- E194-135** Aspiration Tube, curved, 2mm dia, 50cm
- E194-140** Aspiration Tube, curved, 2mm dia, 60cm
- E194-205** Aspiration Tube, curved, 3mm dia, 20cm
- E194-210** Aspiration Tube, curved, 3mm dia, 25cm
- E194-215** Aspiration Tube, curved, 3mm dia, 30cm
- E194-220** Aspiration Tube, curved, 3mm dia, 35cm
- E194-225** Aspiration Tube, curved, 3mm dia, 40cm
- E194-230** Aspiration Tube, curved, 3mm dia, 45cm
- E194-235** Aspiration Tube, curved, 3mm dia, 50cm
- E194-237** Aspiration Tube, curved, 3mm dia, 55cm
- E194-240** Aspiration Tube, curved, 3mm dia, 60cm
- E194-305** Aspiration Tube, curved, 4mm dia, 20cm
- E194-310** Aspiration Tube, curved, 4mm dia, 25cm
- E194-315** Aspiration Tube, curved, 4mm dia, 30cm
- E194-320** Aspiration Tube, curved, 4mm dia, 35cm
- E194-325** Aspiration Tube, curved, 4mm dia, 40cm
- E194-330** Aspiration Tube, curved, 4mm dia, 45cm
- E194-335** Aspiration Tube, curved, 4mm dia, 50cm
- E194-340** Aspiration Tube, curved, 4mm dia, 60cm
- E194-405** Aspiration Tube, curved, 5mm dia, 20cm
- E194-410** Aspiration Tube, curved, 5mm dia, 25cm
- E194-415** Aspiration Tube, curved, 5mm dia, 30cm
- E194-420** Aspiration Tube, curved, 5mm dia, 35cm
- E194-425** Aspiration Tube, curved, 5mm dia, 40cm
- E194-430** Aspiration Tube, curved, 5mm dia, 45cm
- E194-435** Aspiration Tube, curved, 5mm dia, 50cm
- E194-440** Aspiration Tube, curved, 5mm dia, 60cm
- E194-505** Aspiration Tube, curved, 6mm dia, 20cm
- E194-510** Aspiration Tube, curved, 6mm dia, 25cm
- E194-515** Aspiration Tube, curved, 6mm dia, 30cm
- E194-520** Aspiration Tube, curved, 6mm dia, 35cm
- E194-525** Aspiration Tube, curved, 6mm dia, 40cm
- E194-530** Aspiration Tube, curved, 6mm dia, 45cm
- E194-535** Aspiration Tube, curved, 6mm dia, 50cm
- E194-540** Aspiration Tube, curved, 6mm dia, 60cm



Suction Tubes



- E196-083** Barron Suction Tube, 3fg (1mm dia), 13cm
- E196-084** Barron Suction Tube, 4fg (1.3mm dia), 13cm
- E196-085** Barron Suction Tube, 5fg (1.6mm dia), 13cm
- E196-086** Barron Suction Tube, 6fg (2mm dia), 13cm
- E196-087** Barron Suction Tube, 7fg (2.3mm dia), 13cm
- E196-093** Barron Suction Tube, 3fg (1mm dia), 17cm
- E196-094** Barron Suction Tube, 4fg (1.3mm dia), 17cm
- E196-095** Barron Suction Tube, 5fg (1.6mm dia), 17cm
- E196-096** Barron Suction Tube, 6fg (2mm dia), 17cm
- E196-097** Barron Suction Tube, 7fg (2.3mm dia), 17cm
- E196-103** Barron Suction Tube, 3fg (1mm dia), 23cm
- E196-104** Barron Suction Tube, 4fg (1.3mm dia), 23cm
- E196-105** Barron Suction Tube, 5fg (1.6mm dia), 23cm
- E196-106** Barron Suction Tube, 6fg (2mm dia), 23cm
- E196-107** Barron Suction Tube, 7fg (2.3mm dia), 23cm



- E196-205** Bellucci Suction Tube, tube and 3 ends, 20, 22 & 24swg
- E196-210** Bellucci Suction Tube, tube only, (6fg for suction ends), 10cm working length, 14cm overall length




- N150-105** Frazier Suction Tube, 30o angle, 5fg (1.66mm dia), 120mm
- N150-106** Frazier Suction Tube, 30o angle, 6fg (2mm dia), 120mm
- N150-107** Frazier Suction Tube, 30o angle, 7fg (2.3mm dia), 120mm
- N150-108** Frazier Suction Tube, 30o angle, 8fg (2.6mm dia), 120mm
- N150-109** Frazier Suction Tube, 30o angle, 9fg (3mm dia), 120mm
- N150-110** Frazier Suction Tube, 30o angle, 10fg (3.3mm dia), 120mm
- N150-111** Frazier Suction Tube, 30o angle, 11fg (3.6mm dia), 120mm
- N150-112** Frazier Suction Tube, 30o angle, 12fg (4mm dia), 120mm
- N150-115** Frazier Suction Tube, 30o angle, 15fg (5mm dia), 120mm



- N150-205** Fergusson Suction Tube, 90o angle, 5fg (1.66mm dia), 90mm
- N150-206** Fergusson Suction Tube, 90o angle, 6fg (2mm dia), 90mm
- N150-207** Fergusson Suction Tube, 90o angle, 7fg (2.3mm dia), 90mm
- N150-208** Fergusson Suction Tube, 90o angle, 8fg (2.6mm dia), 90mm
- N150-209** Fergusson Suction Tube, 90o angle, 9fg (3mm dia), 90mm
- N150-210** Fergusson Suction Tube, 90o angle, 10fg (3.3mm dia), 90mm
- N150-211** Fergusson Suction Tube, 90o angle, 11fg (3.6mm dia), 90mm
- N150-212** Fergusson Suction Tube, 90o angle, 12fg (4mm dia), 90mm
- N150-215** Fergusson Suction Tube, 90o angle, 15fg (5mm dia), 90mm

E196-275 Formby Suction Tube, angled tapered tube to 10fg (3.3mm)



E196-405 House Suction and Irrigation Tube, 8fg (2.6mm) suction x 6fg (2mm) irrigation


E196-410 House Suction and Irrigation Tube, 10fg (3.3mm) suction x 7fg (2.3mm) irrigation

E196-415 House Suction and Irrigation Tube, 12fg (4mm) suction x 8fg (2.6mm) irrigation

E196-420 House Suction and Irrigation Tube, 16swg (1.63mm) suction x 18swg (1.22mm) irrigation

E196-425 House Suction and Irrigation Tube, 18swg (1.22mm) suction x 20swg (0.9mm) irrigation

E196-430 House Suction and Irrigation Tube, 20swg (0.9mm) suction x 22swg (0.7mm) irrigation



E196-515 Lempert Suction Tube, 6fg, 2mm dia, 120mm working length

E196-520 Lempert Suction Tube, 9fg, 3mm dia, 120mm working length



E196-605 Magill Suction Tube, small, 3mm dia

E196-610 Magill Suction Tube, medium, 3.5mm dia

E196-615 Magill Suction Tube, large, 4mm dia

E196-620 Magill Suction Tube, extra large, 5mm dia

E196-603 Magill Suction Tube, set of 3 sizes, small, medium and large




E196-640 Magill Suction Tube, insulated, small, 3mm dia

E196-645 Magill Suction Tube, insulated, medium, 3.5mm dia

E196-650 Magill Suction Tube, insulated, large, 4mm dia

E196-655 Magill Suction Tube, detachable ends with handle, set/4



E196-672 Micro Laryngeal Suction Tube, with finger control, 2mm dia, effective length 26cm

E196-674 Micro Laryngeal Suction Tube, with finger control, 2.5mm dia, effective length 26cm

E196-676 Micro Laryngeal Suction Tube, with finger control, 3mm dia, effective length 26cm

E196-669 Micro Laryngeal Suction Tube, with finger control, 3.5mm dia, effective length 26cm

E196-670 Micro Laryngeal Suction Tube, with finger control, 4mm dia, effective length 26cm

E196-671 Micro Laryngeal Suction Tube, with finger control, 5mm dia, effective length 26cm



Suction Tubes



E196-710 Rowbotham Suction Tube, 28.5cm



E196-805 Schuknecht Suction Tube, 20swg (0.9mm)

E196-810 Schuknecht Suction Tube, 24swg (0.55mm)



E196-815 Self - Clearing Suction Tube, 1.5mm dia

E196-820 Self - Clearing Suction Tube, 2.0mm dia

E196-825 Self - Clearing Suction Tube, 3.0mm dia



E196-916 Verhoeven Suction Tube, 16swg (1.63mm)


E196-918 Verhoeven Suction Tube, 18swg (1.22mm)

E196-920 Verhoeven Suction Tube, 20swg (0.9mm)

E196-922 Verhoeven Suction Tube, 22swg (0.7mm)

E196-924 Verhoeven Suction Tube, 24swg (0.55mm)

E196-900 Verhoeven Suction Tube, set /5
(16,18, 20, 22 and 24swg)




E196-930 House (Fisch) Suction Tube, cut out adaptor,
5.5cm

Suction Tubes

- E196-990** Suction Tip (For Zoellner And Bellucci Tubes), detachable end, 18swg (1.22mm)
- E196-991** Suction Tip (For Zoellner And Bellucci Tubes), detachable end, 20swg (0.9mm)
- E196-992** Suction Tip (For Zoellner And Bellucci Tubes), detachable end, 22swg (0.7mm)
- E196-993** Suction Tip (For Zoellner And Bellucci Tubes), detachable end, 24swg (0.55mm)
- E196-994** Suction Tip (For Zoellner And Bellucci Tubes), detachable end, 26swg (0.45mm)




- E196-975** Zoellner Suction Tube, (original pattern) 16cm




Syringes


- E184-405** Bacon Aural Syringe, complete



- E184-105** Aural Syringe, 2oz (50ml), complete with 2 pipes
- E184-110** Aural Syringe, 3oz (75ml), complete with 2 pipes
- E184-115** Aural Syringe, 4oz (100ml), complete with 2 pipes
- E184-120** Aural Syringe, 6oz (150ml), complete with 2 pipes




- E186-105** Dental Chip Syringe, complete
- E186-115** Dental Chip Syringe, rubber bulb only



Ear, Nose and Throat


Tongue Depressors



- E024-010** Lack Tongue Depressor, child, 13mm
- E024-015** Lack Tongue Depressor, adolescent, 19mm
- E024-020** Lack Tongue Depressor, adult, 22mm




- E024-003** Flat (Buchwald) Tongue Depressor, 16mm, 14.5cm
- E024-005** Flat (Buchwald) Tongue Depressor, 23mm tapering to 18mm, 20cm




- E024-500** Tongue Depressor, folding

Tongue Forceps



- E080-205** Guy Tongue Holding Forcep, child, 12.5cm
- E080-210** Guy Tongue Holding Forcep, adult, 18cm


Tonsil Holding Forceps



- E082-030** White Tonsil Seizing Forcep, curved, 18cm
- E082-033** White Tonsil Seizing Forcep, curved, 23cm

Tonsil Holding Forceps

E082-105 Denis Browne Tonsil Holding Forcep, angled, 20cm



Tonsil Snare

E180-105 Eve Tonsil Snare, stainless steel, slide action, 28cm



Tracheal Dilators

E056-205 Trousseau Tracheal Dilating Forcep, child, 10.5cm




E056-210 Trousseau Tracheal Dilating Forcep, adult, 14cm




Trocars

E190-100 Tilley Harpoon Trocar, 17cm




Ear, Nose and Throat

Trocars



- E190-105** Tilley Lichtwitz Trocar, luer lock, with tubing mount, 5fg
- E190-106** Tilley Lichtwitz Trocar, luer lock, with tubing mount, 6fg
- E190-107** Tilley Lichtwitz Trocar, luer lock, with tubing mount, 7fg
- E190-108** Tilley Lichtwitz Trocar, luer lock, with tubing mount, 8fg
- E190-109** Tilley Lichtwitz Trocar, luer lock, with tubing mount, 9fg

Turbinate Forceps



- E066-105** Grunwald Nasal Turbinate Punch Forcep, pattern A, small, 11.5cm
- E066-110** Grunwald Nasal Turbinate Punch Forcep, pattern B, large, 11.5cm
- E066-115** Grunwald Nasal Turbinate Punch Forcep, pattern C, large, fenestrated, 11.5cm



- E066-205** Irwin Moore Nasal Turbinate Forcep, 19.5cm



- E066-305** Luc Nasal Turbinate Forcep, 6mm wide heart shaped jaw, small, 19.5cm
- E066-310** Luc Nasal Turbinate Forcep, 7.5mm wide heart shaped jaw, medium, 19.5cm
- E066-315** Luc Nasal Turbinate Forcep, 9mm wide heart shaped jaw, large, 19.5cm


Turbinate Forceps

- E066-405** Luc Nasal Turbinate Forcep, 6mm wide oval shaped jaw, small, 19.5cm
- E066-410** Luc Nasal Turbinate Forcep, 7.5mm wide oval shaped jaw, medium, 19.5cm
- E066-415** Luc Nasal Turbinate Forcep, 9mm wide oval shaped jaw, large, 19.5cm



Wool Carriers


- E150-210** Heath Wool Carrier, 17cm



- E150-605** Musgrove Hook / Wool Carrier, double ended, 18cm



- E150-705** Tumarkin Wool Carrier, twisted wire, pk/10, 15cm




- E150-305** Jobson Horne Probe, double ended, 14cm

- E150-310** Jobson Horne Probe, double ended, 18cm



Various

- E149-250** King Vein Graft Forcep, round flat blades, 12.5cm



- E096-305** Gardiner Brown Hook and Imray Loop, double ended, 12.5cm



- E083-105** Gardiner Brown Tuning Fork, C0, 128dv
- E083-110** Gardiner Brown Tuning Fork, C1, 256dv
- E083-115** Gardiner Brown Tuning Fork, C2, 512dv
- E083-120** Gardiner Brown Tuning Fork, C3, 1024dv
- E083-125** Gardiner Brown Tuning Fork, C4, 2048dv
- E083-130** Gardiner Brown Tuning Fork, C5, 4096dv





Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk

We are an Approved Supplier to the
NHS Supply Chain



Supply Chain
v1