

Judd
medical

General Surgery Instrumentation Catalogue

www.judd-medical.co.uk

REUSABLE

We are an Approved
Supplier to the
NHS Supply Chain

NHS

Supply Chain


Index



	Pages		Pages
Artery Forceps	1-7	Probes	30
Cholecystectomy Forceps	7	Post Mortem	30
Clamps	8-10	Retractors	31-41
Clamps - Atraumatic	10-11	Retractors - Self Retaining	41-56
Clip Forceps	11-12	Saws	56
Decompressors	12	Scalpel Handles	56-57
Diathermy Quivers	12	Scissors	57-64
Dilators	13	Scissors - Tungsten Carbide	65-67
Directors	13	- Super Cut	67-69
Dissecting Forceps	14-17	Sinus Forceps	69
Dissectors	18	Splinter Forceps	69
Dressing Forceps	18	Sponge Holding Forceps	70
Hernia Forceps	19	Steriliser Forceps	70
Instrument Pins	19	Suction Tubes	71-72
Intestinal Forceps	19	Suture Forceps	72
Introducing Forceps	19	Tissue Forceps	72-74
Knives	20	Tongue Forceps	74
Ligature Forceps	7	Towel Clips	74-75
Measures	20	Bags Clips	75
Mouth Gags	20	Vein Strippers	76
Needleholders	21-27	Vasectomy	76
Needles	27-28	Holloware	77
Percussors	28		
Plaster Instruments	28-29		
Ring Cutters	29		

Artery Forceps

- G052-002** Adson Fraser Artery Forcep, straight, 18cm
G052-004 Adson Fraser Artery Forcep, curved, 18cm
G052-006 Adson Fraser Artery Forcep, straight, 1:2 teeth, 18cm
G052-008 Adson Fraser Artery Forcep, curved, 1:2 teeth, 18cm



- G052-062** Dieffenbach Artery Forcep, straight, 3.5cm
G052-064 Dieffenbach Artery Forcep, straight, 5cm
G052-066 Dieffenbach Artery Forcep, curved, 3.5cm
G052-068 Dieffenbach Artery Forcep, curved, 5cm



- G052-072** Dunhill Artery Forcep, straight, 12.5cm
G052-076 Dunhill Artery Forcep, straight, 18cm


- G052-074** Dunhill Artery Forcep, curved, 12.5cm
G052-078 Dunhill Artery Forcep, curved, 18cm



- G052-090** Fraser Kelly Artery Forcep, curved, 18cm




- G052-094** Grey Turner Artery Forcep, curved, 14.5cm
G052-098 Grey Turner Artery Forcep, curved, 18cm



General

Artery Forceps



G052-092 Grey Turner Artery Forcep, straight, 14.5cm

G052-096 Grey Turner Artery Forcep, straight, 18cm



G052-100 Halstead Mosquito Artery Forcep, straight, 10cm

G052-104 Halstead Mosquito Artery Forcep, straight, 12.5cm

G052-102 Halstead Mosquito Artery Forcep, curved, 10cm

G052-106 Halstead Mosquito Artery Forcep, curved, 12.5cm

G052-107 Halstead Mosquito Artery Forcep, straight, 1:2 teeth, 12.5cm

G052-109 Halstead Mosquito Artery Forcep, curved, 1:2 teeth, 12.5cm



G052-099 Halstead Mosquito Artery Forcep, straight, fine, 10cm

G052-103 Halstead Mosquito Artery Forcep, straight, fine, 12cm

G052-101 Halstead Mosquito Artery Forcep, curved, fine, 10cm

G052-105 Halstead Mosquito Artery Forcep, curved, fine, 12cm



G052-112 Heiss Artery Forcep, straight, 20cm

G052-114 Heiss Artery Forcep, curved, 20cm

G052-115 Heiss Artery Forcep, fully curved, 19.5cm

G052-116 Heiss Artery Forcep, straight, 1:2 teeth, 20cm

G052-118 Heiss Artery Forcep, curved, 1:2 teeth, 20cm



G052-122 Kelly Artery Forcep, straight, 14cm


G052-126 Kelly Artery Forcep, straight, 16cm

G052-124 Kelly Artery Forcep, curved, 14cm


G052-128 Kelly Artery Forcep, curved, 16cm

Artery Forceps


- G052-142** Kilner Artery Forcep, straight, 14cm
G052-144 Kilner Artery Forcep, curved, 14cm
G052-147 Kilner Artery Forcep, straight, 1:2 teeth, 14cm
G052-149 Kilner Artery Forcep, curved, 1:2 teeth, 14cm




- G052-152** Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 12.5cm
G052-156 Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 14cm
G052-162 Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 16cm
G052-166 Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 18cm
G052-172 Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 20cm
G052-176 Kocher (Ochsner) Artery Forcep, straight, 1:2 teeth, 22cm




- G052-154** Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 12.5cm
G052-158 Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 14cm
G052-164 Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 16cm
G052-168 Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 18cm
G052-174 Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 20cm
G052-178 Kocher (Ochsner) Artery Forcep, curved, 1:2 teeth, 22cm



- G052-202** Mayo Oschner Artery Forcep, 1:2 teeth, straight, 18.5cm
G052-204 Mayo Oschner Artery Forcep, 1:2 teeth, curved, 18.5cm



- G052-190** Mayo Artery Forcep, curved, 15cm



General

Artery Forceps



G052-500 Halstead Mosquito Artery Forcep, straight, 14cm

G052-510 Halstead Mosquito Artery Forcep, straight, 18cm

G052-502 Halstead Mosquito Artery Forcep, curved, 14cm

G052-515 Halstead Mosquito Artery Forcep, curved, 18cm

G052-498 Halstead Mosquito Artery Forcep, 1:2 teeth, straight, 10cm

G052-504 Halstead Mosquito Artery Forcep, 1:2 teeth, straight, 14cm

G052-499 Halstead Mosquito Artery Forcep, 1:2 teeth, curved, 10cm

G052-506 Halstead Mosquito Artery Forcep, 1:2 teeth, curved, 14cm




G052-213 Mixer Artery Forcep, curved, fine jaw, 14cm

G052-215 Mixer Artery Forcep, curved, fine jaw, 18.5cm




G052-217 Mixer Artery Forcep, curved, fine jaw, 23cm



G052-222 Moynihan Artery Forcep, straight, 14.5cm

G052-226 Moynihan Artery Forcep, straight, 18cm




G052-224 Moynihan Artery Forcep, curved, 14.5cm

G052-228 Moynihan Artery Forcep, curved, 18cm

Artery Forceps

- G052-281** O'Shaughnessy Artery Forcep, curved, 15cm
G052-283 O'Shaughnessy Artery Forcep, curved, 18cm
G052-285 O'Shaughnessy Artery Forcep, curved, 20cm
G052-287 O'Shaughnessy Artery Forcep, curved, 23cm
G052-289 O'Shaughnessy Artery Forcep, curved, 25cm



- G052-110** Harrison Cripp Artery Forcep, 17.5cm



- G052-010** Baby Adson Artery Forcep, curved, 14cm




- G052-302** Rochester Pean Artery Forcep, straight, 14cm
G052-306 Rochester Pean Artery Forcep, straight, 16cm
G052-312 Rochester Pean Artery Forcep, straight, 18cm
G052-322 Rochester Pean Artery Forcep, straight, 20cm
G052-325 Rochester Pean Artery Forcep, straight, 22cm
G052-326 Rochester Pean Artery Forcep, straight, 24cm
G052-332 Rochester Pean Artery Forcep, straight, 26cm
G052-336 Rochester Pean Artery Forcep, straight, 30cm

- G052-304** Rochester Pean Artery Forcep, curved, 14cm
G052-308 Rochester Pean Artery Forcep, curved, 16cm
G052-314 Rochester Pean Artery Forcep, curved, 18cm
G052-324 Rochester Pean Artery Forcep, curved, 20cm
G052-327 Rochester Pean Artery Forcep, curved, 22cm
G052-328 Rochester Pean Artery Forcep, curved, 24cm
G052-334 Rochester Pean Artery Forcep, curved, 26cm
G052-338 Rochester Pean Artery Forcep, curved, 30cm




General


Artery Forceps




M062-220 Ronald Edwards Artery Forcep, screw joint, shape A, 21cm



M062-222 Ronald Edwards Artery Forcep, screw joint, shape B, 21cm




M062-224 Ronald Edwards Artery Forcep, screw joint, shape C, 21cm



- G052-372** Spencer Wells Artery Forcep, straight, 12.5cm
- G052-382** Spencer Wells Artery Forcep, straight, 15cm
- G052-392** Spencer Wells Artery Forcep, straight, 18cm
- G052-402** Spencer Wells Artery Forcep, straight, 20cm
- G052-412** Spencer Wells Artery Forcep, straight, 23cm

Artery Forceps

- G052-374** Spencer Wells Artery Forcep, curved, 12.5cm
- G052-384** Spencer Wells Artery Forcep, curved, 15cm
- G052-394** Spencer Wells Artery Forcep, curved, 18cm
- G052-404** Spencer Wells Artery Forcep, curved, 20cm
- G052-414** Spencer Wells Artery Forcep, curved, 23cm




Ligature Forceps

- G056-575** Wickstrom Ligature Forcep, 20cm



Cholecystectomy Forceps

- G056-117** Geary Grant Cholecystectomy Forcep, large curve, longitudinal serrations, 23cm
- G056-115** Geary Grant Cholecystectomy Forcep, small curve, longitudinal serrations, 23cm



- G056-119** Lahey Cholecystectomy Forcep, curved, 19.5cm
- G056-123** Lahey Cholecystectomy Forcep, curved, 23cm




- G056-423** Moynihan Cholecystectomy Forcep, curved, 23cm



General

Clamps



G017-010 De Martel Anastomosis Clamp, complete




G036-093 Doyen Intestinal Clamp, child, serrated jaw, curved, 16.5cm

G036-097 Doyen Intestinal Clamp, adult, serrated jaw, curved, 24cm




G036-091 Doyen Intestinal Clamp, child, serrated jaw, straight, 16.5cm


G036-095 Doyen Intestinal Clamp, adult, serrated jaw, straight, 24cm



G017-020 Galbraith (Modified Wade) kidney pedicle clamp, curved, 25cm



G017-044 Price Thomas Bronchus Clamp, serrated jaws, double angle, 22cm



G036-110 Kocher Intestinal Clamp, curved, longitudinal serrations, 22cm

G036-114 Kocher Intestinal Clamp, curved, longitudinal serrations, 25.5cm

G036-118 Kocher Intestinal Clamp, curved, longitudinal serrations, 28cm

G036-108 Kocher Intestinal Clamp, straight, longitudinal serrations, 22cm

G036-112 Kocher Intestinal Clamp, straight, longitudinal serrations, 25.5cm


G036-116 Kocher Intestinal Clamp, straight, longitudinal serrations, 28cm

G036-210 Lane Intestinal Clamp, straight, longitudinal serrations, 21cm




G036-555 Lanes Twin Anastomosis Clamp, curved, longitudinal serrations, 30.5cm

G036-550 Lanes Twin Anastomosis Clamp, straight, longitudinal serrations, 30.5cm



G036-310 Lang Stevenson Intestinal Clamp, keyed blades, straight, 20cm




G036-482 Parker Kerr Intestinal Clamp, curved, longitudinal serrations, 23cm

G036-486 Parker Kerr Intestinal Clamp, curved, longitudinal serrations, 25.5cm

G036-480 Parker Kerr Intestinal Clamp, straight, longitudinal serrations, 23cm

G036-484 Parker Kerr Intestinal Clamp, straight, longitudinal serrations, 25.5cm



G017-032 Payr Intestinal Clamp, no pin, 20cm, 7cm blades

G017-034 Payr Intestinal Clamp, no pin, 25.5cm, 10cm blades




General

Clamps



G036-500 Schoemaker-Ogilvie Colectomy Clamp, straight, hollow blades, 17cm


G036-510 Schoemaker-Ogilvie Colectomy Clamp, straight, hollow blades, 21.5cm



G096-105 Tube Occluding Clamp, heavy wedge serrated jaw, 15cm

G096-110 Tube Occluding Clamp, heavy wedge serrated jaw, 18cm

G096-115 Tube Occluding Clamp, heavy wedge serrated jaw, 20cm



G017-050 Zachary Cope Anastomosis Clamp, complete


Clamps - Atraumatic




G036-092 Doyen Intestinal Clamp, child, atraumatic, curved, 17cm

Clamps - Atraumatic


G036-096 Doyen Intestinal Clamp, adult, atraumatic, curved, 24cm



G036-090 Doyen Intestinal Clamp, child, atraumatic, straight, 17cm



G036-094 Doyen Intestinal Clamp, adult, atraumatic, straight, 24cm




Clip Forceps


G058-120 Comyn Berkeley Inserting Clip Forcep, 14cm



G058-175 Michel (Modified) Clip Applying And Removing Forcep, Spring Handle




G058-110 Childes Clip Inserting Forcep, up/rack, 18cm



General

Clip Forceps



G058-150 Michel Clip Extracting Forcep, screw joint

Decompressors



G038-210 Savage Intestinal Decompressor, adult, 48cm

G038-215 Savage Intestinal Decompressor, child, 34cm

Diathermy Quivers



G022-905 Diathermy Quiver, 4.8cm dia x 24.5cm

G022-910 Diathermy Quiver, 6.0cm dia x 23cm



G022-915 Julian Taylor Diathermy Quiver, with drainage holes, 4.4cm x 25cm



G022-950 Diathermy Quiver, 6.0cm dia x 37.5cm




G022-975 Incus Diathermy Quiver, with removable screw bottom, 6cm dia x 24cm


G022-985 Incus Diathermy Quiver, with removable screw bottom, 6cm dia x 37.5cm

Dilators


- G040-391** Bake Bile Duct Dilator, malleable stems, 20cm, 3mm
- G040-392** Bake Bile Duct Dilator, malleable stems, 20cm, 4mm
- G040-393** Bake Bile Duct Dilator, malleable stems, 20cm, 5mm
- G040-394** Bake Bile Duct Dilator, malleable stems, 20cm, 6mm
- G040-395** Bake Bile Duct Dilator, malleable stems, 20cm, 7mm
- G040-396** Bake Bile Duct Dilator, malleable stems, 20cm, 8mm
- G040-397** Bake Bile Duct Dilator, malleable stems, 20cm, 9mm
- G040-398** Bake Bile Duct Dilator, malleable stems, 20cm, 10mm
- G040-399** Bake Bile Duct Dilator, malleable stems, 20cm, 11mm
- G040-390** Bake Bile Duct Dilator, malleable stems, set /9, 3mm to 11mm



- G040-422** Clarke Bile Duct Dilator, stainless steel, double ended, spiral wire stem, 23cm, 3/4mm
- G040-424** Clarke Bile Duct Dilator, stainless steel, double ended, spiral wire stem, 23cm, 5/6mm
- G040-426** Clarke Bile Duct Dilator, stainless steel, double ended, spiral wire stem, 23cm, 7/8mm



- G040-534** Maingot Bile Duct Dilator, double ended, 3mm / 4mm
- G040-556** Maingot Bile Duct Dilator, double ended, 5mm / 6mm
- G040-578** Maingot Bile Duct Dilator, double ended, 7mm / 8mm
- G040-500** Maingot Bile Duct Dilator, double ended, set of 3



Directors

- G042-005** Brodie Director, malleable, 14cm
- G042-010** Brodie Director, malleable, 16.5cm
- G042-015** Brodie Director, malleable, 20cm



- G042-120** Grooved Director, straight, 11.5cm
- G042-122** Grooved Director, straight, 13cm
- G042-124** Grooved Director, straight, 14.5cm
- G042-126** Grooved Director, straight, 16cm
- G042-128** Grooved Director, straight, 18cm
- G042-129** Grooved Director, straight, 20cm



- G042-205** Key Hernia Director, 15cm



General

Dissecting Forceps



- G062-112** Adson Dissecting Forcep, 1:2 teeth, 12cm
- G062-113** Adson Dissecting Forcep, 1:2 teeth, 15cm
- G062-114** Adson Dissecting Forcep, 1:2 teeth, micro, 12cm
- G062-115** Adson Dissecting Forcep, 1:2 teeth, micro, 15cm




- G062-102** Adson Dissecting Forcep, serrated, 12cm
- G062-103** Adson Dissecting Forcep, serrated, 15cm
- G062-104** Adson Dissecting Forcep, serrated, micro, 12cm
- G062-105** Adson Dissecting Forcep, serrated, micro, 15cm



- G062-062** Standard Thumb Dissecting Forcep, serrated, straight, 12.5cm
- G062-064** Standard Thumb Dissecting Forcep, serrated, straight, 14cm
- G062-066** Standard Thumb Dissecting Forcep, serrated, straight, 15cm
- G062-068** Standard Thumb Dissecting Forcep, serrated, straight, 16cm
- G062-070** Standard Thumb Dissecting Forcep, serrated, straight, 18cm
- G062-072** Standard Thumb Dissecting Forcep, serrated, straight, 20cm
- G062-074** Standard Thumb Dissecting Forcep, serrated, straight, 25cm
- G062-076** Standard Thumb Dissecting Forcep, serrated, straight, 30cm



- G062-063** Standard Tissue Dissecting Forcep, 1:2 teeth straight, 12.5cm
- G062-065** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 14cm
- G062-067** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 15cm
- G062-069** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 16cm
- G062-071** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 18cm
- G062-073** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 20cm
- G062-075** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 25cm
- G062-077** Standard Tissue Dissecting Forcep, 1:2 teeth, straight, 30cm




- G062-040** Block End Dissecting Forcep, serrated, 12.5cm
- G062-042** Block End Dissecting Forcep, serrated, 15cm
- G062-044** Block End Dissecting Forcep, serrated, 18cm
- G062-046** Block End Dissecting Forcep, serrated, 20cm
- G062-048** Block End Dissecting Forcep, serrated, 25cm
- G062-049** Block End Dissecting Forcep, serrated, 30cm



- G062-137** Bickford Dissecting Forcep, 1:2 teeth, 23cm


Dissecting Forceps

G062-136 Bickford Dissecting Forcep, serrated, 23cm




G062-135 Bonney Dissecting Forcep, 1:2 teeth, 18cm

G062-130 Bonney Dissecting Forcep, serrated, 18cm




G062-139 Braithwaite Dissecting Forcep, 1:2 teeth, 14cm



G062-154 Canadian Pattern Dissecting Forcep, 2:3 teeth, 14cm

G062-156 Canadian Pattern Dissecting Forcep, 2:3 teeth, 15cm



G062-150 Canadian Pattern Dissecting Forcep, serrated, 14cm

G062-152 Canadian Pattern Dissecting Forcep, serrated, 15cm



G062-922 USA Pattern Dissecting Forcep, 1:2 teeth, 12.5cm

G062-924 USA Pattern Dissecting Forcep, 1:2 teeth, 14cm


G062-926 USA Pattern Dissecting Forcep, 1:2 teeth, 15cm

G062-928 USA Pattern Dissecting Forcep, 1:2 teeth, 18cm

G062-930 USA Pattern Dissecting Forcep, 1:2 teeth, 20cm

G062-932 USA Pattern Dissecting Forcep, 1:2 teeth, 25cm

G062-934 USA Pattern Dissecting Forcep, 1:2 teeth, 30cm




G062-230 Emmett (Moynihan) Dissecting Forcep, serrated, 20cm



General


Dissecting Forceps



G062-268 Gillies Dissecting Forcep, 1:2 teeth, 15cm



G062-266 Gillies Dissecting Forcep, serrated, 15cm




G062-272 Jeans Dissecting Forcep, 1:2 teeth, 14cm

G062-274 Jeans Dissecting Forcep, 1:2 teeth, 18cm

G062-276 Jeans Dissecting Forcep, 2:3 teeth, 18cm

G062-270 Jeans Dissecting Forcep, serrated, 14cm




G062-300 Lane Dissecting Forcep, 1:2 teeth, 12.5cm

G062-302 Lane Dissecting Forcep, 1:2 teeth, 15cm

G062-308 Lane Dissecting Forcep, 1:2 teeth, 18cm


G062-309 Lane Dissecting Forcep, 1:2 teeth, 20cm



G062-304 Lane Dissecting Forcep, 2:3 teeth, 15cm

G062-310 Lane Dissecting Forcep, 2:3 teeth, 18cm

G062-311 Lane Dissecting Forcep, 2:3 teeth, 20cm



G062-320 Leedham Green Dissecting Forcep, 2:3 teeth, 14.5cm



G062-322 Leedham Green Dissecting Forcep, 3:4 teeth, 14.5cm




G062-330 Maingot Dissecting Forcep, serrated, 28cm

Dissecting Forceps

G062-345 McIndoe Dissecting Forcep, 1:2 teeth, 15cm



G062-340 McIndoe Dissecting Forcep, serrated, 15cm



G062-432 Ramsey (Ferris Smith) Dissecting Forcep, 1:2 teeth, 18cm




G062-434 Ramsey (Ferris Smith) Dissecting Forcep, 2:3 teeth, 18cm

G062-430 Ramsey (Ferris Smith) Dissecting Forcep, serrated, 18cm



G062-452 St. Marks Dissecting Forcep, 1:2 teeth, 29cm

G062-450 St. Marks Dissecting Forcep, serrated, 29cm




G062-484 Treves Dissecting Forcep, 1:2 teeth, 12.5cm

G062-486 Treves Dissecting Forcep, 1:2 teeth, 15cm

G062-488 Treves Dissecting Forcep, 1:2 teeth, 18cm

G062-490 Treves Dissecting Forcep, 1:2 teeth, 20cm

G062-492 Treves Dissecting Forcep, 1:2 teeth, 25cm



G062-020 Turned Over End Dissecting forcep, 12.5cm

G062-022 Turned Over End Dissecting forcep, 15cm

G062-024 Turned Over End Dissecting forcep, 18cm



G062-512 Waugh Dissecting Forcep, 1:2 teeth, 15cm

G062-514 Waugh Dissecting Forcep, 1:2 teeth, 18cm

G062-516 Waugh Dissecting Forcep, 1:2 teeth, 20cm

G062-502 Waugh Dissecting Forcep, serrated, 15cm

G062-504 Waugh Dissecting Forcep, serrated, 18cm

G062-506 Waugh Dissecting Forcep, serrated, 20cm




General


Dissectors



E028-205 Durham Tonsil Dissector, double ended, 19cm




G042-405 Kocher Thyroid Dissector, slight curved, without eye, 16cm



G042-410 Kocher Thyroid Dissector, slight 11mm, curved, with eye, 15cm




G042-415 Kocher Thyroid Dissector, heavy 12mm, curved, with eye, 16cm



G042-500 MacDonald Dissector, double ended, 19cm

G042-510 MacDonald Dissector, double ended, 25cm




G042-605 Watson Cheyne Dissector, double ended, one probe end, 12.5cm

G042-610 Watson Cheyne Dissector, double ended, one probe end, 18cm

Dressing Forceps




G064-200 Bryant Dressing Forcep, 12.5cm



G064-310 French Pattern Dressing Forcep, 12.5cm

Hernia Forceps

G068-100 Collingwood-Stewart Hernia Ring Forcep




Instrument Pins

G128-090 Mayo Instrument Pin, 11.5cm

G128-095 Mayo Instrument Pin, 15cm

G128-100 Mayo Instrument Pin, double, 12.5cm



Intestinal Forceps

G062-215 Denis Browne Dissecting (Intestinal) Forcep, fenestrated jaw, 18cm



Introducing Forceps

A090-050 Magill Introducing Forcep, adult, right hand, 25cm

A090-100 Magill Introducing Forcep, child, right hand, 20cm

A090-150 Magill Introducing Forcep, infant, right hand, 15cm



General

Knives



- G114-120** Brain Knife, solid forged, 20cm blade
- G114-125** Brain Knife, solid forged, 25cm blade
- G114-130** Brain Knife, solid forged, 30cm blade
- G114-135** Brain Knife, solid forged, 38cm blade



- G114-207** Cartilage Knife, solid forged, 7cm blade
- G114-209** Cartilage Knife, solid forged, 9cm blade
- G114-211** Cartilage Knife, solid forged, 10cm blade
- G114-213** Cartilage Knife, solid forged, 13cm blade



- G114-086** Liston Amputation Knife, solid forged, 13cm blade
- G114-088** Liston Amputation Knife, solid forged, 16cm blade
- G114-090** Liston Amputation Knife, solid forged, 19cm blade



- G138-080** Solid Forged Scalpel, 2.5cm blade
- G138-082** Solid Forged Scalpel, 3.0cm blade
- G138-084** Solid Forged Scalpel, 3.5cm blade
- G138-086** Solid Forged Scalpel, 4.0cm blade
- G138-088** Solid Forged Scalpel, 4.5cm blade
- G138-090** Solid Forged Scalpel, 5.0cm blade
- G138-092** Solid Forged Scalpel, 5.5cm blade

Measures




- T600-250** Metal Ruler, 15cm
- T600-300** Metal Ruler, 30cm

Mouth Gags




- E086-205** Fergusson Mouth Gag, Ackland jaw, child, 17cm
- E086-210** Fergusson Mouth Gag, Ackland jaw, adult, 20cm



- E086-105** Doyen Mouth Gag with cam rack, child, 11.5cm
- E086-110** Doyen Mouth Gag with cam rack, adult, 12.5cm

Needle Holders


- G124-052** Bozemann Needle Holder, 20cm
G124-054 Bozemann Needle Holder, 23cm
G124-056 Bozemann Needle Holder, 25.5cm



- G124-070** Crile Wood Needle Holder, 15cm
G124-072 Crile Wood Needle Holder, 18cm
G124-074 Crile Wood Needle Holder, 20cm
G124-076 Crile Wood Needle Holder, 23cm
G124-078 Crile Wood Needle Holder, 25cm
G124-079 Crile Wood Needle Holder, 30cm




- G124-066** Crile Murray Needle Holder, 14.5cm



- G124-090** Derf Needle Holder, 12cm




- G124-168** Nievert Needle Holder, 12.5cm




General

Needle Holders

G124-102 Halsey Needle Holder, fine serrated jaw, 12.5cm




G124-122 Masson (Mayo) Needle Holder, cross serrated jaw, 27cm




G124-112 Kilner Needle Holder, 13.5cm

G124-114 Kilner Needle Holder, 16.5cm




G124-120 Lawrance Needle Holder, lightweight, serrated jaws, 15cm




G124-030 Baumgartner Needle Holder, 12.5cm

G124-032 Baumgartner Needle Holder, 14cm




- G126-041** Olsen Hegar Needle Holder & Scissor, 14cm
G126-043 Olsen Hegar Needle Holder & Scissor, 16.5cm
G126-045 Olsen Hegar Needle Holder & Scissor, 19cm




- G124-126** Mathieu Needle Holder, 14cm
G124-128 Mathieu Needle Holder, 17cm
G124-129 Mathieu Needle Holder, 20cm



- G124-104** Heaney Needle Holder, curved jaw, 16cm
G124-106 Heaney Needle Holder, curved jaw, 21cm
G124-108 Heaney Needle Holder, curved jaw, 25.5cm



- G124-140** Mayo Needle Holder, 12.5cm
G124-142 Mayo Needle Holder, 14cm
G124-144 Mayo Needle Holder, 15cm
G124-146 Mayo Needle Holder, 16cm
G124-148 Mayo Needle Holder, 18cm
G124-152 Mayo Needle Holder, 20cm
G124-154 Mayo Needle Holder, 23cm
G124-156 Mayo Needle Holder, 25cm
G124-158 Mayo Needle Holder, 30cm




General

Needle Holders




G124-160 Naunton Morgan Needle Holder, compound action, 18cm

G124-162 Naunton Morgan Needle Holder, compound action, 25cm



G124-190 Thomson Walker Needle Holder, 20cm



G124-660 Ward Needle Holder, for right hand, 18cm

Needle Holders - Tungsten Carbide



G125-070 Crile Wood Needle Holder, tungsten carbide jaw, 15cm

G125-072 Crile Wood Needle Holder, tungsten carbide jaw, 18cm

G125-074 Crile Wood Needle Holder, tungsten carbide jaw, 20cm

G125-076 Crile Wood Needle Holder, tungsten carbide jaw, 23cm

G125-078 Crile Wood Needle Holder, tungsten carbide jaw, 25cm

G125-079 Crile Wood Needle Holder, tungsten carbide jaw, 30cm

Needle Holders - Tungsten Carbide


G125-080 DeBakey Needle Holder, tungsten carbide jaw, 18cm

G125-082 DeBakey Needle Holder, tungsten carbide jaw, 20cm


G125-084 DeBakey Needle Holder, tungsten carbide jaw, 23cm

G125-086 DeBakey Needle Holder, tungsten carbide jaw, 26cm

G125-088 DeBakey Needle Holder, tungsten carbide jaw, 30cm




G125-090 Derf Needle Holder, tungsten carbide jaw, 12cm




G125-092 Finochietto Needle Holder, tungsten carbide jaw, 20cm

G125-094 Finochietto Needle Holder, tungsten carbide jaw, 26cm




G125-102 Halsey Needle Holder, serrated tungsten carbide jaw, 12.5cm




General

Needle Holders - Tungsten Carbide




G125-112 Kilner Needle Holder, tungsten carbide jaw, 13.5cm

G125-114 Kilner Needle Holder, tungsten carbide jaw, 16.5cm



G125-120 Lawrance Needle Holder, tungsten carbide jaw, 15cm



G125-140 Mayo Hegar Needle Holder, tungsten carbide jaw, 12.5cm

G125-142 Mayo Hegar Needle Holder, tungsten carbide jaw, 14cm

G125-144 Mayo Hegar Needle Holder, tungsten carbide jaw, 15cm

G125-146 Mayo Hegar Needle Holder, tungsten carbide jaw, 16cm


G125-148 Mayo Hegar Needle Holder, tungsten carbide jaw, 18cm

G125-152 Mayo Hegar Needle Holder, tungsten carbide jaw, 20cm

G125-154 Mayo Hegar Needle Holder, tungsten carbide jaw, 24cm

G125-156 Mayo Hegar Needle Holder, tungsten carbide jaw, 26cm


G125-158 Mayo Hegar Needle Holder, tungsten carbide jaw, 30cm



G125-168 Neivert (Masing) Needle Holder, tungsten carbide jaw, 12.5cm

Needle Holders - Tungsten Carbide

- G125-170** Ryder Needle Holder, tungsten carbide jaw, 12.5cm
- G125-171** Ryder Needle Holder, tungsten carbide jaw, 15cm
- G125-173** Ryder Needle Holder, tungsten carbide jaw, 18cm
- G125-175** Ryder Needle Holder, tungsten carbide jaw, 20cm
- G125-177** Ryder Needle Holder, tungsten carbide jaw, 23cm
- G125-179** Ryder Needle Holder, tungsten carbide jaw, 25cm
- G125-180** Micro Vascular Ryder Needle Holder, tungsten carbide serrated, jaw 1mm wide, 12.5cm
- G125-181** Micro Vascular Ryder Needle Holder, tungsten carbide serrated, jaw 1mm wide, 15cm
- G125-183** Micro Vascular Ryder Needle Holder, tungsten carbide serrated, jaw 1mm wide, 18cm
- G125-185** Micro Vascular Ryder Needle Holder, tungsten carbide serrated, jaw 1mm wide, 20cm
- G125-187** Micro Vascular Ryder Needle Holder, tungsten carbide serrated, jaw, 1mm wide, 23cm




Ligature Needles

- G118-045** Ligature Needle, small curve, 15cm



- G122-500** Bonney Reverdin Set, comprising, needle holder and 3 needles



- G118-310** Cooper Ligature Needle, 20cm



- G118-020** Deschamp Ligature Needle, blunt, for left hand, 21cm

- G118-030** Deschamp Ligature Needle, sharp, for left hand, 21cm




- G118-025** Deschamp Ligature Needle, blunt, for right hand, 21cm

- G118-035** Deschamp Ligature Needle, sharp, for right hand, 21cm



General

Ligature Needles



- G118-205** Jolls Ligature Needle, 20cm, curve no.1
G118-210 Jolls Ligature Needle, 20cm, curve no.2
G118-215 Jolls Ligature Needle, 20cm, curve no.3
G118-200 Jolls Ligature Needles, set of 3 sizes, 20cm



- G118-305** Syme Ligature Needle, lateral curve, 17cm

Percussors



- G162-130** National Hospital Percussor (Rossier), rubber head & nylon handle, 35cm



- G162-135** Taylor Percussor, 19.5cm

Plaster Instruments




- G094-105** Bohler (Wolf) Plaster Cast Bending Forcep, 18cm
G094-110 Bohler (Wolf) Plaster Cast Bending Forcep, 24cm



- G094-205** Bohler Plaster Scissors / Shear, saw edge, 24cm

Plaster Instruments

G094-505 Daws Plaster Cast Spreader, 45cm



G094-605 Guys (Bruns) Plaster Cast Scissors / Shear, smooth edge, 24cm



G094-705 Hennig Plaster Cast Spreader, 28cm



G094-905 Stille-Lorenz Plaster Cast Shear, 23cm

G094-910 Stille-Lorenz Plaster Cast Shear, 26cm

G094-915 Stille-Lorenz Plaster Cast Shear, 30cm

G094-920 Stille-Lorenz Plaster Cast Shear, 37cm



Ring Cutter

G018-105 Ring Cutter, 16.5cm

G018-110 Blade For Ring Cutter



General

Probes




- G130-112** Nickel Silver Probe, s/e, with eye, 12cm
- G130-115** Nickel Silver Probe, s/e, with eye, 15cm
- G130-118** Nickel Silver Probe, s/e, with eye, 18cm
- G130-120** Nickel Silver Probe, s/e, with eye, 20cm
- G130-123** Nickel Silver Probe, s/e, with eye, 23cm
- G130-125** Nickel Silver Probe, s/e, with eye, 25cm
- G130-130** Nickel Silver Probe, s/e, with eye, 30cm
- G130-212** Silver Probe With Eye, s/e, 12.5cm
- G130-215** Silver Probe With Eye, s/e, 15cm
- G130-218** Silver Probe With Eye, s/e, 18cm
- G130-220** Silver Probe With Eye, s/e, 20cm
- G130-223** Silver Probe With Eye, s/e, 23cm
- G130-225** Silver Probe With Eye, s/e, 25cm
- G130-230** Silver Probe With Eye, s/e, 30cm




- G130-500** Sinus Probe, silver, d/e, 15cm


Post Mortem




- G400-200** Post Mortem - Chisel, cross handle, detachable, 25mm wide



- G400-300** Post Mortem - Spine Chisel and Wrench, with locating point




- G400-410** Post Mortem - Eslander Rib Shear, 21.5cm



- G400-413** Post Mortem - Pollocks Rib Shear, 31.1cm

G134-130 Canny Ryall Retractor, double ended, large blunt hook, small blunt double hook, 16.5cm

G134-135 Canny Ryall Retractor, double ended, small blunt hook, large blunt double hook, 16.5cm



G134-160 Canny Ryall Retractor, double ended, 16mm wide fenestrated blade, 5mm wide blunt hook, 21cm

G134-165 Canny Ryall Retractor, double ended, 20mm wide fenestrated blade, 5mm wide blunt hook, 21cm

G134-170 Canny Ryall Retractor, double ended, 25mm wide fenestrated blade, 5mm wide blunt hook, 21cm

G134-175 Canny Ryall Retractor, double ended, 30mm wide fenestrated blade, 5mm wide blunt hook, 21cm



G134-235 Czerny Retractor, d/e, 4.0cm x 2.2cm blades, adult 18Cm



G134-242 Deaver Retractor, flat handle, 25mm x 23cm blade

G134-244 Deaver Retractor, flat handle, 25mm x 25.5cm blade

G134-246 Deaver Retractor, flat handle, 25mm x 30cm blade

G134-239 Deaver Retractor, flat handle, 25mm x 33cm blade

G134-247 Deaver Retractor, flat handle, 40mm x 30cm blade

G134-248 Deaver Retractor, flat handle, 50mm x 30cm blade

G134-249 Deaver Retractor, flat handle, 75mm x 30cm blade



General

Retractors



G134-250 Deaver Retractor, hollow handle, 25mm x 23cm blade

G134-252 Deaver Retractor, hollow handle, 25mm x 25.5cm blade

G134-254 Deaver Retractor, hollow handle, 25mm x 30cm blade

G134-255 Deaver Retractor, hollow handle, 40mm x 30cm blade

G134-256 Deaver Retractor, hollow handle, 50mm x 30cm blade

G134-257 Deaver Retractor, hollow handle, 75mm x 30cm blade



G134-275 Doyen Retractor, retrograde curve, hollow handle, 9cm wide, small

G134-277 Doyen Retractor, retrograde curve, hollow handle, 10cm wide, medium


G134-279 Doyen Retractor, retrograde curve, hollow handle, 11.5cm wide, large



G134-274 Doyen Retractor, retrograde curve, solid handle, 9cm wide, small

G134-276 Doyen Retractor, retrograde curve, solid handle, 10cm wide, medium

G134-278 Doyen Retractor, retrograde curve, solid handle, 11.5cm wide, large



G134-281 Durham Retractor, fenestrated blade, 10mm blade, 21.5cm

G134-282 Durham Retractor, fenestrated blade, 16mm blade, 21.5cm

G134-283 Durham Retractor, fenestrated blade, 22mm blade, 21.5cm




G134-180 Cushing Wound & Vein Retractor, blade 14mm x 11mm, 25cm

G134-182 Cushing Wound & Vein Retractor, blade 18mm x 13mm, 25cm

G134-184 Cushing Wound & Vein Retractor, blade 22mm x 15mm, 25cm

G134-290 Dyball And Lock liver Retractor, 24cm, 6.3cm x 12cm blade



G134-292 Doyen (Fritsch) Retractor, 45mm x 40mm blade

G134-293 Doyen (Fritsch) Retractor, 45mm x 50mm blade

G134-294 Doyen (Fritsch) Retractor, 45mm x 60mm blade


G134-295 Doyen (Fritsch) Retractor, 45mm x 75mm blade

G134-296 Doyen (Fritsch) Retractor, 65mm x 85mm blade



G134-550 Jolls Thyroid Retractor, 30mm x 14mm, 18cm

G134-552 Jolls Thyroid Retractor, 40mm x 20mm, 18cm



G134-311 Hook Retractor, 1 prong, blunt, 16cm

G134-312 Hook Retractor, 2 prong, blunt, 16cm

G134-313 Hook Retractor, 3 prong, blunt, 16cm

G134-314 Hook Retractor, 4 prong, blunt, 16cm

G134-321 Hook Retractor, 1 prong, sharp, 16cm

G134-322 Hook Retractor, 2 prong, sharp, 16cm


G134-323 Hook Retractor, 3 prong, sharp, 16cm

G134-324 Hook Retractor, 4 prong, sharp, 16cm



General

Retractors




G134-382 Kelly Retractor, ring handle, 160mm x 38mm blade , 27cm

G134-384 Kelly Retractor, ring handle, 160mm x 57mm blade, 27cm


G134-386 Kelly Retractor, ring handle, 190mm x 38mm blade, 27cm

G134-388 Kelly Retractor, ring handle, 190mm x 57mm blade, 27cm




G134-375 Kelly Retractor, 190mm x 40mm wide, 27cm


G134-380 Kelly Retractor, 190mm x 60mm wide, 27cm



G134-511 Kocher Thyroid Retractor, 11mm wide shallow blade, 16cm




G134-514 Kocher Thyroid Retractor, 14mm wide deep blade, 23cm



G134-510 Lane Retractor, d/e, 19mm tapering, 13mm Blade, 23cm

G134-515 Lane Retractor, d/e, 32mm tapering, 25mm Blade, 23cm

- G134-430** Langenbeck Retractor, 6mm wide x 22mm, 22cm
- G134-432** Langenbeck Retractor, 13mm wide x 45mm, 22cm
- G134-434** Langenbeck Retractor, 25mm Wide x 65mm, 22cm
- G134-436** Langenbeck Retractor, 20mm wide x 90mm, 22cm



- G134-412** Langenbeck Retractor, ring handle, 30mm x 11mm blade, 21cm
- G134-413** Langenbeck Retractor, ring handle, 30mm x 14mm blade, 21cm
- G134-414** Langenbeck Retractor, ring handle, 33mm x 14mm blade, 21cm
- G134-415** Langenbeck Retractor, ring handle, 30mm x 16mm blade, 21cm
- G134-416** Langenbeck Retractor, ring handle, 40mm x 10mm blade, 21cm
- G134-417** Langenbeck Retractor, ring handle, 42mm x 13mm blade, 21cm
- G134-418** Langenbeck Retractor, ring handle, 50mm x 11mm blade, 21cm
- G134-419** Langenbeck Retractor, ring handle, 60mm x 20mm blade, 21cm



- G134-460** Morris Retractor, single ended, 40mm x 70mm, 25cm
- G134-464** Morris Retractor, single ended, 50mm x 70mm, 25cm
- G134-466** Morris Retractor, single ended, 65mm x 70mm, 25cm




- G134-400** Kocher-Langenbeck Retractor, 25mm x 6mm blade, 21cm
- G134-401** Kocher-Langenbeck Retractor, 25mm x 8mm blade, 21cm
- G134-402** Kocher-Langenbeck Retractor, 35mm x 8mm blade, 21cm
- G134-403** Kocher-Langenbeck Retractor, 35mm x 11mm blade, 21cm
- G134-404** Kocher-Langenbeck Retractor, 35mm x 15mm blade, 21cm
- G134-405** Kocher-Langenbeck Retractor, 41mm x 11mm blade, 21cm
- G134-406** Kocher-Langenbeck Retractor, 55mm x 11mm blade, 21cm
- G134-407** Kocher-Langenbeck Retractor, 70mm x 14mm blade, 21cm
- G134-408** Kocher-Langenbeck Retractor, 80mm x 12mm blade, 21cm
- G134-409** Kocher-Langenbeck Retractor, 80mm x 16mm blade, 21cm



General

Retractors



G134-480 Moynihan Retractor, single ended, 21cm, 13mm x 22mm blade



G134-420 Mini Langenbeck Retractor, 6mm x 10mm blade, 16cm

G134-422 Mini Langenbeck Retractor, 6mm x 17mm blade, 16cm

G134-424 Mini Langenbeck Retractor, 6mm x 20mm blade, 16cm

G134-426 Mini Langenbeck Retractor, 8mm x 22mm blade, 16cm




G134-502 Kocher Retractor, 40mm x 18mm blade, 23cm

G134-503 Kocher Retractor, 60mm x 20mm blade, 23cm

G134-504 Kocher Retractor, 60mm x 25mm blade, 23cm

G134-505 Kocher Retractor, 75mm x 30mm blade, 23cm

G134-506 Kocher Retractor, 75mm x 40mm blade, 23cm



G134-600 Ollier Retractor, 2 prong blunt, 25mm x 45mm, 23cm

G134-602 Ollier Retractor, 3 prong blunt, 36mm x 30mm, 23cm

G134-603 Ollier Retractor, 4 prong blunt, 36mm x 60mm, 23cm




G134-820 Parker-Langenbeck Army-Navy Retractor, set /2, 21cm



G134-485 Paton Higgs Retractor, ring handle, 30mm x 10mm blade, 21.5cm

G134-487 Paton Higgs Retractor, ring handle, 42mm x 13mm blade, 21.5cm

G134-490 Proust Retractor, hollow handle, 65mm x 45mm blade




G134-530 Pyrah Retractor, hollow handle, 75mm x 50mm, 23cm

G134-532 Pyrah Retractor, hollow handle, 100mm x 50mm, 23cm

G134-534 Pyrah Retractor, hollow handle, 100mm x 65mm, 23cm

G134-536 Pyrah Retractor, hollow handle, 125mm x 70mm, 23cm

G134-538 Pyrah Retractor, hollow handle, 150mm x 80mm, 23cm



G134-842 Baby Roux Retractor, 13cm

G134-844 Roux Retractor, small, 14cm

G134-846 Roux Retractor, medium, 15cm

G134-848 Roux Retractor, large, 17cm




G134-921 Middeldorf Retractor, ring handle, fenestrated blade, 15mm x 15mm, 22cm

G134-922 Middeldorf Retractor, ring handle, fenestrated blade, 20mm x 22mm, 22cm


G134-923 Middeldorf Retractor, ring handle, fenestrated blade, 28mm x 28mm, 23cm





G134-587 St. Marks Retractor, hollow handle, 125mm x 55mm blade, with lip, complete with twin fibre light

G134-588 St. Marks Retractor, hollow handle, 175mm x 60mm blade, with lip, complete with twin fibre light




G134-906 Richardson Retractor, loop handle, 28mm x 20mm blade, 26cm

G134-908 Richardson Retractor, loop handle, 30mm x 28mm blade, 26cm

G134-910 Richardson Retractor, loop handle, 38mm x 37mm blade, 26cm

G134-912 Richardson Retractor, loop handle, 52mm x 22mm blade, 26cm



G134-585 St. Marks Retractor, hollow handle, 125mm x 55mm blade, without lip, complete with twin fibre light

G134-586 St. Marks Retractor, hollow handle, 175mm x 65mm blade, without lip, complete with twin fibre light



G134-920 Richardson-Eastman Retractors, set /2, double ended


- E162-115** Tracheal Hook Retractor, double, blunt, 16cm
E162-120 Tracheal Hook Retractor, double, sharp, 16cm




- E162-105** Tracheal Hook Retractor, single, blunt, 16cm
E162-110 Tracheal Hook Retractor, single, sharp, 16cm



- G134-644** Volkmann Retractor, blunt, 22cm, 1 prong
G134-646 Volkmann Retractor, blunt, 22cm, 2 prong
G134-648 Volkmann Retractor, blunt, 22cm, 3 prong
G134-650 Volkmann Retractor, blunt, 22cm, 4 prong
G134-654 Volkmann Retractor, blunt, 22cm, 6 prong
G134-655 Volkmann Retractor, blunt, 22cm, 8 prong
G134-669 Volkmann Retractor, semi-sharp, 22cm, 2 prong
G134-670 Volkmann Retractor, semi-sharp, 22cm, 4 prong
G134-656 Volkmann Retractor, sharp, 22cm, 1 prong
G134-658 Volkmann Retractor, sharp, 22cm, 2 prong
G134-660 Volkmann Retractor, sharp, 22cm, 3 prong
G134-662 Volkmann Retractor, sharp, 22cm, 4 prong
G134-666 Volkmann Retractor, sharp, 22cm, 6 prong
G134-667 Volkmann Retractor, sharp, 22cm, 8 prong




- G134-671** Kocher Retractor, sharp, 2 prong, 16mm x 9mm, 22cm
G134-672 Kocher Retractor, sharp, 3 prong, 16mm x 14mm, 22cm
G134-673 Kocher Retractor, sharp, 4 prong, 16mm x 20mm, 22cm
G134-674 Kocher Retractor, sharp, 6 prong, 16mm x 30mm, 22cm
G134-676 Kocher Retractor, blunt, 2 prong, 16mm x 9mm, 22cm
G134-678 Kocher Retractor, blunt, 4 prong, 16mm x 20mm, 22cm
G134-679 Kocher Retractor, blunt, 6 prong, 16mm x 30mm, 22cm




General

Retractors



- G134-760** Walton Retractor, copper, malleable, 18cm x 13mm
- G134-761** Walton Retractor, copper, malleable, 18cm x 20mm
- G134-762** Walton Retractor, copper, malleable, 18cm x 25mm
- G134-763** Walton Retractor, copper, malleable, 18cm x 30mm
- G134-764** Walton Retractor, copper, malleable, 18cm x 40mm
- G134-765** Walton Retractor, copper, malleable, 18cm x 50mm
- G134-766** Walton Retractor, copper, malleable, 18cm x 65mm
- G134-767** Walton Retractor, copper, malleable, 25.5cm x 13mm
- G134-768** Walton Retractor, copper, malleable, 25.5cm x 20mm
- G134-769** Walton Retractor, copper, malleable, 25.5cm x 25mm
- G134-770** Walton Retractor, copper, malleable, 25.5cm x 30mm
- G134-771** Walton Retractor, copper, malleable, 25.5cm x 40mm
- G134-772** Walton Retractor, copper, malleable, 25.5cm x 50mm
- G134-773** Walton Retractor, copper, malleable, 25.5cm x 65mm
- G134-774** Walton Retractor, copper, malleable, 30cm x 13mm
- G134-775** Walton Retractor, copper, malleable, 30cm x 20mm
- G134-776** Walton Retractor, copper, malleable, 30cm x 25mm
- G134-777** Walton Retractor, copper, malleable, 30cm x 30mm
- G134-778** Walton Retractor, copper, malleable, 30cm x 40mm
- G134-779** Walton Retractor, copper, malleable, 30cm x 50mm
- G134-780** Walton Retractor, copper, malleable, 30cm x 65mm



- G134-751** Ribbon Retractor, malleable, 200mm x 6mm
- G134-752** Ribbon Retractor, malleable, 200mm x 12mm
- G134-754** Ribbon Retractor, malleable, 200mm x 17mm
- G134-755** Ribbon Retractor, malleable, 200mm x 20mm
- G134-756** Ribbon Retractor, malleable, 250mm x 25mm
- G134-758** Ribbon Retractor, malleable, 325mm x 25mm
- G134-797** Ribbon Retractor, malleable, 330mm x 30mm
- G134-798** Ribbon Retractor, malleable, 330mm x 40mm
- G134-799** Ribbon Retractor, malleable, 330mm x 50mm

Retractors

- G134-783** Haberer (Tuffier) Abdominal and Intestinal Retractor, malleable, 17mm x 25mm wide, 20cm
- G134-785** Haberer (Tuffier) Abdominal and Intestinal Retractor, malleable, 25mm x 30mm wide, 28cm
- G134-787** Haberer (Tuffier) Abdominal and Intestinal Retractor, malleable, 37mm x 45mm wide, 30cm
- G134-789** Haberer (Tuffier) Abdominal and Intestinal Retractor, malleable, 40mm x 50mm wide, 30cm




Retractors - Self Retaining


- G136-118** Balfour Self Retaining Retractor, baby, 30mm deep outer blades, 27mm x 27mm centre blade 95mm spread
- G136-120** Balfour Self Retaining Retractor, standard, 63mm x 35mm deep outer blades, 45mm x 80mm centre blade 180mm spread
- G136-122** Balfour Self Retaining Retractor, large, 100mm x 35mm deep outer blades, 70mm x 85mm centre blade 250mm spread
- G136-124** Balfour Self Retaining Retractor, deep, 100mm x 35mm deep outer blades, 100mm x 70mm centre blade 200mm spread



- G136-100** Blake Amputation Retractor, self retaining, adjusting circular shield




- G136-102** Percy Amputation Retractor, with detachable handles



Retractors - Self Retaining



G136-235 Comyns Berkeley Abdominal Self Retaining Retractor, 2 fixed and adjustable blades



G136-151 Clelland Abdominal Self Retaining Retractor, deep angled blade only, 115mm x 65mm

G136-153 Clelland Abdominal Self Retaining Retractor, frame only, 25cm inside diameter

G136-152 Clelland Abdominal Self Retaining Retractor, shallow curved blade only, 50mm x 40mm

G136-150 Clelland Abdominal Retractor Set, complete set, 1 x frame 25cm, 4 x shallow curved blades, 50mm x 40mm and 4 x deep angled blades, 115mm x 65mm



G136-112 Beckmann-Adson Retractor, self retaining, hinged arms, 25mm blade depth, 4:4 sharp teeth, 31cm

G136-113 Beckmann-Adson Retractor, self retaining, hinged arms, 44mm blade depth, 4:4 sharp teeth, 31cm


G136-114 Beckmann-Adson Retractor, self retaining, hinged arms, 25mm blade depth, 4:4 blunt teeth, 31cm

G136-115 Beckmann-Adson Retractor, self retaining, hinged arms, 44mm blade depth, 4:4 blunt teeth, 31cm

Retractors - Self Retaining

G136-525 Collin (Pozzi) Tri-valve Self Retaining Retractor, complete with detachable 3.8cm x 6.5cm blades and centre blade 5.5cm x 5.5cm

G136-520 Collin (Pozzi) Bi-valve Self Retaining Retractor, complete with detachable 3.8cm x 6.5cm blades



G136-080 Adson Retractor, self retaining, straight arms, 4:4 sharp prongs, 20cm

G136-085 Adson Retractor, self retaining, straight arms, 4:4 blunt prongs, 20cm



G136-103 Adson Retractor, self retaining, 3:4 teeth, sharp, 27cm

G136-104 Adson Retractor, self retaining, 3:4 teeth, blunt, 27cm


G136-105 Adson Retractor, self retaining, 5:6 teeth, sharp, 27cm

G136-107 Adson Retractor, self retaining, 5:6 teeth, blunt, 27cm




G136-161 Denis Browne Oval Sprocket Frame, 18cm x 14cm

G136-160 Denis Browne Abdominal Retractor set, set /1




General


Retractors - Self Retaining




G136-162 Denis Browne Curved Claw Blade, part A, set /4




G136-163 Denis Browne Angled Claw Blade, part B, set /4



G136-164 Denis Browne Curved Blade, part C, set /4




G136-165 Denis Browne Sharp Hook Blade, part D, set /2




G136-166 Denis Browne Lung Retractor Blade (To Fit Part H), part E, set /1

Retractors - Self Retaining


- G136-167** Denis Browne Curved Blade End (To Fit Part H),
part F, set /1




- G136-168** Denis Browne Regressive Blade End (To Fit Part H),
part G, set /1



- G136-169** Denis Browne Retractor Shaft (To Fit Parts E, F & G),
part H, set /1



- G136-170** Denis Browne Spanner (To Fit Part H),
part J, set /1




- G136-731** Weitlaner Retractor, self retaining, straight arms, 3:4
sharp prongs, 14cm
- G136-732** Weitlaner Retractor, self retaining, straight arms, 3:4
sharp prongs, 16.5cm
- G136-733** Weitlaner Retractor, self retaining, straight arms, 3:4
sharp prongs, 20cm
- G136-734** Weitlaner Retractor, self retaining, straight arms, 3:4
sharp prongs, 25cm
- G136-735** Weitlaner Retractor, self retaining, straight arms, 3:4
blunt prongs, 14cm
- G136-736** Weitlaner Retractor, self retaining, straight arms, 3:4
blunt prongs, 16.5cm
- G136-737** Weitlaner Retractor, self retaining, straight arms, 3:4
blunt prongs, 20cm
- G136-738** Weitlaner Retractor, self retaining, straight arms, 3:4
blunt prongs, 25cm




General


Retractors - Self Retaining



- G136-330** Goligher Abdominal Retractor, self retaining frame, complete with 1 x centre blade 51mm deep, 1 x centre blade 102mm deep, 1 x pair lateral blades 51mm deep, 1 x pair lateral blades 102mm deep




- G136-345** Gosset Abdominal Self Retaining Retractor, with 6.5cm x 3.5cm fenestrated blades and centre blade



- G136-342** Gosset Abdominal Retractor, self retaining, large, 6.5cm x 3.5cm fenestrated blades
- G136-340** Gosset Abdominal Retractor, self retaining, small, 4.0cm x 2.8cm fenestrated blades

Retractors - Self Retaining


- G136-600** Irwin Self Retaining Retractor, 4:4 sharp prongs, 17.5cm



- G136-555** Jolls Thyroid Self Retaining Retractor, sharp prong and threaded adjustment bar, 12.5cm



- G136-192** Lace Abdominal Retractor, self retaining, large



- G136-317** Gelpi Retractor, self retaining, single sharp hooks, 13cm

- G136-318** Gelpi Retractor, self retaining, single sharp hooks, blade depth 30mm, 18cm

- G136-320** Gelpi Retractor, self retaining, single sharp hooks, blade depth 50mm, 19.5cm




General


Retractors - Self Retaining




G136-585 Norfolk And Norwich Self Retaining Retractor,
4:5 blunt prongs, 21.5cm



G136-275 Ricard-Begouin Abdominal Retractor,
self retaining, 20cm spread, adjustable
centre blade and two sets of 2 blades,
35mm x 65mm and 65mm x 65mm




G136-293 Swift Joly Retractor Blade, curved, 25mm x 50mm,
part A




G136-294 Swift Joly Retractor Blade, dishd, 25mm x 70mm,
part B

Retractors - Self Retaining

G136-295 Swift Joly Retractor Blade, dished, 25mm x 90mm, part C




G136-296 Swift Joly Retractor Blade, angled, 65mm x 65mm, part D




G136-292 Swift Joly Oval Frame Only, 180mm X 140mm


G136-290 Swift Joly Abdominal Retractor, self retaining, set of 6 blades & oval frame



Retractors - Self Retaining



- G136-860** Stevens Abdominal Self Retaining Retractor, lightweight, 60mm deep x 45mm wide blades, 190g (7oz)




- G136-872** St. Thomas Hospital Pattern Self Retaining Retractor Blade A Only, right angled 145mm deep x 50mm wide, 9 pegs, 23cm

- G136-874** St. Thomas Hospital Pattern Self Retaining Retractor Blade B Only, right angled 125mm deep x 38mm wide, 9 pegs, 23cm

- G136-876** St. Thomas Hospital Pattern Self Retaining Retractor Blade C Only, right angled 120mm deep x 38mm wide, 6 pegs, 18cm

- G136-878** St. Thomas Hospital Pattern Self Retaining Retractor Blade D Only, right angled 95mm deep x 38mm wide, 6 pegs, 18cm



- G136-880** St. Thomas Hospital Pattern Self Retaining Retractor Blade E Only, curved, 45mm deep x 38mm wide, 9 pegs, 21.5cm

- G136-882** St. Thomas Hospital Pattern Self Retaining Retractor Blade F Only, curved, 45mm deep x 38mm wide, 7 pegs, 18cm

- G136-884** St. Thomas Hospital Pattern Self Retaining Retractor Blade G Only, curved, 45mm deep x 38mm wide, 4 pegs, 12.5cm

- G136-886** St. Thomas Hospital Pattern Self Retaining Retractor Blade K Only, curved, 65mm deep x 38mm wide, 9 pegs, 21.5cm

- G136-888** St. Thomas Hospital Pattern Self Retaining Retractor Blade L Only, curved, 65mm deep x 38mm wide, 7 pegs, 18cm


- G136-890** St. Thomas Hospital Pattern Self Retaining Retractor Blade M Only, curved, 65mm deep x 38mm wide, 4 pegs, 12.5cm

- G136-892** St. Thomas Hospital Pattern Self Retaining Retractor Blade N Only, curved, 85mm deep x 38mm wide, 9 pegs, 21.5cm

Retractors - Self Retaining

G136-894 St. Thomas Hospital Pattern Self Retaining Retractor Blade H Only, deaver type, 115mm x 70mm tapering to 64mm, 6 pegs, 23cm


G136-896 St. Thomas Hospital Pattern Self Retaining Retractor Blade J Only, deaver type, 150mm x 70mm tapering to 64mm, 6 pegs, 25.5cm



G136-864 St. Thomas Hospital Pattern Self Retaining Retractor Frame A Only, 255mm x 345mm

G136-866 St. Thomas Hospital Pattern Self Retaining Retractor Frame B Only, 330mm x 445mm

G136-870 St. Thomas Hospital Pattern Self Retaining Retractor Complete With Frame A, and blades A, B, C, D, E, F, G, H And J



G136-580 Travers Self Retaining Retractor, 4:5 blunt prongs, straight, 20cm



G136-808 Turner Warwick Retractor, blade A, 19mm x 76mm



General

Retractors - Self Retaining



G136-810 Turner Warwick Retractor, blade B,
38mm x 102mm




G136-812 Turner Warwick Retractor, blade C,
57mm x 89mm




G136-814 Turner Warwick Retractor, blade D,
57mm x 114mm

Retractors - Self Retaining

G136-816 Turner Warwick Retractor, blade E,
51mm x 38mm



G136-818 Turner Warwick Retractor, blade F,
76mm x 60mm



G136-820 Turner Warwick Retractor, blade G,
89mm x 60mm



G136-822 Turner Warwick Retractor, blade H,
114mm x 60mm




General


Retractors - Self Retaining




G136-824 Turner Warwick Retractor, blade J,
13mm



G136-825 Turner Warwick Retractor, blade K,
55mm x 90mm




G136-826 Turner Warwick Retractor, blade L,
55mm x 120mm




G136-827 Turner Warwick Retractor, blade M,
31mm x 55mm

Retractors - Self Retaining

- G136-828** Turner Warwick Retractor, blade N,
60mm x 55mm



- G136-829** Turner Warwick Retractor, blade P,
28mm x 30mm



- G136-800** Turner Warwick Circular Retractor, Frame Only,
18cm Dia


- G136-802** Turner Warwick Circular Retractor, Frame Only,
23cm Dia

- G136-804** Turner Warwick Circular Retractor, Frame Only,
30cm Dia

- G136-805** Turner Warwick Oval Retractor, 457mm x 330mm
I.D. Frame Only

- G136-806** Turner Warwick Oval Retractor, 370mm x 300mm
I.D. Frame Only

- G136-799** Turner Warwick Abdominal Retractor Set,
consisting of 1 x ea round frame 18cm, 23cm
and 30cm, 2 x ea blades A, B, C, D, E, F, G, H
and J




- G136-710** Weislander Self Retaining Retractor, 3:4 blunt prongs,
14cm




General

Retractors - Self Retaining



G136-720 West Mastoid Self Retaining Retractor, 3:4 sharp prongs, 14cm

Saws



G138-010 Satterlee Amputation Saw, 31cm



G094-305 Engel Plaster Saw

Scalpel Handles



G102-043 Swann Morton Surgical Handle, no.3, graduated, stainless steel, for no.3 fitment blade



G102-048 Swann Morton Surgical Handle, no.3L, stainless steel, for no.3 fitment blade



G102-044 Swann Morton Surgical Handle, no.4, graduated, stainless steel, for no.4 fitment blade



G102-045 Swann Morton Surgical Handle no.4L, stainless steel, for no.4 fitment blade

Scalpel Handles


G102-047 Swann Morton Surgical Handle no.7, stainless steel, for no.3 fitment blade



Scissors


G140-022 Aufricht Dissecting Scissor, curved, fine sharp points, 14cm

G140-020 Aufricht Dissecting Scissor, straight, fine sharp points, 14cm




G140-026 Aufricht Dissecting Scissor, curved, fine triangular points, 14cm

G140-024 Aufricht Dissecting Scissor, straight, fine triangular points, 14cm




G140-076 Carless Scissor, straight, lightweight, blunt / blunt, 16.5cm

G140-078 Carless Scissor, straight, lightweight, blunt / blunt, 20cm



G140-150 Doyen Abdominal Scissor, heavy, straight, blunt / blunt, 18cm



G140-201 Dressing / Operating Scissor, straight, blunt / blunt, 10.5cm

G140-202 Dressing / Operating Scissor, straight, blunt / blunt, 12.5cm

G140-203 Dressing / Operating Scissor, straight, blunt / blunt, 14cm

G140-204 Dressing / Operating Scissor, straight, blunt / blunt, 15cm

G140-205 Dressing / Operating Scissor, straight, blunt / blunt, 18cm

G140-206 Dressing / Operating Scissor, straight, blunt / blunt, 20cm



General

Scissors



- G140-221** Dressing / Operating Scissor, sharp / blunt, straight, 10.5cm
- G140-222** Dressing / Operating Scissor, sharp / blunt, straight, 12.5cm
- G140-223** Dressing / Operating Scissor, sharp / blunt, straight, 14cm
- G140-224** Dressing / Operating Scissor, sharp / blunt, straight, 15cm
- G140-225** Dressing / Operating Scissor, sharp / blunt, straight, 18cm
- G140-226** Dressing / Operating Scissor, sharp / blunt, straight, 20cm



- G140-241** Dressing / Operating Scissor, sharp / sharp, straight, 10.5cm
- G140-242** Dressing / Operating Scissor, sharp / sharp, straight, 12.5cm
- G140-243** Dressing / Operating Scissor, sharp / sharp, straight, 14cm
- G140-244** Dressing / Operating Scissor, sharp / sharp, straight, 15cm
- G140-245** Dressing / Operating Scissor, sharp / sharp, straight, 18cm
- G140-246** Dressing / Operating Scissor, sharp / sharp, straight, 20cm



- G140-211** Dressing / Operating Scissor, blunt / blunt, curved, 10.5cm
- G140-212** Dressing / Operating Scissor, blunt / blunt, curved, 12.5cm
- G140-213** Dressing / Operating Scissor, blunt / blunt, curved, 14cm
- G140-214** Dressing / Operating Scissor, blunt / blunt, curved, 15cm
- G140-215** Dressing / Operating Scissor, blunt / blunt, curved, 18cm
- G140-216** Dressing / Operating Scissor, blunt / blunt, curved, 20cm




- G140-231** Dressing / Operating Scissor, sharp / blunt, curved, 10.5cm
- G140-232** Dressing / Operating Scissor, sharp / blunt, curved, 12.5cm
- G140-233** Dressing / Operating Scissor, sharp / blunt, curved, 14cm
- G140-234** Dressing / Operating Scissor, sharp / blunt, curved, 15cm
- G140-235** Dressing / Operating Scissor, sharp / blunt, curved, 18cm
- G140-236** Dressing / Operating Scissor, sharp / blunt, curved, 20cm



- G140-251** Dressing / Operating Scissor, sharp / sharp, curved, 10.5cm
- G140-252** Dressing / Operating Scissor, sharp / sharp, curved, 12.5cm
- G140-253** Dressing / Operating Scissor, sharp / sharp, curved, 14cm
- G140-254** Dressing / Operating Scissor, sharp / sharp, curved, 15cm
- G140-255** Dressing / Operating Scissor, sharp / sharp, curved, 18cm
- G140-256** Dressing / Operating Scissor, sharp / sharp, curved, 20cm

G140-041 Blood Vessel Scissor, probe ended blade, straight, 11.5cm



G140-250 Iris Scissor, straight, sharp / sharp, 9cm


G140-260 Iris Scissor, curved, sharp / sharp, 9cm

G140-262 Iris Scissor, straight, sharp / sharp, 10.5cm

G140-264 Iris Scissor, curved, sharp / sharp, 10.5cm

G140-270 Iris Scissor, straight, sharp / sharp, 11.5cm

G140-275 Iris Scissor, curved, sharp / sharp, 11.5cm




G140-282 Fergusson Dressing / Operating Scissor, blunt / blunt, angled up, 14cm

G140-284 Fergusson Dressing / Operating Scissor, blunt / blunt, angled up, 15.5cm

G140-286 Fergusson dressing / operating Scissor, blunt / blunt, angled up, 18cm

G140-287 Fergusson dressing / operating Scissor, blunt / blunt, angled up, 20cm




G140-280 Joseph Scissor, straight, sharp / sharp, 14cm

G140-281 Joseph Scissor, curved, sharp / sharp, 14cm



G140-315 Heath Stitch Scissor, sharp / sharp, serrated blades, angled laterally, 15cm



General

Scissors




G140-031 Lister Bandage Scissor, probe end, 14cm

G140-033 Lister Bandage Scissor, probe end, 18cm

G140-035 Lister (Esmarch) Bandage Scissor, probe end, 20cm

G140-037 Lister (Bergmann) Bandage / Plaster Scissor, probe end, 23cm



G140-027 Lister Excentric Bandage Scissor, probe end, saw edge, 16cm

G140-029 Lister Bandage Scissor, probe end, 9cm



G140-058 Bowel Scissor, hook ended lower blade, straight, 21cm

G140-060 Bowel Scissor, probe ended lower blade, straight, 21cm




G140-774 Littauer Stitch Scissor, hook end, 14cm

G140-470 Mayo Scissor, conical blades, round points, curved, 14cm

G140-478 Mayo Scissor, conical blades, round points, curved, 17cm

G140-482 Mayo Scissor, conical blades, round points, curved, 19cm

G140-486 Mayo Scissor, conical blades, round points, curved, 23cm




G140-495 Mayo-Stille Scissor, curved, 14cm

G140-498 Mayo-Stille Scissor, curved, 17cm

G140-502 Mayo-Stille Scissor, curved, 19cm

G140-504 Mayo-Stille Scissor, curved, 21cm




G140-468 Mayo Scissor, conical blades, round points, straight, 14cm

G140-476 Mayo Scissor, conical blades, round points, straight, 17cm

G140-480 Mayo Scissor, conical blades, round points, straight, 19cm

G140-484 Mayo Scissor, conical blades, round points, straight, 23cm




G140-494 Mayo-Stille Scissor, straight, 14cm

G140-496 Mayo-Stille Scissor, straight, 17cm

G140-500 Mayo-Stille Scissor, straight, 19cm

G140-503 Mayo-Stille Scissor, straight, 21cm



G140-562 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 14.5cm

G140-564 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 18cm


G140-566 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 20cm

G140-570 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 23cm

G140-574 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 25cm


G140-576 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 28cm

G140-580 Metzenbaum Dissecting Scissor, standard blades, blunt / blunt, curved, 30cm




G140-561 Metzenbaum Dissecting Scissor, standard blades, curved, sharp / sharp, 14.5cm

G140-602 Metzenbaum Dissecting Scissor, standard blades, curved, sharp / sharp, 18cm



Scissors



G140-642 Metzenbaum Dissecting Scissor, fine blades, straight, blunt / blunt, 14.5cm

G140-644 Metzenbaum Dissecting Scissor, fine blades, straight, blunt / blunt, 18cm

G140-646 Metzenbaum Dissecting Scissor, fine blades, straight, blunt / blunt, 20cm


G140-650 Metzenbaum Dissecting Scissor, fine blades, straight, blunt / blunt, 23cm

G140-662 Metzenbaum Dissecting Scissor, fine blades, curved, blunt / blunt, 14.5cm

G140-664 Metzenbaum Dissecting Scissor, fine blades, curved, blunt / blunt, 18cm


G140-666 Metzenbaum Dissecting Scissor, fine blades, curved, blunt / blunt, 20cm

G140-670 Metzenbaum Dissecting Scissor, fine blades, curved, blunt / blunt, 23cm



G140-541 Metzenbaum Dissecting Scissor, standard blades, straight, sharp / sharp, 14.5cm

G140-590 Metzenbaum Dissecting Scissor, standard blades, straight, sharp / sharp, 18cm



G140-542 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 14.5cm

G140-544 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 18cm

G140-546 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 20cm

G140-550 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 23cm


G140-554 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 25cm

G140-556 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 28cm

G140-560 Metzenbaum Dissecting Scissor, standard blades, Blunt / Blunt, Straight, 30cm




G140-365 Nelson Lobectomy Scissor, curved, blunt / blunt, 23cm




G140-363 Nelson Lobectomy Scissor, straight, blunt / blunt, 23cm

G400-047 Post Mortem - Scissor, ball end, for ducts, 14cm




G400-045 Post Mortem - Scissor, open shank, sharp / probe end, for coronary arteries, 10cm




G140-509 Spencer Stitch Scissor, hook end, 9cm
G140-511 Spencer Stitch Scissor, hook end, 11cm



G140-826 Stevens Tenotomy Scissor, curved, blunt / blunt, 11.5cm




G140-836 Stevens Tenotomy Scissor, curved, sharp / sharp, 11.5cm




General

Scissors




G140-824 Stevens Tenotomy Scissor, straight, blunt / blunt, 11.5cm



G140-834 Stevens Tenotomy Scissor, straight, sharp / sharp, 11.5cm



G140-920 Universal Bandage Scissor, saw edge, black plastic handles, autoclavable, angled, 18cm



E172-805 Universal Wire Cutting Scissor, angled serrated blade, to cut max 0.6mm hard wire and 1.2mm soft wire, 12cm



G140-800 Strabismus Scissor, straight, conical blades, round points, 9cm

G140-802 Strabismus Scissor, straight, conical blades, round points, 11cm

G140-804 Strabismus Scissor, curved, conical blades, round points, 9cm


G140-805 Strabismus Scissor, curved, conical blades, round points, 11cm

Scissors - Tungsten Carbide

G141-372 Operating Scissor, tungsten carbide blades, straight, blunt / blunt, 14.5cm

G141-374 Operating Scissor, tungsten carbide blades, straight, blunt / sharp, 14.5cm


G141-376 Operating Scissor, tungsten carbide blades, straight, sharp / sharp, 14.5cm



G141-378 Operating Scissor, tungsten carbide blades, curved, blunt / blunt, 14.5cm


G141-380 Operating Scissor, tungsten carbide blades, curved, blunt / sharp, 14.5cm

G141-382 Operating Scissor, tungsten carbide blades, curved, sharp / sharp, 14.5cm



G141-270 Iris Scissor, tungsten carbide blades, straight, sharp / sharp, 11.5cm

G141-275 Iris Scissor, tungsten carbide blades, curved, sharp / sharp, 11.5cm



G141-468 Mayo Scissor, tungsten carbide blades, conical blades, round points, straight, 14cm


G141-476 Mayo Scissor, tungsten carbide blades, conical blades, round points, straight, 17cm

G141-480 Mayo Scissor, tungsten carbide blades, conical blades, round points, straight, 19cm

G141-470 Mayo Scissor, tungsten carbide blades, conical blades, round points, curved, 14cm


G141-478 Mayo Scissor, tungsten carbide blades, conical blades, round points, curved, 17cm

G141-482 Mayo Scissor, tungsten carbide blades, conical blades, round points, curved, 19cm




General


Scissors - Tungsten Carbide




- G141-494** Mayo-Stille Scissor, tungsten carbide blades, straight, 14cm
- G141-496** Mayo-Stille Scissor, tungsten carbide blades, straight, 17cm
- G141-495** Mayo-Stille Scissor, tungsten carbide blades, curved, 14cm
- G141-498** Mayo-Stille Scissor, tungsten carbide blades, curved, 17cm




- G141-542** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, straight, 14.5cm
- G141-544** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, straight, 18cm
- G141-546** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, straight, 20cm
- G141-550** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, straight, 23cm



- G141-562** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, curved, 14.5cm
- G141-564** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, curved, 18cm
- G141-566** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, curved, 20cm
- G141-570** Metzenbaum Dissecting Scissor, tungsten carbide blade, standard blade, blunt / blunt, curved, 23cm



- G141-365** Nelson Lobectomy Scissor, tungsten carbide, blades, blunt / blunt, curved, 23cm




- G141-363** Nelson Lobectomy Scissor, tungsten carbide, blades, blunt / blunt, straight, 23cm



- G141-802** Strabismus Scissor, tungsten carbide blades, conical blades, round points, straight, 11cm
- G141-805** Strabismus Scissor, tungsten carbide blades, conical blades, round points, curved, 11cm

Scissors - Tungsten Carbide


- E173-805** Universal Wire Cutting Scissor, tungsten carbide edge, angled serrated blade, to cut max 0.8mm hard wire and 1.2mm soft wire, 12cm



Scissors - Super Cut


- G142-374** Micro / Super Cut, operating scissor, precision cutting blade / micro serrated blade, sharp / blunt, straight, 14.5cm

- G142-380** Micro / Super Cut, operating scissor, precision cutting blade / micro serrated blade, sharp / blunt, curved, 14.5cm



- G142-270** Micro / Super Cut, iris Scissor, precision cutting blade / micro serrated blade, flat blade, straight, 11.5cm


- G142-275** Micro / Super Cut, iris Scissor, precision cutting blade / micro serrated blade, flat blade, curved, 11.5cm



- G142-780** Micro / Super Cut, sims uterine scissor, precision cutting blade / micro serrated blade, straight, 20cm


- G142-783** Micro / Super Cut, sims uterine scissor, precision cutting blade / micro serrated blade, curved, 20cm

- G142-786** Micro / Super Cut, sims uterine scissor, precision cutting blade / micro serrated blade, curved, 23cm




General


Scissors - Super Cut




- G142-494** Mayo Scissor, micro / super cut, straight, 14cm
- G142-496** Mayo Scissor, micro / super cut, straight, 17cm
- G142-503** Mayo Scissor, micro / super cut, straight, 23cm
- G142-495** Mayo Scissor, micro / super cut, curved, 14cm
- G142-498** Mayo Scissor, micro / super cut, curved, 17cm
- G142-504** Mayo Scissor, micro / super cut, curved, 23cm




- G142-468** Mayo-Stille Scissor, micro / super cut, straight, 15cm
- G142-476** Mayo-Stille Scissor, micro / super cut, straight, 17cm
- G142-470** Mayo-Stille Scissor, micro / super cut, curved, 15cm
- G142-478** Mayo-Stille Scissor, micro / super cut, curved, 17cm



- G142-542** Metzenbaum Scissor, micro / super cut, straight, 14.5cm
- G142-544** Metzenbaum Scissor, micro / super cut, straight, 18cm
- G142-562** Metzenbaum Scissor, micro / super cut, curved, 14.5cm
- G142-564** Metzenbaum Scissor, micro / super cut, curved, 18cm
- G142-566** Metzenbaum Scissor, micro / super cut, curved, 20cm
- G142-570** Metzenbaum Scissor, micro / super cut, curved, 23cm
- G142-574** Metzenbaum Scissor, micro / super cut, curved, 25cm



- G142-824** Micro / Super Cut, Stevens tenotomy scissor, precision cutting blade / micro serrated blade, straight, 11.5cm
- G142-826** Micro / Super Cut, Stevens tenotomy scissor, precision cutting blade / micro serrated blade, curved, 11.5cm




- G142-725** Potts Smith Scissor, micro / super cut, 25o, 19cm
- G142-745** Potts Smith Scissor, micro / super cut, 45o, 19cm
- G142-760** Potts Smith Scissor, micro / super cut, 60o, 19cm

Scissors - Super Cut

G142-802 Micro / Super Cut, strabismus scissor, precision cutting blade / micro serrated blade, conical blade, round points, straight, 11cm

G142-805 Micro / Super Cut, strabismus scissor, precision cutting blade / micro serrated blade, conical blade, round points, curved, 11cm




Sinus Forceps

E078-105 Lister Sinus Forcep, serrated jaw, straight, 12.5cm

E078-110 Lister Sinus Forcep, serrated jaw, straight, 15cm

E078-115 Lister Sinus Forcep, serrated jaw, straight, 18cm

E078-120 Lister Sinus Forcep, serrated jaw, straight, 20cm



Splinter Forceps

G082-102 Arthur (Stieglitz) Splinter Forcep, straight, 14cm



G082-110 Hunter Splinter Forcep, straight, 11.5cm



G082-100 Martin (Feilchenfeld) Splinter Forcep, straight, 7.5cm

G082-105 Martin (Feilchenfeld) Splinter Forcep, straight, 9cm

G082-112 Martin (Feilchenfeld) Splinter Forcep, straight, 11.5cm



G082-216 Splinter / Delicate Tweezer Forcep, smooth, straight, 9cm

G082-218 Splinter / Delicate Tweezer Forcep, smooth, straight, 10.5cm


G082-220 Splinter / Delicate Tweezer Forcep, smooth, straight, 11.5cm

G082-222 Splinter / Delicate Tweezer Forcep, smooth, straight, 13cm




General

Sponge Holding Forceps



G084-126 Foerster Sponge Holder, curved, serrated jaw, 18cm

G084-128 Foerster Sponge Holder, curved, serrated jaw, 24cm



G084-118 Foerster Sponge Holder, straight, serrated jaw, 18cm

G084-124 Foerster Sponge Holder, straight, serrated jaw, 24cm



G084-218 Rampley Sponge Holder, straight, 18cm

G084-224 Rampley Sponge Holder, straight, 24cm

Steriliser Forceps




G084-305 Cheattle Sterilizer Forcep, fine

Suction Tubes

G150-010 American Pattern Suction Tube, chrome finish




G150-030 Millin Suction Tube with thumb control, chrome




G150-045 Poole Suction Tube, stainless steel, curved
23fg (7.66mm) diameter



G150-040 Poole Suction Tube, stainless steel, straight
30fg (10mm) diameter



G150-070 Simpson Smith Suction Tube, with locking
thumb control, stainless steel



N150-105 Frazier Suction Tube, 30o angle, 120mm
working length, 5fg (1.66mm dia)

N150-106 Frazier Suction Tube, 30o angle, 120mm
working length, 6fg (2mm dia)

N150-107 Frazier Suction Tube, 30o angle, 120mm
working length, 7fg (2.3mm dia)

N150-108 Frazier Suction Tube, 30o angle, 120mm
working length, 8fg (2.6mm dia)

N150-109 Frazier Suction Tube, 30o angle, 120mm
working length, 9fg (3mm dia)

N150-110 Frazier Suction Tube, 30o angle, 120mm
working length, 10fg (3.3mm dia)

N150-111 Frazier Suction Tube, 30o angle, 120mm
working length, 11fg (3.6mm dia)


N150-112 Frazier Suction Tube, 30o angle, 120mm
working length, 12fg (4mm dia)

N150-115 Frazier Suction Tube, 30o angle, 120mm
working length, 15fg (5mm dia)



General

Suction Tubes




- N150-205** Fergusson Suction Tube, 90o angle, 90mm working length, 5fg (1.66mm dia)
- N150-206** Fergusson Suction Tube, 90o angle, 90mm working length, 6fg (2mm dia)
- N150-207** Fergusson Suction Tube, 90o angle, 90mm working length, 7fg (2.3mm dia)
- N150-208** Fergusson Suction Tube, 90o angle, 90mm working length, 8fg (2.6mm dia)
- N150-209** Fergusson Suction Tube, 90o angle, 90mm working length, 9fg (3mm dia)
- N150-210** Fergusson Suction Tube, 90o angle, 90mm working length, 10fg (3.3mm dia)
- N150-211** Fergusson Suction Tube, 90o angle, 90mm working length, 11fg (3.6mm dia)
- N150-212** Fergusson Suction Tube, 90o angle, 90mm working length, 12fg (4mm dia)
- N150-215** Fergusson Suction Tube, 90o angle, 90mm working length, 15fg (5mm dia)




- G150-160** Yankauer Suction Tube, double angle, stainless steel
- G150-162** Yankauer Suction Tube, paediatric, double angle, stainless steel, 21.5cm

Suture Forceps



- G086-655** Guys Suture Forcep, 16cm
- G086-650** Guys Suture Forcep, 12.5cm

Tissue Forceps



- G088-070** American Pattern Allis (Judd Allis) Tissue Forcep, 3:4 teeth, 15cm




- G088-204** Poirier Tissue Forcep, 2:3 teeth, 12.5cm

Tissue Forceps


- G088-023** Allis Tissue Forcep, 15cm, 2:3 teeth
G088-034 Allis Tissue Forcep, 15cm, 3:4 teeth
G088-045 Allis Tissue Forcep, 15cm, 4:5 teeth
G088-056 Allis Tissue Forcep, 15cm, 5:6 teeth
G088-067 Allis Tissue Forcep, 15cm, 6:7 teeth




- G088-116** Babcock Tissue Forcep, 16cm
G088-120 Babcock Tissue Forcep, 20cm
G088-123 Babcock Tissue Forcep, 24cm




- G088-155** Duval Tissue Forcep, 19cm, 1.3cm wide
G088-160 Duval Lung / Tissue Forcep, 22cm, 2.5cm wide



- G088-170** Judd-Allis Tissue Forcep, 19cm, 3:4 teeth




- G088-206** Rutherford Morrison Tissue Forcep, 4:5 teeth, 16cm




General

Tissue Forceps




- G088-182** Lane Tissue Forcep, 1:2 teeth, 12.5cm
- G088-184** Lane Tissue Forcep, 1:2 teeth, 15cm
- G088-186** Lane Tissue Forcep, 1:2 teeth, 19.5cm



- G088-200** Littlewood Tissue Forcep, 2:3 teeth, 18.5cm

Tongue Forceps



- E080-205** Guy Tongue Holding Forcep, child, 12.5cm
- E080-210** Guy Tongue Holding Forcep, adult, 18cm

Towel Clips



- G090-202** Backhaus Towel Clip, 9cm
- G090-203** Backhaus Towel Clip, 10cm
- G090-204** Backhaus Towel Clip, 12.5cm
- G090-206** Backhaus Towel Clip, 15cm




- G090-290** Ball And Socket Towel Clip, 11cm
- G090-295** Ball And Socket Towel Clip, 14cm

Towel Clips


- G090-601** Lorna (Edna) Towel Clip, 9.5cm
G090-602 Lorna (Edna) Towel Clip, 13cm



- G090-605** Peers Towel Clip, 14.5cm




- G090-620** Tohoku Towel Clip, 10cm
G090-623 Tohoku Towel Clip, 13cm



- G090-538** Shardle Towel Clip, cross action, 9cm




- G090-600** Bernhard Towel Clip, 16cm



Bag Clips

- G090-100** Bag / X-Ray Clip, bulldog type, 2.5cm wide,
 stainless steel




- G090-150** Sling / Zimmer Clip, bulldog type, 6.5cm wide,
 serrated jaw




General

Vein Strippers




- G147-110** Graves Myers Vein Stripper Set, consisting of 1 each 6mm, 10mm, and 13mm interchangeable heads, t-handle and 900mm cable.
(does not include rack or case)




- G147-100** Nabatoff Vein Stripper Set, consisting of 1 x 15mm, 12mm, 9mm and 6mm interchangeable heads, 2 x 3mm probe tips, t-handle and 1 mtr cable complete in case

Vasectomy



- G068-030** Vasectomy Clamp, ring type, 3mm jaw, 14cm
G068-035 Vasectomy Clamp, ring type, 3.5mm jaw, 14cm
G068-040 Vasectomy Clamp, ring type, 4mm jaw, 14cm



- G068-050** Vasectomy (Quinsy) Forceps, sharp pointed tip, curved, 14cm

- S002-802** Gallipot, stainless steel, 55ml, 2oz, rimless, 50mm x 35mm
S002-804 Gallipot, stainless steel, 110ml, 4oz, rimless, 65mm x 40mm
S002-806 Gallipot, stainless steel, 170ml, 6oz, rimless, 80mm x 45mm
S002-810 Gallipot, stainless steel, 280ml, 10oz, rimless, 90mm x 55mm
S002-812 Gallipot, stainless steel, 340ml, 12oz, rimless, 105mm x 60mm



- S008-517** Kidney Dish, stainless steel, rimless, 16cm
S008-520 Kidney Dish, stainless steel, rimless, 20cm
S008-525 Kidney Dish, stainless steel, rimless, 25cm
S008-530 Kidney Dish, stainless steel, rimless, 30cm




- S010-404** Lotion Bowl, stainless steel, rimless, 240ml, 100mm
S010-414 Lotion Bowl, stainless steel, rimless, 700ml, 150mm
S010-424 Lotion Bowl, stainless steel, rimless, 2000ml, 205mm
S010-442 Lotion Bowl, stainless steel, rimless, 4500ml, 255mm
S010-452 Lotion Bowl, stainless steel, rimless, 6500ml, 305mm



- S006-600** Instrument Tray Only, stainless steel, 21cm x 15cm x 4cm
S006-605 Instrument Tray Only, stainless steel, 25cm x 20cm x 4cm
S006-610 Instrument Tray Only, stainless steel, 30cm x 23cm x 4cm
S006-615 Instrument Tray Only, stainless steel, 30cm x 25cm x 4cm
S006-620 Instrument Tray Only, stainless steel, 36cm x 30cm x 4cm



- S007-600** Instrument Tray Lid Only, stainless steel, for S006-600
S007-605 Instrument Tray Lid Only, stainless steel, for S006-605
S007-610 Instrument Tray Lid Only, stainless steel, for S006-610
S007-615 Instrument Tray Lid Only, stainless steel, for S006-615
S007-620 Instrument Tray Lid Only, stainless steel, for S006-620





Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk

We are an Approved Supplier to the
NHS Supply Chain



Supply Chain
v1