

Judd
medical

Neurosurgery Instrumentation Catalogue

www.judd-medical.co.uk

REUSABLE

We are an Approved
Supplier to the
NHS Supply Chain


NHS

Supply Chain

Index




	Pages
Artery Forceps	236-237
Braces	238
Burs	238
Curettes	238-240
Dissecting Forceps	240
Dissectors	241
Dural Hooks	241
Elevators	242-243
Gouges	243
Micro Instruments	253-254
Needleholders - Tungsten Carbide	244
Retractors	244-245
Retractors - Self Retaining	245-246
Rongeurs	247-252
Saws	252



G052-036 Cairns Artery Forcep, curved to side, 14.5cm



G052-034 Cairns (Providence) Artery Forcep, curved, 14cm



G052-032 Cairns (Providence) Artery Forcep, straight, 14cm



G052-041 Crile Artery Forcep, curved, 14cm

G052-043 Crile Artery Forcep, curved, 16cm

G052-045 Crile Artery Forcep, curved, 18cm


G052-040 Crile Artery Forcep, straight, 14cm

G052-042 Crile Artery Forcep, straight, 16cm


G052-044 Crile Artery Forcep, straight, 18cm




G052-054 Cushing Artery Forcep, curved to side, 14.5cm



G052-052 Cushing Artery Forcep, curved, 14.5cm



G052-050 Cushing Artery Forcep, straight, 14.5cm



Neurology

Braces



N040-013 Hudson (English Pattern) Brace, extension piece, 10cm



N040-010 Hudson (English Pattern) Brace, standard, 25cm

Burrs



N040-032 Hudson (English Pattern) Conical Burr, 9mm

N040-034 Hudson (English Pattern) Conical Burr, 11mm

N040-036 Hudson (English Pattern) Conical Burr, 13mm

N040-038 Hudson (English Pattern) Conical Burr, 16mm

N040-040 Hudson (English Pattern) Conical Burr, 19mm



N040-022 Hudson (English Pattern) Perforator Burr, 10mm

N040-024 Hudson (English Pattern) Perforator Burr, 12mm

N040-026 Hudson (English Pattern) Perforator Burr, 15mm




N040-042 Hudson (English Pattern) Spherical Burr, 11mm

N040-044 Hudson (English Pattern) Spherical Burr, 13mm

N040-046 Hudson (English Pattern) Spherical Burr, 16mm

N040-048 Hudson (English Pattern) Spherical Burr, 19mm

N040-050 Hudson (English Pattern) Spherical Burr, 22mm



N040-052 Spiral Burr, up cutting, 100mm x 12mm

Curettes



T090-140 Halle Flexible Curette, hollow handle, single ended, oval cup, 4.4mm wide, no.1, 21cm

T090-141 Halle Flexible Curette, hollow handle, single ended, oval cup, 5.2mm wide, no.2, 21cm

T090-142 Halle Flexible Curette, hollow handle, single ended, oval cup, 6.8mm wide, no.3, 21cm

- T090-120** Volkmann tufnol handle Curette, oval cup
4.4mm wide, no.0, 27cm
- T090-121** Volkmann tufnol handle Curette, oval cup
6.8mm wide, no.1, 27cm
- T090-122** Volkmann tufnol handle Curette, oval cup
8.5mm wide, no.2, 27cm
- T090-123** Volkmann tufnol handle Curette, oval cup
10mm wide, no.3, 27cm



- N014-020** Bruns Spinal Curette, hollow handle, straight,
oval cup 2mm wide, no.000, 23cm
- N014-022** Bruns Spinal Curette, hollow handle, straight,
oval cup 2.5mm wide, no.00, 23cm
- N014-024** Bruns Spinal Curette, hollow handle, straight,
oval cup 3.3mm wide, no.0, 23cm
- N014-026** Bruns Spinal Curette, hollow handle, straight,
oval cup 3.8mm wide, no.1, 23cm
- N014-028** Bruns Spinal Curette, hollow handle, straight,
oval cup 4.8mm wide, no.2, 23cm
- N014-030** Bruns Spinal Curette, hollow handle, straight,
oval cup 5.3mm wide, no.3, 23cm
- N014-032** Bruns Spinal Curette, hollow handle, straight,
oval cup 6mm wide, no.4, 23cm
- N014-034** Bruns Spinal Curette, hollow handle, straight,
oval cup 6.7mm wide, no.5, 23cm
- N014-036** Bruns Spinal Curette, hollow handle, straight,
oval cup 7.3mm wide, no.6, 23cm



- N014-110** Epstein Spinal Curette, hollow handle, reverse angle,
no.000, 20cm
- N014-111** Epstein Spinal Curette, hollow handle, reverse angle,
no.00, 20cm
- N014-112** Epstein Spinal Curette, hollow handle, reverse angle,
no.0, 20cm



- N014-120** Spinal Curette, hollow handle, straight, no.000,
20cm
- N014-121** Spinal Curette, hollow handle, straight, no.00,
20cm
- N014-122** Spinal Curette, hollow handle, straight, no.0,
20cm
- N014-123** Spinal Curette, hollow handle, straight, no.1,
20cm
- N014-124** Spinal Curette, hollow handle, straight, no.2,
20cm
- N014-125** Spinal Curette, hollow handle, straight, no.3,
20cm
- N014-126** Spinal Curette, hollow handle, straight, no.4,
20cm
- N014-127** Spinal Curette, hollow handle, straight, no.5,
20cm
- N014-128** Spinal Curette, hollow handle, straight, no.6,
20cm



Curettes



N014-130 Spinal Curette, hollow handle, angled, no.000, 20cm

N014-131 Spinal Curette, hollow handle, angled, no.00, 20cm

N014-132 Spinal Curette, hollow handle, angled, no.0, 20cm

N014-133 Spinal Curette, hollow handle, angled, no.1, 20cm

N014-134 Spinal Curette, hollow handle, angled, no.2, 20cm

N014-135 Spinal Curette, hollow handle, angled, no.3, 20cm

N014-136 Spinal Curette, hollow handle, angled, no.4, 20cm

N014-137 Spinal Curette, hollow handle, angled, no.5, 20cm

N014-138 Spinal Curette, hollow handle, angled, no.6, 20cm

Dissecting Forceps



G062-118 Adson Bayonet Dissecting Forcep, 1:2 teeth, 19.5cm



G062-108 Adson Bayonet Dissecting Forcep, serrated, 19.5cm



G062-106 Adson Dissecting Forcep, serrated, 18cm



G062-116 Adson Dissecting Forcep, 1:2 teeth, 18cm



G062-170 Cushing Dissecting Forcep, serrated, straight, 18cm



G062-172 Cushing Dissecting Forcep, 1:2 teeth, straight, 18cm



G062-416 Potts Smith Dissecting Forcep, straight, serrated, 18cm

G062-418 Potts Smith Dissecting Forcep, straight, serrated, 21cm

G062-420 Potts Smith Dissecting Forcep, straight, serrated, 25cm




G062-417 Potts Smith Dissecting Forcep, straight, 1:2 teeth, 18cm

G062-419 Potts Smith Dissecting Forcep, straight, 1:2 teeth, 21cm

G062-421 Potts Smith Dissecting Forcep, straight, 1:2 teeth, 25cm

Dissectors

N045-700 Pennybacker Dissector Elevator and Probe, double ended, 22.5cm



N045-602 Penfield Dissector, double ended, broad curved disscetor and sharp round spoon, no.1, 18cm



N045-604 Penfield Dissector, double ended, slightly curved dissector and wax packer, no.2, 19.5cm



N045-606 Penfield Dissector, double ended, fully curved dissector and wax packer, no.3, 19.5cm



N045-608 Penfield Dissector, single ended, slightly curved dissector, no.4, 21cm



N045-610 Penfield Dissector, double ended, slightly curved blunt dissectors, no.5, 29cm



Dural Hooks

N104-310 Adson Nerve / Vein Hook, blunt, 19cm


N104-312 Adson Nerve / Vein Hook, sharp, 19cm



N104-400 Cairns Dural Hook, fine, 12.5cm



N104-412 Cairns Dural Hook, 90o, sharp, with flat forged handle, 12.5cm



M135-250 Blalock Hook, 90o blunt hook, 24cm

M135-252 Blalock Hook, 90o sharp hook, 24cm



M135-350 Nelson Nerve / Vein Hook, retrograde curve, blunt, 24cm



Elevators



N060-014 Adson Periosteal Elevator, curved, 5mm, 16.5cm

N060-015 Adson Periosteal Elevator, curved, 7mm, 16.5cm



N060-010 Adson Periosteal Elevator, straight, 5mm, 16.5cm

N060-012 Adson Periosteal Elevator, straight, 7mm, 16.5cm



N050-210 Cobb Spinal Elevator, 10mm, 28cm

N050-213 Cobb Spinal Elevator, 13mm, 28cm

N050-219 Cobb Spinal Elevator, 19mm, 28cm



N050-220 Cobb Spinal Elevator, tufnol handle, 10mm wide, 25cm

N050-223 Cobb Spinal Elevator, tufnol handle, 13mm wide, 25cm

N050-229 Cobb Spinal Elevator, tufnol handle, 19mm wide, 25cm

N050-235 Cobb Spinal Elevator, tufnol handle, 25mm wide, 25cm



N050-404 Wagner Periosteal Raspator, tufnol handle, 8mm wide, 33cm

N050-408 Wagner Periosteal Raspator, tufnol handle, 17mm wide, 33cm



N060-250 Pennybacker Elevator, curved, sharp, 6mm, 17cm



Elevators

N350-502 Periosteal Raspatory, tufnol handle, straight shaft, 6mm wide, straight edge, 19cm

N350-504 Periosteal Raspatory, tufnol handle, curved, 6mm straight edge, 18.5cm

N350-506 Periosteal Raspatory, tufnol handle, curved, 6mm round edge, 18.5cm

N350-513 Periosteal Raspatory, tufnol handle, curved, 3mm straight edge, 18.5cm

N350-522 Periosteal Raspatory, tufnol handle, curved shaft, 13mm wide straight edge, 18.5cm

N350-523 Periosteal Raspatory, tufnol handle, straight shaft, 13mm wide straight edge, 18.5cm

N350-524 Periosteal Raspatory, tufnol handle, straight shaft, 14mm round edge, 19.5cm

N350-526 Periosteal Raspatory, tufnol handle, curved shaft, 14mm round edge, 18.5cm

N350-530 Periosteal Raspatory, tufnol handle, straight shaft, 20mm round edge, 20cm



Gouges

N350-086 Capener Laminar Gouge, tufnol handle, cranked stem, 6mm wide, 25cm

N350-088 Capener Laminar Gouge, tufnol handle, cranked stem, 8mm wide, 25cm

N350-091 Capener Laminar Gouge, tufnol handle, cranked stem, 11mm wide, 25cm

N350-096 Capener Laminar Gouge, tufnol handle, cranked stem, 16mm wide, 25cm



N350-046 Capener Gouge, straight, tufnol handle, 6mm wide, 25cm

N350-048 Capener Gouge, straight, tufnol handle, 8mm wide, 25cm

N350-050 Capener Gouge, straight, tufnol handle, 10mm wide, 25cm

N350-052 Capener Gouge, straight, tufnol handle, 12mm wide, 25cm

N350-055 Capener Gouge, straight, tufnol handle, 15mm wide, 25cm


N350-070 Capener Gouge, straight, tufnol handle, 20mm wide, 25cm

N350-075 Capener Gouge, straight, tufnol handle, 25mm wide, 25cm

N350-077 Capener Gouge, straight, tufnol handle, 30mm wide, 25cm




Needle Holders - Tungsten Carbide




- G126-141** Olsen Hegar Needle Holder & Scissor, tungsten carbide inserts, 14cm
- G126-145** Olsen Hegar Needle Holder & Scissor, tungsten carbide inserts, 19cm
- G126-143** Olsen Hegar Needle Holder & Scissor, tungsten carbide inserts, 16.5cm


Retractors



- N134-041** Cairns Strip Retractor, malleable, 6mm
- N134-042** Cairns Strip Retractor, malleable, 8mm
- N134-043** Cairns Strip Retractor, malleable, 10mm
- N134-044** Cairns Strip Retractor, malleable, 11mm
- N134-045** Cairns Strip Retractor, malleable, 13mm
- N134-046** Cairns Strip Retractor, malleable, 16mm
- N134-047** Cairns Strip Retractor, malleable, 19mm
- N134-048** Cairns Strip Retractor, malleable, 25mm
- N134-040** Cairn Strip Retractors, set/ 8, malleable, sizes 6, 8, 10, 11, 13, 16, 19 and 25mm wide, 20cm



- N134-202** Charnley Laminectomy Retractor, single sharp blade only, 70mm deep
- N134-204** Charnley Laminectomy Retractor, double blunt blade only, 70mm deep
- N134-200** Charnley Laminectomy Retractor set, 1 x sharp and blunt blade, 70mm deep, complete with weight and chain



- N104-250** Cushing Nerve / Vein Hook, right angle, 19cm
- N104-252** Cushing Nerve / Vein Hook, right angle, 23cm
- N104-255** Cushing Nerve / Vein Hook, right angle, 28cm




- N134-240** Hibbs Retractor, 25mm x 75mm claw blade, 24cm
- N134-242** Hibbs Retractor, 25mm x 80mm claw blade, 24cm
- N134-244** Hibbs Retractor, toothed, set /2, 24cm

N134-750 Love Nerve Root Retractor, straight, 20cm

N134-753 Love Nerve Root Retractor, 45o, 20cm

N134-755 Love Nerve Root Retractor, 90o, 20cm



N134-820 O'Connell Nerve Root Retractor, 20cm



Retractors - Self Retaining

N136-090 Cone (Adson - Baby) hinged retractor, 3:4 teeth, sharp, 14cm

N136-092 Cone (Adson - Baby) hinged retractor, 3:4 teeth, sharp, 16cm

N136-094 Cone Hinged Retractor, 3:4 teeth, sharp, 20cm

N136-096 Cone Hinged Retractor, 3:4 teeth, sharp, 26cm

N136-080 Cone (Adson - Baby) hinged retractor, 3:4 teeth, blunt, 14cm

N136-082 Cone (Adson - Baby) hinged retractor, 3:4 teeth, blunt, 16cm

N136-084 Cone Hinged Retractor, 3:4 teeth, blunt, 20cm

N136-086 Cone Hinged Retractor, 3:4 teeth, blunt, 26cm



G136-750 Weitlaner Retractor, 2:3 sharp prongs, 11cm


G136-755 Weitlaner Retractor, 2:3 blunt prongs, 11cm

G136-731 Weitlaner Retractor, self retaining, straight arms, 3:4 sharp prongs, 14cm

G136-732 Weitlaner Retractor, self retaining, straight arms, 3:4 sharp prongs, 16.5cm


G136-735 Weitlaner Retractor, self retaining, straight arms, 3:4 blunt prongs, 14cm

G136-736 Weitlaner Retractor, self retaining, straight arms, 3:4 blunt prongs, 16.5cm



G136-313 Hardman Laminectomy Retractor, short left 30mm, long right 50mm, 21.5cm


G136-312 Hardman Laminectomy Retractor, short right 30mm, long left 50mm, 21.5cm




Retractors - Self Retaining



G136-310 Harvey Jackson Laminectomy Self Retaining Retractor, 30.5cm



N136-250 Henderson Self Retaining Retractor, 13cm



G136-595 Jefferson (Scalp - Contour) Self Retaining Retractor, 3:4 blunt teeth, 14cm



G136-110 Beckmann-Eaton Retractor, self retaining, hinged arms, 7:7 sharp prongs, 33cm




G136-112 Beckmann-Adson Self Retaining Retractor, hinged arms, 25mm blade depth, 4:4 sharp teeth, 31cm

G136-113 Beckmann-Adson Self Retaining Retractor, hinged arms, 44mm blade depth, 4:4 sharp teeth, 31cm

G136-114 Beckmann-Adson Self Retaining Retractor, hinged arms, 25mm blade depth, 4:4 blunt teeth, 31cm

G136-115 Beckmann-Adson Self Retaining Retractor, hinged arms, 44mm blade depth, 4:4 blunt teeth, 31cm

Intervertebral Disc Rongeurs

LENGTH OF SHAFT 1.5MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

Straight Jaw

15cm (6")	N074-197	N074-200	N074-230	N074-260	N073-060	N073-065
18cm (7")	N074-207	N074-210	N074-240	N074-270	N073-070	N073-075
20cm (8")	N073-007	N073-020	N073-025	N074-275	N073-080	N073-085
23cm (9")	N073-010	N073-030	N073-035	N074-280	N073-090	N073-095
25cm (10")	-----	N073-040	N073-045	N074-285	N073-100	N073-105
28cm (11")	-----	N073-050	N073-055	N074-290	N073-110	N073-115
30cm (12")	-----	N074-217	N074-247	N074-297	N073-120	N073-125

LENGTH OF SHAFT 1.5MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

Angled Upwards

15cm (6")	N074-198	N074-202	N074-232	N074-262	N073-062	N073-067
18cm (7")	N074-208	N074-212	N074-242	N074-272	N073-072	N073-077
20cm (8")	N073-008	N073-022	N073-027	N074-276	N073-082	N073-087
23cm (9")	N073-011	N073-032	N073-037	N074-281	N073-092	N073-097
25cm (10")	-----	N073-042	N073-047	N074-286	N073-102	N073-107
28cm (11")	-----	N073-052	N073-057	N074-291	N073-112	N073-117
30cm (12")	-----	N074-218	N074-248	N074-298	N073-122	N073-127

LENGTH OF SHAFT 1.5MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

Angled Downwards

15cm (6")	N074-199	N074-203	N074-233	N074-263	N073-064	N073-069
18cm (7")	N074-209	N074-213	N074-243	N074-273	N073-074	N073-079
20cm (8")	N073-009	N073-024	N073-029	N074-277	N073-084	N073-089
23cm (9")	N073-012	N073-034	N073-039	N074-282	N073-094	N073-099
25cm (10")	-----	N073-044	N073-049	N074-287	N073-104	N073-109
28cm (11")	-----	N073-054	N073-059	N074-292	N073-114	N073-119
30cm (12")	-----	N074-219	N074-249	N074-299	N073-124	N073-129



Ferris Smith Kerrison Rongeur



Upward Cutting, standard handle, 10mm jaw opening.

LENGTH OF SHAFT 1MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

40° UP BITE




18cm (7")	N074-601	N074-602	N074-603	N074-604	N074-605	N074-606
20cm (8")	N074-721	N074-722	N074-723	N074-724	N074-725	N074-726
23cm (9")	N074-731	N074-732	N074-733	N074-734	N074-735	N074-736
25cm (10")	N074-741	N074-742	N074-743	N074-744	N074-745	N074-746
28cm (11")	N074-771	N074-772	N074-773	N074-774	N074-775	N074-776
30cm (12")	N074-761	N074-762	N074-763	N074-764	N074-765	N074-766

LENGTH OF SHAFT 1MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

90° UP BITE



18cm (7")	N074-501	N074-502	N074-503	N074-504	N074-505	N074-506
20cm (8")	N074-701	N074-702	N074-703	N074-704	N074-705	N074-706
23cm (9")	N074-781	N074-782	N074-783	N074-784	N074-785	N074-786
25cm (10")	N074-801	N074-802	N074-803	N074-804	N074-805	N074-806
28cm (11")	N074-831	N074-832	N074-833	N074-834	N074-835	N074-836
30cm (12")	N074-841	N074-842	N074-843	N074-844	N074-845	N074-846



Rongeurs also available with Bone Ejector. Catalogue number as above ending with an 'E'

Ferris Smith Kerrison Rongeur**Downward Cutting, standard handle, 10mm jaw opening.**

LENGTH OF SHAFT 1MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE


40° DOWN BITE

18cm (7")	N074-621	N074-622	N074-623	N074-624	N074-625	N074-626
20cm (8")	N074-631	N074-632	N074-633	N074-634	N074-635	N074-636
23cm (9")	N074-671	N074-672	N074-673	N074-674	N074-675	N074-676
25cm (10")	N074-681	N074-682	N074-683	N074-684	N074-685	N074-686
28cm (11")	N074-851	N074-852	N074-853	N074-854	N074-855	N074-856
30cm (12")	N074-861	N074-862	N074-863	N074-864	N074-865	N074-866

LENGTH OF SHAFT 1MM BITE 2MM BITE 3MM BITE 4MM BITE 5MM BITE 6MM BITE

90° DOWN BITE

18cm (7")	N074-511	N074-512	N074-513	N074-514	N074-515	N074-516
20cm (8")	N074-711	N074-712	N074-713	N074-714	N074-715	N074-716
23cm (9")	N074-811	N074-812	N074-813	N074-814	N074-815	N074-816
25cm (10")	N074-871	N074-872	N074-873	N074-874	N074-875	N074-876
28cm (11")	N074-881	N074-882	N074-883	N074-884	N074-885	N074-886
30cm (12")	N074-891	N074-892	N074-893	N074-894	N074-895	N074-896



Rongeurs also available with Bone Ejector. Catalogue number as above ending with an 'E'





N076-110 Glasgow Rongeur, curved, 4mm, 14.5cm



N076-130 Jansen (Heanley) Bone Rongeur, straight, 17.5cm


N076-133 Jansen Bone Rongeur, curved, 17.5cm




N076-350 Luer Rongeur, heavy pattern, straight, 17cm

N076-355 Luer Rongeur, heavy pattern, curved, 17cm

N076-004 Adson Bone Rongeur, straight, 20cm



N076-470 Northfield Rongeur, straight, 20cm



N076-600 Pennybacker Rongeur, straight, 4mm, 20cm



N076-760 Ruskin Rongeur, compound action, curved, 5mm bite, 18cm

N076-764 Ruskin Rongeur, compound action, straight, 5mm bite, 18cm



Rongeurs



N077-320 Stille Rongeur, angled, compound action, 6mm bite, 23cm




N077-330 Stille Ruskin Rongeur, curved, compound action, 23cm

Saws




T140-130 De Martel Gigli Saw Guide, 30cm



T140-120 Gigli Saw Handle, pair

- T140-122** Gigli Saw Wires, 30cm, pk/5
T140-124 Gigli Saw Wires, 50cm, pk/5
T140-125 Gigli Saw Wires, 60cm, pk/5



Various - Micro Instruments

- F275-105** Micro Forcep, wide flat handle, 0.3mm smooth straight jaw, 15cm
F275-106 Micro Forcep, wide flat handle, 0.3mm smooth curved jaw, 15cm
F275-110 Micro Forcep, wide flat handle, 0.3mm smooth straight jaw, 18cm
F275-111 Micro Forcep, wide flat handle, 0.3mm smooth curved jaw, 18cm
F275-115 Micro Forcep, wide flat handle, 0.7mm smooth straight jaw, 15cm
F275-116 Micro Forcep, wide flat handle, 0.7mm smooth curved jaw, 15cm
F275-120 Micro Forcep, wide flat handle, 0.7mm smooth straight jaw, 18cm
F275-121 Micro Forcep, wide flat handle, 0.7mm smooth curved jaw, 18cm
F275-130 Micro Forcep, wide flat handle, 0.3mm smooth straight jaw, 21cm
F275-131 Micro Forcep, wide flat handle, 0.3mm smooth curved jaw, 21cm



- F275-155** Micro Needle Holder, wide flat handle, smooth jaw, without catch, straight, 14cm
F275-160 Micro Needle Holder, wide flat handle, smooth jaw, without catch, curved, 14cm
F275-165 Micro Needle Holder, wide flat handle, smooth jaw, without catch, straight, 18cm
F275-170 Micro Needle Holder, wide flat handle, smooth jaw, without catch, curved, 18cm
F275-171 Micro Needle Holder, wide flat handle, smooth jaw, without catch, straight, 21.5cm
F275-172 Micro Needle Holder, wide flat handle, smooth jaw, without catch, curved, 21.5cm
F275-173 Micro Needle Holder, wide flat handle, smooth jaw, without catch, straight, 23.5cm
F275-174 Micro Needle Holder, wide flat handle, smooth jaw, without catch, curved, 23.5cm





- F275-175** Micro Needle Holder, wide flat handle, smooth jaw, with catch, straight, 14cm
- F275-180** Micro Needle Holder, wide flat handle, smooth jaw, with catch, curved, 14cm
- F275-185** Micro Needle Holder, wide flat handle, smooth jaw, with catch, straight, 18cm
- F275-190** Micro Needle Holder, wide flat handle, smooth jaw, with catch, curved, 18cm
- F275-191** Micro Needle Holder, wide flat handle, smooth jaw, with catch, straight, 21.5cm
- F275-192** Micro Needle Holder, wide flat handle, smooth jaw, with catch, curved, 21.5cm
- F275-193** Micro Needle Holder, wide flat handle, smooth jaw, with catch, straight, 23.5cm
- F275-194** Micro Needle Holder, wide flat handle, smooth jaw, with catch, curved, 23.5cm



- F275-202** Micro Scissor, spring pattern, flat handle, straight, 12cm
- F275-203** Micro Scissor, spring pattern, flat handle, curved, 12cm
- F275-205** Micro Scissor, spring pattern, flat handle, straight, 15cm
- F275-210** Micro Scissor, spring pattern, flat handle, curved, 15cm
- F275-215** Micro Scissor, spring pattern, flat handle, straight, 18cm
- F275-220** Micro Scissor, spring pattern, flat handle, curved, 18cm



Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk

We are an Approved Supplier to the
NHS Supply Chain



Supply Chain
v1