

Judd
medical

Orthopaedic Consumables Catalogue

www.judd-medical.co.uk

NHS

Supply Chain

APPROVED SUPPLIER

TABLE OF CONTENTS

3	Who We Are / How Can We Help
4	What We Do / Training & Education
5	How To Order / Delivery Options
6-7	Kirschner Wires
8-9	Arthrodesis Wires
10	Steinmann & Denham Pins
11	Cerclage Wire
12	Quick Coupling Twist Drills
13	Standard Twist Drills
14	Pulse Lavage System
15	Cement Mixer Systems
16	Hip Revision Extraction Set
17-19	Headlights, Light Sources, Cables & Adaptors
20-26	Instrumentation (Reusable)
27	Plaster Room Instruments (Reusable)
28	Orthopaedic Instrumentation Catalogue (Reusable)
29	Saw Blades
30	Smoke Evacuation Pencil
31	Electrosurgery: Energy & Diathermy Catalogue
32	Catalogue Resources & Product Videos
33	Quality & Accreditations
34	Notes
35	How To Reach Us

Disclaimer: Product specifications in this catalogue are correct at time of publication but may change from time-to-time without notice. It is the customer's responsibility to validate that a particular product with the properties described in the product specification meets the requirements of their application.

WHO WE ARE



Judd Medical Limited is recognised as being amongst one of the most progressive medical companies in the UK and during the last 25 years have seen numerous developments that have enabled us to grow into a leading medical device supplier in the industry.

With the expansion of our product portfolio to include many innovative products and services, it has allowed Judd Medical to make a significant impact on the healthcare market to NHS hospitals, private groups and independent clinics throughout the UK.

HOW CAN WE HELP



CUSTOMER SERVICE LEVELS

We ensure that should you need any help and advice or encounter any problems we have a fully trained team ready to support you. Our specifically trained Customer Service Team offer friendly, knowledgeable advice and support for every requirement.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any products proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts, standing and call off orders dependent on your requirements.

JUDD MEDICAL

“YOUR PARTNER IN PRECISION HEALTHCARE”

WHAT WE DO.....



Our company is proud to offer a comprehensive product portfolio of
Single Use & Reusable Surgical Instruments and Medical Devices.

-  Arthrodesis & Kirschner Wires, Twist Drills, Steinmann Pins & Denham Pins
-  Electrosurgery: Diathermy Products – Single Use & Reusable
-  Light Sources: Xenon & LED, Surgical Headlights, Fibre Light Cables & Adaptors
-  Pulse Lavage Systems
-  Replacement Surgical Saw Blades & Burs
-  Single-Use & Reusable Surgical Instrumentation

At Judd Medical we are committed to our partnership with the NHS Supply Chain, NHS Scotland and NHS Wales, ensuring optimal value for our NHS Customers.

We are an approved supplier to the NHS Supply Chain.



TRAINING & EDUCATION



We are dedicated to providing comprehensive training and education on regulatory compliance, organising international seminars and hosting healthcare roundtable events. We facilitate peer-to-peer clinical discussions and offer training on the care and handling of surgical instrumentation and diathermy safety.

Through these initiatives, we aim to enhance professional knowledge, ensure adherence to industry standards and promote the highest levels of safety and efficacy in healthcare practices.

HOW TO ORDER:

You can place an order with us via:



01527 559010



sales@judd-medical.co.uk



NHS Framework Agreements



SCAN TO VIEW OUR
NHS FRAMEWORK
AGREEMENTS

DELIVERY OPTIONS:



Shipping & Despatch – Where stock is available this will be despatched on the same day.
(Order must be processed before 3pm)

Delivery options available:

- Standard Next Day Delivery
- Pre 9am Next Day Delivery*
- Pre 10.30am Next Day Delivery*
- Pre 12 Noon Next Day Delivery*
- Same Day Delivery – Must be processed before 12 noon on date of request.*

*Price on Application

We deliver to:

- Mainland UK
- Northern Ireland
- Channel Islands
- Scottish Highlands & Islands



**Centrally located in Bromsgrove,
Worcestershire**



Kirschner Wires, Single Trocar

STERILE packed disposables for use in
Orthopaedic & Trauma Surgery.



Product Code	Supply Chain Code	Product Description	Pk Size
0.9mm diameter			
362.100S	FQD5968	0.9mm x 10cm K Wire, Single Trocar	10
362.125S	FQD5969	0.9mm x 12.5cm K Wire, Single Trocar	10
362.150S	FQD5970	0.9mm x 15cm K Wire, Single Trocar	10
1.0mm diameter			
320.070S	FQD5971	1.0mm x 7cm K Wire, Single Trocar	10
320.100S	FQD5972	1.0mm x 10cm K Wire, Single Trocar	10
320.150S	FQD5973	1.0mm x 15cm K Wire, Single Trocar	10
1.1mm diameter			
363.100S	FQD5975	1.1mm x 10cm K Wire, Single Trocar	10
363.125S	FQD5976	1.1mm x 12.5cm K Wire, Single Trocar	10
363.150S	FQD5977	1.1mm x 15cm K Wire, Single Trocar	10
1.2mm diameter			
321.225S	FQD5978	1.2mm x 22.5cm K Wire, Single Trocar	10
321.300S	FQX88834	1.2mm x 30cm K Wire, Single Trocar	10
1.25mm diameter			
346.150S	FQD5980	1.25mm x 15cm K Wire, Single Trocar	10
1.4mm diameter			
322.100S	FQD5981	1.4mm x 10cm K Wire, Single Trocar	10
322.125S	FQD5982	1.4mm x 12.5cm K Wire, Single Trocar	10
322.150S	FQD5983	1.4mm x 15cm K Wire, Single Trocar	10
1.5mm diameter			
364.150S	FQD5984	1.5mm x 15cm K Wire, Single Trocar	10


Kirschner Wires, Single Trocar

All wires are:


- Individually **STERILE** packed in tough Tyvek® pouches
- Available in additional sizes upon request



Product Code	Supply Chain Code	Product Description	Pk Size
1.6mm diameter			
323.100S	FQX88825	1.6mm x 10cm K Wire, Single Trocar	10
323.125S	FQX88827	1.6mm x 12.5cm K Wire, Single Trocar	10
323.150S	FQD5988	1.6mm x 15cm K Wire, Single Trocar	10
323.200S	FQD5989	1.6mm x 20cm K Wire, Single Trocar	10
323.230S	FQD5921	1.6mm x 23cm K Wire, Single Trocar	10
1.8mm diameter			
325.150S	FQD5922	1.8mm x 15cm K Wire, Single Trocar	10
2.0mm diameter			
324.100S	Direct Only	2.0mm x 10cm K Wire, Single Trocar	10
324.125S	FQD5923	2.0mm x 12.5cm K Wire, Single Trocar	10
324.150S	FQD5924	2.0mm x 15cm K Wire, Single Trocar	10
324.200S	FQD5925	2.0mm x 20cm K Wire, Single Trocar	10
324.230S	FQD5926	2.0mm x 23cm K Wire, Single Trocar	10
2.5mm diameter			
335.150S	FQD5927	2.5mm x 15cm K Wire, Single Trocar	10
335.175S	FQX88826	2.5mm x 17.5cm K Wire, Single Trocar	10
335.250S	FQX88824	2.5mm x 25cm K Wire, Single Trocar	10



Arthrodesis Wires, Double Trocar



STERILE packed disposables for use in
Orthopaedic & Trauma Surgery.

Product Code	Supply Chain Code	Product Description	Pk Size
338.100S	FQX88832	0.8mm x 10cm A Wire, Double Trocar	10
370.075S	FQD5992	0.9mm x 7.5cm A Wire, Double Trocar	10
370.100S	FQD5993	0.9mm x 10cm A Wire, Double Trocar	10
370.125S	FQD5994	0.9mm x 12.5cm A Wire, Double Trocar	10
370.150S	FQD5996	0.9mm x 15cm A Wire, Double Trocar	10
339.070S	FQD5997	1.0mm x 7cm A Wire, Double Trocar	10
339.100S	FQD5998	1.0mm x 10cm A Wire, Double Trocar	10
339.125S	Direct Only	1.0mm x 12.5cm A Wire, Double Trocar	10
339.150S	FQD5999	1.0mm x 15cm A Wire, Double Trocar	10
371.070S	FQD6001	1.1mm x 7cm A Wire, Double Trocar	10
371.075S	FQD6000	1.1mm x 7.5cm A Wire, Double Trocar	10
371.100S	FQX88833	1.1mm x 10cm A Wire, Double Trocar	10
371.125S	FQD6002	1.1mm x 12.5cm A Wire, Double Trocar	10
371.130S	FQD6003	1.1mm x 13cm A Wire, Double Trocar	10
371.150S	FQD6004	1.1mm x 15cm A Wire, Double Trocar	10
371.200S	FQD6005	1.1mm x 20cm A Wire, Double Trocar	10
340.100S	FQD6006	1.2mm x 10cm A Wire, Double Trocar	10
340.125S	FQD6007	1.2mm x 12.5cm A Wire, Double Trocar	10
340.150S	FQD6008	1.2mm x 15cm A Wire, Double Trocar	10
345.075S	FQD6009	1.25mm x 7.5cm A Wire, Double Trocar	10
345.100S	FQD6010	1.25mm x 10cm A Wire, Double Trocar	10

Arthrodesis Wires, Double Trocar

All wires are:

- Individually **STERILE** packed in tough Tyvek® pouches
- Available in additional sizes upon request



Product Code	Supply Chain Code	Product Description	Pk Size
345.125S	FQX88829	1.25mm x 12.5cm A Wire, Double Trocar	10
345.150S	FQD6011	1.25mm x 15cm A Wire, Double Trocar	10
341.100S	FQD6013	1.4mm x 10cm A Wire, Double Trocar	10
341.125S	FQD6014	1.4mm x 12.5cm A Wire, Double Trocar	10
341.150S	FQD6015	1.4mm x 15cm A Wire, Double Trocar	10
342.075S	FQD5913	1.6mm x 7.5cm A Wire, Double Trocar	10
342.100S	FQD5914	1.6mm x 10cm A Wire, Double Trocar	10
342.125S	FQD5915	1.6mm x 12.5cm A Wire, Double Trocar	10
342.130S	FQD5916	1.6mm x 13cm A Wire, Double Trocar	10
342.150S	FQD8225	1.6mm x 15cm A Wire, Double Trocar	10
342.200S	FQD5919	1.6mm x 20cm A Wire, Double Trocar	10
342.230S	FQD5920	1.6mm x 23cm A Wire, Double Trocar	10
343.125S	FQD5955	1.8mm x 12.5cm A Wire, Double Trocar	10
343.150S	FQD5956	1.8mm x 15cm A Wire, Double Trocar	10
344.100S	FQD5958	2.0mm x 10cm A Wire, Double Trocar	10
344.125S	FQD5959	2.0mm x 12.5cm A Wire, Double Trocar	10
344.150S	FQD5961	2.0mm x 15cm A Wire, Double Trocar	10
344.200S	FQD5963	2.0mm x 20cm A Wire, Double Trocar	10
344.230S	FQD5964	2.0mm x 23cm A Wire, Double Trocar	10



The **Steinmann Pins** listed here are the most popular sizes in our range.

All our pins are:

- Individually **STERILE** packed in tough Tyvek® pouches
- Available in additional sizes upon request



Product Code	Supply Chain Code	Product Description	Pk Size
Steinmann Pins, Plain/Trocar, STERILE			
334.100S	Direct Only	2.4mm x 100mm Steinmann Pin, Pl/Tr	10
334.200S	FQX85007	2.4mm x 200mm Steinmann Pin, Pl/Tr	10
309.230S	FQX85010	2.8mm x 230mm Steinmann Pin, Pl/Tr	10
312.200S	FQX85020	3.2mm x 200mm Steinmann Pin, Pl/Tr	10
312.230S	FQX88808	3.2mm x 230mm Steinmann Pin, Pl/Tr	10
336.200S	FQX88811	3.5mm x 200mm Steinmann Pin, Pl/Tr	10
336.230S	FQX85006	3.5mm x 230mm Steinmann Pin, Pl/Tr	10
500.200S	FQX85026	4.8mm x 200mm Steinmann Pin, Pl/Tr	10
500.250S	FQX85021	4.8mm x 250mm Steinmann Pin, Pl/Tr	10

Steinmann Pins, Plain/Triangle, Fully Threaded, STERILE			
380.150S	FQX85001	1.6mm x 150mm Steinmann Pin, Pl/Tr, Fully Threaded	10
381.150S	FQX85023	2.0mm x 150mm Steinmann Pin, Pl/Tr, Fully Threaded	10
383.200S	FQX85014	3.0mm x 200mm Steinmann Pin, Pl/Tr, Fully Threaded	10
383.230S	Direct Only	3.0mm x 230mm Steinmann Pin, Pl/Tr, Fully Threaded	10

Product Code	Supply Chain Code	Product Description	Ends	Diameter	Length	Pk Size
Denham Pins, STERILE						
390.200S	FQX85019	Mid Thread	Trocar/3sq	4.0mm	20cm	10
390.230S	FQX85008	Mid Thread	Trocar/3sq	4.0mm	23cm	10
391.200S	Direct Only	Mid Thread	Trocar/3sq	4.8mm	20cm	10
391.230S	FQX85012	Mid Thread	Trocar/3sq	4.8mm	23cm	10
397.200S	FQX85016	Mid Thread	Trocar/3sq	5.0mm	20cm	10
397.230S	FQX85005	Mid Thread	Trocar/3sq	5.0mm	23cm	10

T520-010	Direct Only	Steinmann Pin Cap, Reusable					1
----------	-------------	-----------------------------	--	--	--	--	---



CERCLAGE WIRE


Sterile

Pre-Cut **STERILE** Cerclage Wire

- Ready to Use
- No Wastage
- No need to Re-Sterilise
- Eliminates the need to use Reel Wire
- Manufactured from annealed 316lvm grade wire
- Individually double pouched in tough Tyvek® pouches
- Three peel off adhesive patient labels per pouch allowing for full product traceability
- Sterilised by gamma irradiation to ISO 11137



Product Code	Supply Chain Code	Product Description	Pk Size
720.400S	FQX88813	Cerclage Bone Wire 0.40mm x 400mm	10
730.500S	FQX85000	Cerclage Bone Wire 0.8mm x 500mm	10
731.400S	FQX88812	Cerclage Bone Wire 0.9mm x 400mm	10
735.500S	Direct Only	Cerclage Bone Wire 1.0mm x 500mm	10
741.400S	FQX88815	Cerclage Bone Wire 1.5mm x 400mm	10
741.500S	FQX88817	Cerclage Bone Wire 1.5mm x 500mm	10
750.090S	Direct Only	Cerclage Bone Wire 1.2mm x 900mm (Looped)	10
750.400S	FQX88807	Cerclage Bone Wire 1.2mm x 400mm	10
750.500S	FQX85024	Cerclage Bone Wire 1.2mm x 500mm	10
750.600S	FQX88816	Cerclage Bone Wire 1.2mm x 600mm	10
761.400S	Direct Only	Cerclage Bone Wire 0.7mm x 400mm	10
790.400S	FQX88806	Cerclage Bone Wire 2.0mm x 400mm	10



Drill bits with AO-type connection

All our drills are:

- Individually **STERILE** packed
- Available in additional sizes and fittings
- Cannulated option also available



Product Code	Supply Chain Code	Product Description	Pk Size
912.060S	Direct Only	1.1mm x 60mm drill bit, quick fit	5
913.085S	FBX85081	1.5mm x 85mm drill bit, quick fit	5
920.100S	Direct Only	1.8mm x 100mm drill bit, quick fit	5
914.100S	FBX85079	2.0mm x 100mm drill bit, quick fit	5
914.125S	FBX85084	2.0mm x 125mm drill bit, quick fit	5
921.100S	Direct Only	2.4mm x 100mm drill bit, quick fit	5
915.110S	FBX85083	2.5mm x 110mm drill bit, quick fit	5
915.125S	FBX85076	2.5mm x 125mm drill bit, quick fit	5
915.170S	FBX85161	2.5mm x 170mm drill bit, quick fit	5
916.100S	FBX85075	2.7mm x 100mm drill bit, quick fit	5
916.125S	FBX85078	2.7mm x 125mm drill bit, quick fit	5
917.125S	FBX85077	3.2mm x 125mm drill bit, quick fit	5
917.145S	FBX85074	3.2mm x 145mm drill bit, quick fit	5
917.195S	Direct Only	3.2mm x 195mm drill bit, quick fit	5
918.110S	FBX85086	3.5mm x 110mm drill bit, quick fit	5
918.145S	FBX85082	3.5mm x 145mm drill bit, quick fit	5
918.195S	Direct Only	3.5mm x 195mm drill bit, quick fit	5
919.145S	FBX85080	4.5mm x 145mm drill bit, quick fit	5
919.195S	FBX85085	4.5mm x 195mm drill bit, quick fit	5

Drill bits with Plain Shank

All our drills are:

- Individually **STERILE** packed
- Available in additional sizes and fittings



Product Code	Supply Chain Code	Product Description	Pk Size
930.045S	Direct Only	1.1mm x 45mm drill bit, plain shank	5
932.085S	FBX21807	2.0mm x 85mm drill bit, plain shank	5
931.070S	Direct Only	1.5mm x 45mm drill bit, plain shank	5
933.095S	FBX21811	2.5mm x 95mm drill bit, plain shank	5
934.100S	Direct Only	2.7mm x 100mm drill bit, plain shank	5
934.085S	Direct Only	2.7mm x 85mm drill bit, plain shank	5
935.106S	Direct Only	3.2mm x 106mm drill bit, plain shank	5
935.130S	Direct Only	3.2mm x 130mm drill bit, plain shank	5
936.180S	Direct Only	3.5mm x 180mm drill bit, plain shank	5
939.120S	FBX21754	4.0mm x 120mm drill bit, plain shank	5
936.112S	Direct Only	3.5mm x 112mm drill bit, plain shank	5
939.180S	Direct Only	4.0mm x 180mm drill bit, plain shank	5
937.126S	Direct Only	4.5mm x 126mm drill bit, plain shank	5
937.180S	Direct Only	4.5mm x 180mm drill bit, plain shank	5
938.150S	Direct Only	6.5mm x 150mm drill bit, plain shank	5

PULSE LAVAGE SYSTEM

Disposable

The Pulse Lavage System has a **BUILT-IN Battery Pack** and is a **FULLY DISPOSABLE, Single-Use** medical device designed for use in Orthopaedic Surgical Procedures, Trauma and General Wound Care.

FLEXIBILITY FOR HIP & KNEE ARTHROPLASTY PROCEDURES IN ONE SYSTEM!

Specifications:

- High speed flow rate up to 1000ml/min
- Low speed flow rate 700ml/min
- Pressure up to 15PSI
- 10.5V Alkaline Battery
- Suction tube, 3 metre length

Orthopaedic Surgery

Pulsed jets of irrigated saline penetrate into the cancellous bone, removing blood and debris, providing an improved bone to cement interface and reducing the risk of fat embolism.

Wound Management

Facilitates quick, efficient removal of necrotic tissue, bacteria and foreign tissue from the wound.



System includes:

- 1 x Hand Piece
- 1 x Suction Tube (3M)
- 1 x Small Cone Tip
- 1 x Coaxial Canal Tip



Coaxial Canal Tip

Product Code	Supply Chain Code	Product Description	Pk Size
L41100	FJT31333	Pulse Lavage System	1
BT04Z	FJT31335	Coaxial Femoral Brush Tip (Sold separately - not illustrated)	1
RT23012Z	FBU86325	Hip Revision Tip (Sold separately - not illustrated)	1



The **SOLOMAX Bone Cement Mixer Systems** enable **bone cement for Orthopaedic Surgery to be mixed in a sealed manner** without human contact or exposure to the outer environment.

This is vital for two main purposes:

- **Decreased Cement Porosity**
- **Prevents Exposure to Harmful Fumes**

Bone Cement Mixer Bowl

Bone Cement Mixer Bowl Kit Includes:

- Mixing bowl
- Filter / Vacuum tubing set
- Spatula
- Curette
- Ampoule breaker



Product Code	Supply Chain Code	Product Description	Pk Size
SMX-5000-9	Direct Only	Bone Cement Mixer Bowl Kit	10
SMX-5000-16	Direct Only	Ampoule Filter (Sold Separately)	10
SMX-5000-6	Direct Only	Vacuum Foot Switch (Sold Separately)	10

Bone Cement Mixer Syringe

Bone Cement Mixer Syringe Standard Kit Includes:

- Syringe cap
- Suction tube with filter
- Mixing stand
- Funnel
- Syringe body
- Mixing rod
- Stopper
- Long nozzle
- Push rod for nozzle

Additional accessories available upon request.



Product Code	Supply Chain Code	Product Description	Pk Size
SMX-5000-12	Direct Only	Bone Cement Mixer Dispensing Syringe Kit	10

HIP REVISION EXTRACTION SET

The **Gomina Extraction Set** allows a hip revision to be done with minimal **bone loss**. This innovative, pioneering extraction set consists of **Single-Use Chisels/Osteotomes** and a **Reusable Handpiece**.

- Dimensionally rigid chisels/osteotomes with **different radii**.
- Precision ground, sharp cutting edges of chisels/osteotomes make a **precise separation of bone & prosthesis** possible.
- **Ergonomically designed Reusable Handpiece** for optimal performance and maximum power transmission.
- **Quick & secure** clamping mechanism.
- **Various fixation positions** for chisels/osteotomes offering adaptability.



Product Code	Supply Chain Code	Product Description	Pk Size
400.001	FBU85969	Handpiece (Non-Sterile, Reusable)	1
401.101S	Direct Only	Osteotome/Straight long (Sterile, Single-Use)	1
401.102S	Direct Only	Osteotome/Straight short (Sterile, Single-Use)	1
401.201S	Direct Only	Osteotome/Curved lateral (Sterile, Single-Use)	1
401.202S	Direct Only	Osteotome/Curved medial (Sterile, Single-Use)	1

The Gomina Extraction Set is available to PURCHASE or HIRE

Contact our Customer Service Team for further information: **01527 559010**

Also Available...

The Gomina Extraction Set for Hip Revision with a **Jig Guided System fitting Corail® (DePuy Synthes) Hip Stems**.

A key feature of this Extraction Set is the various **chisels/osteotomes match to the guide slots and grooves of the jig**.



HEADLIGHTS & LIGHT SOURCES

LLS-9500 Battery Powered LED Headlight System


CUDA® Surgical's LLS-9500 LED headlight produces an **unparalleled bright, white light with superior battery life** so you can work without interruption. The revolutionary **Rear Cranial Support design** will create a new wearing experience which is **lightweight, comfortable, flexible and durable!**



CUDA® Surgical Battery and Charger:

Each battery lasts **2 hours 30 minutes at full intensity, 5 hours at half power (output)**. The charger will fully charge an empty four cell battery in 3 ½ hours.

The **LLS-9500 system comes with two batteries for easy replacement** when one battery becomes low during procedures and is now available with a single bay or a dual bay charger.



LED Light Source – LED4000

The **LED4000** is a state-of-the-art LED Light Source from **CUDA® Surgical** designed to provide illumination for **Headlights & Endoscopes**.

Expected LED Life of 30,000 hours eliminates the need to purchase replacement lamps, thus saving money on operational costs.



- Compact, efficient, durable source of brilliant “daylight” quality light for headlights & endoscopes.
- LED automatically shuts off when cable is removed from the active port.
- Costly periodic lamp changes eliminated.
- Lumen output comparable to many Xenon light sources but more stable and consistent.
- **DISCOUNTS** available for **Bundle Options**:
(Bundle consists of: Light Source, Headlight, Headlight Cable, Floorstand)

Xenon Light Sources also available

ENDOSCOPIC CABLES

CUDA®
SURGICAL

CUDA® Surgical's range of Endoscopic cables are designed to help deliver optimal light for surgical procedures.

High quality, **versatility**, performance & reliability go hand-in-hand with these **Universal** Fibre-Optic cables.



Compatible with:

- KeyMed®
 - Stryker®
 - Eschmann®
 - Wolf®
 - Dyonics®
 - Karl Storz®
 - Olympus® & GU®
- Light Sources**



BLUE

PINK

RED

GREEN

Features & Benefits:

- Bundle Diameter: 5.0mm
- Cable Length: 7.5ft
- Latex Free Silicone Sheathing
- Fully Autoclavable
- Interchangeable Adaptors to suit your requirements
- Available in different colours for ease of identification

**Other colours & size options
available upon request**

FIBRE LIGHT CABLES & ADAPTORS

Interchangeable Cable Adaptors for use with our range of **CUDA® Surgical Universal Endoscopic Cables**.

A Light Source Adaptor and an Instrument Adaptor are required with each cable. Simply choose your required light source and instrument adaptor and screw to the cable ends.

Cable Adaptors (Light Source End)

Karl Storz®
Code C0100



Wolf®/Dyonics®
Code C0101



ACMI®
Code C0102



Olympus®
Code C0103



Fujinon®
Code C0104



Machida®
Code C0105



Pentax®
Code C0106



ACMI®/Olympus®
Code C0107

Not illustrated

Wolf® Extended
Code C0108

Not illustrated



Cable Adaptors (Instrument End)

**Storz®/
Olympus®**
Code C0200



**Wolf®/Dyonics®
Female**
Code C0201



**ACMI®
Female**
Code C0202



**ACMI®
Male**
Code C0204



**ACMI®
Female Clip**
Code C0205



**Wolf®/Dyonics®
External Clip**
Code EX-0203



**Wolf®/Dyonics®
Male**
Code C0203

Not illustrated

**Olympus®
Retractable**
C0206

Not illustrated

Our range of adaptors provide users with the **versatility** of creating a custom Fibre Light Cable, in a **cost effective** manner, by simply changing the adaptors on each end.



REUSABLE INSTRUMENTATION



Chisels, Gouges & Osteotomes

Capener Osteotome

Tufnol Handle, curved, 250mm

N350-206 (6mm)

N350-212 (12mm)

N350-220 (20mm)



Capener Gouge

Tufnol Handle, straight, 250mm

N350-046 (6mm)

N350-055 (15mm)

N350-075 (25mm)



Capener Osteotome

Tufnol Handle, straight, 250mm

N350-106 (6mm)

N350-112 (12mm)

N350-120 (20mm)

N350-130 (30mm)



Capener Chisel

Tufnol Handle, straight, 250mm

N350-006 (6mm)

N350-012 (12mm)

N350-020 (20mm)

N350-030 (30mm)



Capener Gouge

Tufnol Handle, cranked, 235mm

N350-086 (6mm)

N350-088 (8mm)

N350-096 (16mm)





Bone Rongeurs

Lempert

E068-300 (16cm) Straight

N076-275 (19cm)


E068-305 (16cm) Curve

N076-277 (19cm)




Pennybacker

N076-600 (4mm, 20cm) Straight



Northfield


N076-470 (10mm, 20cm) Straight



Olivercrona

N076-582 (5mm, 20cm)


N076-584 (6mm, 23cm)




Ferris Smith Kerrison Rongeur - Detachable



- **Upward Cutting**, standard handle, 10mm jaw opening. With pin ejector (apart from the 1mm bite)




40° UP BITE



90° UP BITE



- **Downward Cutting**, standard handle, 10mm jaw opening. With pin ejector (apart from the 1mm bite)




90° DOWNWARD BITE

SIZE OPTIONS:


- **Length of Shaft:** 18cm (7"), 20cm (8"), 23cm (9"), 25cm (10")
- **Bite:** 1mm, 2mm, 3mm, 4mm, 5mm, 6mm

Ferris Smith Kerrison Rongeur - Bayonet


- **Upward Cutting**, standard handle, 10mm jaw opening.



40° UP BITE



90° UP BITE



SIZE OPTIONS:

- **Length of Shaft:** 13cm
- **Bite:** 2mm, 3mm, 4mm, 5mm


Love Kerrison Rongeur - Upward Cutting



- Upward Cutting,
extra-large handle,
16mm jaw



40° UP BITE



90° UP BITE



SIZE OPTIONS:


- **Length of Shaft:** 18cm (7"), 20cm (8"), 23cm (9"), 25cm (10")
- **Bite:** 2mm, 3mm, 4mm, 5mm, 6mm

Love Kerrison Rongeur - Downward Cutting

- Downward Cutting,
extra-large handle,
14mm jaw opening.



40° DOWN BITE



90° DOWN BITE



SIZE OPTIONS:

- **Length of Shaft:** 18cm (7"), 20cm (8"), 23cm (9"), 25cm (10")
- **Bite:** 2mm, 3mm, 4mm, 5mm, 6mm

For further information, please contact our
Customer Service Team on: **01527 559010**



Wire Cutters


T800-130

(to cut pins up to 5.55mm diameter)




T800-055

Top Cutter, Hard Edge, 21.5cm
(to cut wires up to 2.5mm diameter)



T800-060

Side Cutter, Hard Edge, 23cm
(to cut wires up to 2.2mm diameter)



M129-302

(15cm) straight edge


M129-304

(20cm) straight edge



M129-306

(15cm) curved edge





Rugines & Raspatories

Farabeuf



T120-040

(Straight, Hollow Handle, 16cm)

Farabeuf

T120-050

(Curved, Hollow Handle, 16cm)



Bone Holding Forceps

Farabeuf-Lambotte

T014-550 (26cm)

Langenbeck

T014-320 (21cm)



Bone Levers

Hohmann

T040-075 (22cm/8 3/4in)



Self Retaining Retractors

Travers

G136-580 (4:5 blunt prongs, 20cm)





Mallets

Toffee Mallet

E132-105 (340g)



Heath Mallet


Stainless steel head and handle

G116-118 (510g)

G116-128 (795g)

G116-132 (905g)

G116-136 (1020g)



Partsch Mallet

T050-180 (120g)



Hajek Mallet

E132-375 (140g)



Tufnol Handle Mallet

Stainless steel head

T050-102 (350g)

T050-104 (500g)

T050-106 (700g)



PLASTER ROOM INSTRUMENTS



We offer an **extensive range of Plaster Room Instruments** for use in **Orthopaedic** procedures. Featured are a selection of our most popular products:

HENNIG G094-705

Plaster Spreader,
28cm



ESMARCH G094-405

Plaster Knife,
Solid forged, 18cm



GUYS G094-605

Plaster Shears
Heavy, 29cm



STILLE LORENZ G094-905

Plaster Shears, 23cm
Other sizes available



WIRE CUTTER T800-060

Wire side cutter,
hard edge, T/C edge,
compound action, 23cm



PLIERS T800-045

Orthopaedic Pliers
S/S, 17cm
Cuts up to 1.6mm
soft wires




TOWEL CLIP G090-290 - 10cm G090-295 - 14cm

Ball & Socket
Towel Clip



DAWS G094-505

Plaster Spreader,
45cm



RING CUTTER G018-105 - 16.5cm G018-110

Spare blade



BOHLER SHEARS G094-205

Plaster Shears,
24cm



MICK CAST SCISSORS Female Pattern 48-010-17F



MICK CAST SCISSORS Male Pattern 48-010-17M



ORTHOPAEDIC INSTRUMENTATION



BONE AWLS
BONE CURETTES
BONE HOOKS
BONE LEVERS
BURRS
CHISELS
COMPRESSION CLAMPS
DISLOCATING HOOKS
ELEVATORS
EXTRACTORS
FILES
FORCEPS - VARIOUS
GOUGES
HAND SPLINTS
KNIVES
MALLETS
OSTEOTOMES
PUNCHES
RASPATORIES
RASPS
RETRACTORS
RONGEURS
RUGINES
SCREWDRIVERS
STIRRUPS
TWIST DRILLS
WIRE INSTRUMENTS

Our Comprehensive **Reusable Orthopaedic Instrumentation** Catalogue is available.

To request your copy, contact our Customer Service Team on:




01527 559010



sales@judd-medical.co.uk



www.judd-medical.co.uk



SAW BLADES

Our extensive portfolio includes a wide range of **surgical saw blades** designed for both small & large bone surgeries. These blades feature **hub fittings that are compatible with all leading manufacturers' power tool systems**, tailored specifically for use by Orthopaedic Surgeons.


Our Replacement Surgical Saw Blades are available in a variety of shapes and sizes, fitting to:

- Aesculap®
- Aesculap® Acculan 3Ti
- Arthrex®
- De Soutter®
- Linvatec® / Hall®
- MicroAire®
- Stryker®
- Synthes® (AO)
- Tuke Saw™
- Zimmer®



All our Replacement Surgical Saw Blades are supplied:

- **Sterile packaged**
- **Complete with batch traceability labelling**
- **Sterilised via Gamma Irradiation**
- **Sold in Packs of 5**



Our dedicated Blades Catalogue is now available, contact us:



T: 01527 559010


E: sales@judd-medical.co.uk

W: www.judd-medical.co.uk

Disposable Hand Switch Pencil housing 2.38mm (3/32") with cut-coag double button and Smoke Evacuator System.

This Diathermy Pencil with Smoke Evacuation System is designed to **fulgurate, coagulate, cut or dissect biological tissue and also provide the effective removal of surgical smoke** produced during electrosurgery.

- Available with either **2 standard blade or 2 teflon (non-stick) blade electrodes**
- Generator Connector: Universal tripolar moulded **Valleylab**



- 2 x Cannula with screw connections. (Lengths 37mm & 122mm)
- Pencil Length - 170mm
- Material - High resistance plastic, Latex Free
- Corrugated tubing, highly flexible, 300cm, + -5%

Product Code	Supply Chain Code	Product Description	Pk Size
F4797SE	FDJ2771	Smoke Evacuation Pencil (includes 2 Standard Electrodes)	10
F4797SE/TEF	FDJ2736	Smoke Evacuation Pencil (includes 2 Teflon Electrodes)	10

ELECTROSURGERY

Our product portfolio offers a trusted range of Electrosurgery products for both Surgeons and Clinicians.

Updated Energy & Diathermy Catalogue Now Available!

We are excited to announce the release of our **updated Energy & Diathermy Catalogue**. This latest edition features a **comprehensive range of both disposable and reusable products** designed to meet all your surgical needs.

Catalogue Highlights:

- **Disposable Products:** Explore our extensive collection of disposable diathermy instruments, designed for precision, safety and convenience.
- **Reusable Products:** Discover our durable and reliable reusable diathermy instruments, crafted to withstand rigorous sterilisation processes & repeated use.

How to Access the Catalogue:


You can access the updated catalogue through our website:

www.judd-medical.co.uk

For more information or to request a physical copy of the catalogue, please reach out to us on:

01527 559010

Stay tuned for more updates and product launches as we continue to innovate and provide you with the best surgical solutions.











CATALOGUE RESOURCES & PRODUCT VIDEOS

Our **product catalogues** are available for **download directly from our website**, providing you with instant access to our comprehensive range of products at your convenience. Alternatively, if you prefer a **physical copy** you can easily **request one** and we'll ensure it is promptly **delivered** to you.

Explore our extensive selection of products in the format that suits you best!



SCAN TO VIEW
OUR CATALOGUES



SCAN TO VIEW OUR
VIDEO LIBRARY



If you're looking for some informative content, be sure to check out our **YouTube channel**. Our videos cover a wide range of specialities that provide product features & benefits through to useful tips.

Dive into our videos and join our growing community today!



QUALITY & ACCREDITATIONS

Judd Medical is certified to **ISO 13485: 2016 and ISO 14001: 2015** across our two locations with all relevant policies and procedures being robustly implemented.

These certifications collectively highlight our commitment to excellence and our leadership in both quality assurance and environmental sustainability.



ISO 13485: 2016

Demonstrating our commitment to the safety & quality of our medical devices.



ISO 14001: 2015

Demonstrating our commitment to identifying, managing, monitoring and controlling our environmental processes.

LSAS: Level 2

In addition, we have achieved **Labour Standards Assurance System (LSAS) Level 2.**

NOTES



HOW TO REACH US



01527 559010



sales@judd-medical.co.uk



www.judd-medical.co.uk



**Unity House, Buntsford Park Road,
Bromsgrove, Worcestershire, B60 3DX**



Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk



We are an Approved Supplier to the
NHS Supply Chain



Supply Chain