

Judd
medical

ENT Theatre & Outpatients Catalogue

www.judd-medical.co.uk



NHS

Supply Chain

APPROVED SUPPLIER

TABLE OF CONTENTS

3	Who We Are / How Can We Help
4	What We Do / Training & Education
5	How To Order / Delivery Options
6	Diathermy Pencils (Disposable & Reusable)
7	Diathermy Smoke Evacuation Pencil (Disposable)
8	Bipolar Forceps (Reusable)
9-10	Bipolar Forceps (Disposable)
11-12	Monopolar Forceps (Reusable)
13	Monopolar Forceps (Disposable)
14	Grounding Pads & Cables (Disposable & Reusable)
15	Abbey Needles & Suction Coagulator (Disposable)
16	Monopolar & Bipolar Cables (Reusable)
17-19	Diathermy Electrodes (Disposable & Reusable)
20	Electrocauteries (Disposable)
21	Diathermy Quivers (Disposable)
22	Diathermy Quivers (Reusable)
23-38	Instrumentation (Reusable)
39	Oesophageal Bougies (Disposable)
40	Tie Wires & Snare Wires (Disposable)
41-43	Headlights, Light Sources, Cables & Adaptors
44	Catalogue Resources & Product Videos
45	Quality & Accreditations
46	Notes
47	How To Reach Us

Disclaimer: Product specifications in this catalogue are correct at time of publication but may change from time-to-time without notice. It is the customer's responsibility to validate that a particular product with the properties described in the product specification meets the requirements of their application.

WHO WE ARE



Judd Medical Limited is recognised as being amongst one of the most progressive medical companies in the UK and during the last 25 years have seen numerous developments that have enabled us to grow into a leading medical device supplier in the industry.

With the expansion of our product portfolio to include many innovative products and services, it has allowed Judd Medical to make a significant impact on the healthcare market to NHS hospitals, private groups and independent clinics throughout the UK.

HOW CAN WE HELP



CUSTOMER SERVICE LEVELS

We ensure that should you need any help and advice or encounter any problems we have a fully trained team ready to support you. Our specifically trained Customer Service Team offer friendly, knowledgeable advice and support for every requirement.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any products proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts, standing and call off orders dependent on your requirements.

JUDD MEDICAL

“YOUR PARTNER IN PRECISION HEALTHCARE”

WHAT WE DO.....



Our company is proud to offer a comprehensive product portfolio of
Single Use & Reusable Surgical Instruments and Medical Devices.

-  Arthrodesis & Kirschner Wires, Twist Drills, Steinmann Pins & Denham Pins
-  Electrosurgery: Diathermy Products – Single Use & Reusable
-  Light Sources: Xenon & LED, Surgical Headlights, Fibre Light Cables & Adaptors
-  Pulse Lavage Systems
-  Replacement Surgical Saw Blades & Burs
-  Single-Use & Reusable Surgical Instrumentation

At Judd Medical we are committed to our partnership with the NHS Supply Chain, NHS Scotland and NHS Wales, ensuring optimal value for our NHS Customers.

We are an approved supplier to the NHS Supply Chain.



TRAINING & EDUCATION



We are dedicated to providing comprehensive training and education on regulatory compliance, organising international seminars and hosting healthcare roundtable events. We facilitate peer-to-peer clinical discussions and offer training on the care and handling of surgical instrumentation and diathermy safety.

Through these initiatives, we aim to enhance professional knowledge, ensure adherence to industry standards and promote the highest levels of safety and efficacy in healthcare practices.

HOW TO ORDER:

You can place an order with us via:



01527 559010



sales@judd-medical.co.uk



NHS Framework Agreements



SCAN TO VIEW OUR
NHS FRAMEWORK
AGREEMENTS

DELIVERY OPTIONS:



Shipping & Despatch – Where stock is available this will be despatched on the same day.
(Order must be processed before 3pm)

Delivery options available:

- Standard Next Day Delivery
- Pre 9am Next Day Delivery*
- Pre 10.30am Next Day Delivery*
- Pre 12 Noon Next Day Delivery*
- Same Day Delivery – Must be processed before 12 noon on date of request.*


*Price on Application

We deliver to:

- Mainland UK
- Northern Ireland
- Channel Islands
- Scottish Highlands & Islands



**Centrally located in Bromsgrove,
Worcestershire**



DIATHERMY PENCILS

DISPOSABLE PENCIL with PUSH-BUTTON CONTROLS

Designed to fulgurate, coagulate, cut or dissect biological tissue during electrosurgery.

- Push-button controls for **Cut & Coagulation** functions
- Available with either a **standard** blade or **non-stick** blade electrode
- Individually Sterile Packed

Disposable Diathermy Pencil - Standard

Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4797	FDJ2684	320mm	2.38mm	Valleylab	Blade F4046	25



Disposable Diathermy Pencil - Non-Stick

Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4797TEF	Direct Only	320cm	2.38mm	Valleylab	Blade F4046/TEF	25

REUSABLE PENCIL with PUSH-BUTTON CONTROLS

These Reusable Diathermy Pencils, with **Push-Button Controls for Cut & Coagulation functions**, are designed for optimal performance during electrosurgery and are complete with **standard blade electrodes**.

REUSABLE (can be sterilised up to **100 times**)




Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4244	FDJ2835	300cm	2.38mm	Valleylab	Blade F4060	1

Disposable Hand Switch Pencil housing 2.38mm (3/32") with cut-coag double button and Smoke Evacuator System.

This Diathermy Pencil with Smoke Evacuation System is designed to **fulgurate, coagulate, cut or dissect biological tissue and also provide the effective removal of surgical smoke** produced during electrosurgery.

- Available with either **2 standard blade or 2 teflon (non-stick) blade electrodes**
- Generator Connector: Universal tripolar moulded **Valleylab**



- 2 x Cannula with screw connections. (Lengths 37mm & 122mm)
- Pencil Length - 170mm
- Material - High resistance plastic, Latex Free
- Corrugated tubing, highly flexible, 300cm, + -5%

Product Code	Supply Chain Code	Product Description	Pk Size
F4797SE	FDJ2771	Smoke Evacuation Pencil (includes 2 Standard Electrodes)	10
F4797SE/TEF	FDJ2736	Smoke Evacuation Pencil (includes 2 Teflon Electrodes)	10

BIPOLAR FORCEPS

Reusable

Precision Reusable Bipolar Forceps designed for grasping & coagulation of tissue during electrosurgery.



- Types available: **Adson, Bayonet, Bayonet (angled down), Jeweller#5, McPherson & Straight options**
- Tip size options: **0.5mm or 1.0mm**

50 times Autoclavable



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-103/R1.0	Direct Only	Straight Forcep	12cm	1.0mm	1
PH786-105/R1.0	FDJ2209	Straight Forcep	15cm	1.0mm	1
PH786-109/R1.0	Direct Only	Straight Forcep	18cm	1.0mm	1
PH786-111/R1.0	Direct Only	Straight Forcep	20cm	1.0mm	1
PH786-113/R1.0	Direct Only	Straight Forcep	22cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-127/R1.0	Direct Only	Bayonet Forcep	18cm	1.0mm	1
PH786-130/R1.0	Direct Only	Bayonet Forcep	20cm	1.0mm	1
PH786-131/R1.0	Direct Only	Bayonet Forcep (angled down)	20cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-141/R1.0	Direct Only	Adson Forcep	12cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-151/R0.5	Direct Only	Jeweller Forcep #5 (straight point)	11.5cm	0.5mm	1
PH786-155/R0.5	Direct Only	McPherson Forcep	10.2cm	0.5mm	1

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

BIPOLAR FORCEPS

(with Attached Cable)

Disposable

STANDARD

- Individually sterile packaged & boxed in 25's
- Tip Sizes: **0.5mm & 1.0mm**
- Various Lengths
- Flexible bridge connector with 2 x 4mm plugs
- Suitable for the most popular diathermy machines including:
Valleylab, Conmed, ERBE & more
- Pre-attached 3M cable
- All STERILE bipolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-105/DWC0.5	FDJ2195	Straight Forcep	15cm	0.5mm	STERILE	25
PH786-105/DWC1.0	FDJ2196	Straight Forcep	15cm	1.0mm	STERILE	25
PH786-107/DWC0.5	FDJ2680	Straight Forcep	16cm	0.5mm	STERILE	25
PH786-107/DWC1.0	FDJ2724	Straight Forcep	16cm	1.0mm	STERILE	25
PH786-109/DWC1.0	FDJ2197	Straight Forcep	18cm	1.0mm	STERILE	25
PH786-111/DWC1.0	FDJ2198	Straight Forcep	20cm	1.0mm	STERILE	25
PH786-112/DWC1.0	FDJ2199	Straight Forcep (angled tip)	20cm	1.0mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-124/DWC1.0	FDJ2200	Bayonet Forcep	16cm	1.0mm	STERILE	25
PH786-130/DWC1.0	FDJ2201	Bayonet Forcep	20cm	1.0mm	STERILE	25
PH786-133/DWC1.0	FDJ2202	Bayonet Forcep	22cm	1.0mm	STERILE	25
PH786-136/DWC1.0	FDJ2203	Bayonet Forcep	24cm	1.0mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-141/DWC1.0	FDJ2194	Adson Forcep	12cm	1.0mm	STERILE	25
PH786-151/DWC0.5	FDJ2193	Jeweller Forcep	11.5cm	0.5mm	STERILE	25
PH786-155/DWC0.5	FDJ2192	McPherson Forcep	10.2cm	0.5mm	STERILE	25

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

BIPOLAR FORCEPS

(with Attached Cable)

Disposable

NON-STICK

- Individually sterile packaged & boxed in 25's
- Tip Sizes: **0.5mm & 1.0mm**
- Various Lengths
- Gold Coloured Tips
- Flexible bridge connector with 2 x 4mm plugs
- Suitable for the most popular diathermy machines including:

Valleylab, Conmed, ERBE & more

- Pre-attached 3M cable
- All STERILE bipolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-105/NS0.5	FDJ2733	Non-Stick Straight Bipolar Forcep	15cm	0.5mm	STERILE	25
PH786-105/NS1.0	FDJ2754	Non-Stick Straight Bipolar Forcep	15cm	1mm	STERILE	25
PH786-111/NS1.0	FDJ2700	Non-Stick Straight Bipolar Forcep	20cm	1mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-124/NS1.0	FDJ2755	Non-Stick Bayonet Bipolar Forcep	16cm	1mm	STERILE	25
PH786-130/NS1.0	FDJ2653	Non-Stick Bayonet Bipolar Forcep	20cm	1mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-141/NS1.0	FDJ2751	Non-Stick Adson Bipolar Forcep	12cm	1mm	STERILE	25
PH786-151/NS0.5	FDJ2648	Non-Stick Jeweller Bipolar Forcep	11.5cm	0.5mm	STERILE	25



Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

MONOPOLAR FORCEPS

Reusable

Reusable Monopolar Forceps for grasping & coagulation of tissue during electrosurgical diathermy procedures.



- Types available: **Adson, Gillie (Serrated), McIndoe, Waugh, Gillies 1:2 Toothed, DeBakey & Straight options**
- Tip Size: **2.0mm**
- Various Lengths available



McIndoe Forcep



Straight Forcep



Adson Forcep



Gillie (Serrated) Forcep

50 times Autoclavable

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-301/R2.0	FDJ2666	Straight Standard Forcep	12cm	2.0mm	1
PH786-303/R2.0	FDJ2738	Straight Standard Forcep	15cm	2.0mm	1
PH786-305/R2.0	FDJ2732	Straight Standard Forcep	18cm	2.0mm	1
PH786-307/R2.0	FDJ2670	Straight Standard Forcep	20cm	2.0mm	1
PH786-309/R2.0	FDJ2656	Straight Standard Forcep	22cm	2.0mm	1
PH786-311/R2.0	FDJ2725	Straight Standard Forcep	25cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-313/R2.0	FDJ2645	McIndoe Forcep	15cm	2.0mm	1
PH786-315/R2.0	FDJ2762	Gillie (Serrated) Forcep	15cm	2.0mm	1
PH786-317/R2.0	FDJ2675	Adson Forcep	12cm	2.0mm	1
PH786-319/R2.0	FDJ2187	Adson Forcep	15cm	2.0mm	1


Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

MONOPOLAR FORCEPS



Reusable

- Tip Size: **2.0mm**
- Various Lengths available



Waugh Forcep



De Bakey Forcep



Gillies 1:2 Forcep

50 times Autoclavable

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-321/R2.0	FDJ2756	Waugh Forcep	15cm	2.0mm	1
PH786-323/R2.0	FDJ2763	Waugh Forcep	18cm	2.0mm	1
PH786-325/R2.0	FDJ2721	Waugh Forcep	20cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-335/R2.0	FDJ2635	Gillies 1:2 Toothed Forcep	15cm	2.0mm	1
PH786-337/R2.0	FDJ2765	Gillies 1:2 Toothed Forcep	18cm	2.0mm	1
PH786-339/R2.0	FDJ2746	Gillies 1:2 Toothed Forcep	20cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-345/R2.0	FDJ2695	De Bakey Forcep	15cm	2.0mm	1
PH786-347/R2.0	FDJ2848	De Bakey Forcep	18cm	2.0mm	1
PH786-349/R2.0	FDJ2649	De Bakey Forcep	20cm	2.0mm	1

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.


MONOPOLAR FORCEPS

(with Attached Cable)

Disposable

- Suitable for the most popular diathermy machines including:
 - **Valleylab**
 - **Conmed**
 - **ERBE & more**
- Pre-attached 3M cable with **8mm generator pin**
- All STERILE monopolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue

Straight
Forcep



McIndoe
Forcep



1:2 Teeth
Gillies Forcep



Waugh
Forcep



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-305/DWC2.0	FDJ2191	Serrated Straight Dissecting Monopolar Forcep	18cm	2.0mm	STERILE	25
PH786-313/DWC2.0	FDJ2189	Serrated McIndoe Monopolar Forcep	15cm	2.0mm	STERILE	25
PH786-339/DWC2.0	FDJ2184	1:2 Teeth Gillies Monopolar Forcep	20cm	2.0mm	STERILE	25
PH786-523/DWC2.0	FDJ2190	Serrated Waugh Monopolar Forcep	18cm	2.0mm	STERILE	25



Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

GROUNDING PADS & CABLES

Electrosurgical Grounding Pads & Cables – used during electrosurgery to protect the patient from burn injuries & harmful effects of the electric current.

Disposable & Reusable options available.

Single Use, Disposable Grounding Pads available in **Adult, Paediatric & Neonatal** sizes.

Both disposable & reusable grounding pads **offer secure & safer surgery**, however, a proper evaluation of both equipment and patient should be undertaken when deciding which grounding pads to use.

Connection Types Available:

- **Tab Connection**
- **Valleylab Connector**
- **Ø6.3mm Bipolar Connector**
- **4mm Socket**

- Adult, Paediatric & Neonatal Sizes
- Foam, Non-Woven, Silicone Backings
- Connection Types: Tab Connection, Valleylab Connector, Ø6.3mm Bipolar Connector, 4mm Socket





Also available: **Connection Cables and Reusable Adaptors.**

Product Code	Product Description	Cable Length	Machine Fitting
F7922	Universal 6.3mm Mono Jack	500cm	
F7922/WB	Universal 6.3mm Mono Jack (without Derivation Box)	500cm	
F7922/24	Double 4mm Plug	500cm	
F7923	Valleylab Type (ESU with REM)	500cm	
F7924	6.3mm Mono Jack (ESU without REM)	500cm	

ABBAY NEEDLES

Disposable



Product Code	Supply Chain Code	Product Description	Size	Pk Size
PH786-389/DWC	FDJ2786	Abbay Needle	17.5cm	25

SUCTION COAGULATOR

Disposable

Single Use, Cannula Suction Coagulator ideally suited for ENT surgery provides **increased efficiency** during procedures requiring **controlled fluid suction**.

Features:

- Suction & Coagulation can be performed **independently or simultaneously**
- **Vented handle** for suction control
- **Malleable cannula** allows surgeons to match instruments to specific surgical need
- Provide suction corrugated tubing & cleaning stylet
- Sold in **Packs of 10**

Specifications:

- Control Type: **Button Control**
- Cannula: 3.8mm diameter, 150mm length
- Suction Corrugated Tubing Length: 300mm
- Connection Cable & Plug: 3m with 3 pin plug

No greater than standard coagulation power of 30 watts is recommended for ENT operation.



Product Code	Product Description	Cannula Length	Tubing Length	Connector	Pk Size
OBS0020S	Suction Coagulator	150mm	300mm	3 pin plug	10

MONOPOLAR CABLES




Reusable

These Monopolar Cables are compatible with most standard electrosurgical devices and compatible with the following generators and forceps:

Valleylab, Conmed, Megadyne, Bovie, Dremed and all others that are similar to the aforementioned brands.

50 times Autoclavable

Product Code	Product Description	No. of Uses	Machine Fitting	Instrument Fitting	Pk Size
PHF685-3M	Monopolar Cable 4mm pin to 4mm	50			1
PHF690-3M	Monopolar Cable 4mm pin to 4.8mm	50			1
PHF675-3M	Monopolar Cable 8mm pin to 4mm	50			1
PHF680-3M	Monopolar Cable 8mm pin to 4.8mm	50			1

BIPOLAR CABLES



Reusable

Product Code	Product Description	No. of Uses	Machine Fitting	Instrument Fitting	Pk Size
OBS581F	Reusable Bipolar Cable Blue with 4mm Pins & Flexible Bridge Connector European Block End Connector 3M Length	50			1
OBS580TF	Reusable Bipolar Cable Blue with 4mm Pins & Flexible Bridge Connector US Two Pin Forceps Connector 3M Length	50			1

Universal Bipolar Cables available upon request.

ELECTRODES

Disposable

The **Teflon (PTFE) coating** makes these electrodes “Non-Stick” and is only present on the active part of the electrode. Suitable for general diathermy.

STAINLESS STEEL - Sterile

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044	FDJ2652	70mm	BALL 4mm	24
F4046	FDJ2715	70mm	BLADE	10
F4046/S	Direct Only	50mm (44mm sleeve)	BLADE	24
F4048	FDJ2661	70mm	NEEDLE	10
F4048/S	Direct Only	69mm (63mm sleeve)	NEEDLE	24
F4049	FDJ2636	54mm	BALL 2mm	24
F4050	Direct Only	153mm	BLADE	24
F4050/S	Direct Only	153mm (147mm sleeve)	BLADE	24
F4052	Direct Only	152mm	NEEDLE	24
F4052/S	Direct Only	152mm (146mm sleeve)	NEEDLE	24
F4053	Direct Only	100mm	BLADE	24
F4067	FDJ2787	156mm	BALL 4mm	24
F4067/S	FDJ2792	157mm	BALL 5mm	24



TEFLON COATED - Sterile - NON STICK

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044/TEF	FDJ2836	70mm	BALL 4mm	24
F4046/TEF	FDJ2647	70mm	BLADE	10
F4046S/TEF	FDJ2658	70mm (65mm sleeve)	BLADE	10
F4048/TEF	FDJ2663	70mm	NEEDLE	10
F4048S/TEF	FDJ2729	70mm (65mm sleeve)	NEEDLE	10
F4050/TEF	FDJ2757	150mm	BLADE	10
F4050S/TEF	FDJ2639	150mm (145mm sleeve)	BLADE	10
F4052/TEF	Direct Only	150mm	NEEDLE	24
F4053/TEF	FDJ2846	100mm	BLADE	10
F4053S/TEF	FDJ2791	100mm (95mm sleeve)	BLADE	10
F4067/TEF	FDJ2828	150mm	BALL 4mm	10
OBS0312C64	Direct Only	64mm (shielded tip)	BLADE	10
OBS0324C120	Direct Only	120mm	BALL 3mm	10
OBS0323C120	Direct Only	120mm	BALL 5mm	10
OBS0314C130	Direct Only	130mm	BLADE	10




These devices are Latex Free.

Images shown do not reflect the size or entire range of Electrodes available.

The **Titanium coating** makes these electrodes “Non-Stick” and is only present on the active part of the electrode.

TITANIUM COATED - Sterile - NON STICK

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044T	FDJ2651	70mm	BALL 4mm	10
F4046T	Direct Only	70mm	BLADE	24
F4048T	Direct Only	70mm	NEEDLE	24
F4050T	FDJ2671	153mm	BLADE	24
F4052T	Direct Only	152mm	NEEDLE	24
F4067T	FDJ2840	156mm	BALL 4mm	24




Disposable **Micro Dissection Needles** suitable for a wide variety of diathermy procedures and allowing for greater accuracy. The **TUNGSTEN Dissection Needles** offer **precise tissue dissection** with improved cutting & coagulation performance.

MICRO DISSECTION NEEDLES - Sterile



Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4051/TU108	FDJ2837	108mm (104mm sleeve)	NEEDLE	10
F4051/TU75	Direct Only	75mm (71mm sleeve)	NEEDLE	24
F4048/TU	Direct Only	65mm (61mm sleeve)	NEEDLE	24
F4051/TU	FDJ2808	55mm (51mm sleeve)	NEEDLE	10
F4051/45TU	FDJ2854	55mm (51mm sleeve)	NEEDLE 45° Angled	10
F4042/TU	FDJ2676	55mm (51mm sleeve)	NEEDLE 30° Angled	10
F4047/TU	Direct Only	45mm (41mm sleeve)	NEEDLE	24



Images shown do not reflect the size or entire range of Electrodes available.

ELECTRODES

Reusable

STAINLESS STEEL 2.38MM SHAFT

Product Code	Supply Chain Code	Length	Shape	Pk Size
F4045	FDJ2816	156mm	SPHERE 4mm	5
F4045/10	FDJ2655	156mm	SPHERE 4mm, 10° angled	5
F4045/20	FDJ2735	156mm	SPHERE 4mm, 20° angled	5
F4045/30	FDJ2654	156mm	SPHERE 4mm, 30° angled	5
F4045/40	FDJ2728	156mm	SPHERE 4mm, 40° angled	5
F4053R	FDJ2822	100mm	BLADE	5
F4060	FDJ2826	69mm	BLADE	5
F4061	FDJ2820	100mm	ROUND NEEDLE, 20° angled	5
F4062	FDJ2827	150mm	CURVED BLADE	5
F4063	FDJ2779	152mm	NEEDLE 0.78mm	5
F4064	FDJ2857	150mm	SPHERE 2.5mm	5
F4065	FDJ2782	153mm	BLADE	5
F4066	FDJ2781	152mm	NEEDLE 1.30mm	5
F4068	FDJ2780	70mm	SPHERE 4mm	5
F4084	FDJ2803	69mm	NEEDLE 0.75mm	5
F4086	FDJ2818	69mm	NEEDLE 1.30mm	5



STAINLESS STEEL 4MM SHAFT

Product Code	Supply Chain Code	Length	Shape	Pk Size
F4064/4	Direct Only	55mm (51mm sleeve)	BLADE	5
F4065/4	FDJ2793	160mm	BLADE	5
F4068/4	FDJ2847	50mm	SPHERE 4mm	5
F4084/4	FDJ2807	55mm	NEEDLE 0.8mm	5

EXTENSIONS

Product Code	Length	Pk Size
EXTENSIONS DISPOSABLE		
F4902	102mm	5
F4903	150mm	5
F4904	182mm	5
EXTENSIONS REUSABLE		
F4952	102mm	5
F4953	150mm	5
F4954	182mm	5



FIAB Sterile Cautery pens are used in haemostasis, granting a **fast and precise cauterisation of small blood vessels**. They are extremely **easy to handle**, thanks to an **internal long life battery**; several models are available, different from each other for **working temperatures (1200c, 800c and 600c)**, **shape, thickness and length** of the cautery tip.

The available range allows the user to select a cautery pen that perfectly matches the clinical demands of the procedure.

Supplied in Packs of 10



Product Code	Supply Chain Code	Tip Shape	Tip Length (mm)	Temp at activation (°C/°F)	Power Supply (V)	Diameter (mm)	Weight (GR)	Series
F7244	FDJ2694	FINE	28	1200/2200	3,0	16	65	HTC*
F7266	FDJ2770	LARGE	28	1200/2200	3,0	16	65	HTC*
F7234	FDJ2679	FINE	75	1200/2200	3,0	16	78	HTC*
F7277	FDJ2764	FINE	125	1200/2200	3,0	16	78	HTC*
F7288	FDJ2683	THICK	28	800/1400	3,0	16	65	MTC*
F7255	FDJ2711	FINE	28	600/1100	1,5	16	45	LTC*
F7255UF	FDJ2772	ULTRA FINE	28	600/1100	1,5	16	45	LTC*

HTC* High Temperature Cautery
MTC* Med. Temperature Cautery
LTC* Low Temperature Cautery



DIATHERMY QUIVER with Laparoscopic Extension

This Disposable Diathermy Quiver is a **cost effective** means of offering protection to patients and medical staff by **safely securing instrumentation during surgery**.

An additional benefit of this product is the availability of a **Laparoscopic Extension** which is designed to be used in conjunction with the aforementioned quiver. This open ended extension fits into the quiver providing the end user with an **extended overall working length of 400mm**.



Features & Benefits:

- Single Use, STERILE Packaged
- Lightweight, Strong & Durable
- Integral Towel Clip ring
- Clear in colour
- Double Wrapped
- COMPETITIVELY PRICED
- Ease of accessibility



Product Code	Supply Chain Code	Product Description	Pk Size
AVM4047	Direct Only	Single Use Diathermy Quiver - 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	65
AVM4085	Direct Only	Single Use Laparoscopic Extension - 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	40
AVM4086	Direct Only	Laparoscopic Quiver Combined - 20cm + 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	40

With **3 different options available**, these **Reusable Diathermy Quivers** provide a secure and hygienic method of storing diathermy instruments while not in use during surgical procedures.

Options Available:

- **Standard Diathermy Quiver**
- **Diathermy Quiver with Removable Screw Bottom**
- **Julian Taylor Diathermy Quiver with Drainage Holes**



Standard Diathermy Quiver



Diathermy Quiver with Removable Screw Bottom



Julian Taylor Diathermy Quiver with Drainage Holes



Features & Benefits:

- REUSABLE & SUSTAINABLE
- Robust & Long Lasting
- Easy to Clean
- Available in different sizes
- UK Manufactured
- Material: Polypropylene
- DISCOUNTS available for quantity purchases
- COMPETITIVELY PRICED, excellent value for money

Product Code	Product Description	Size	Pk Size
G022-987	Diathermy Quiver with Removable Screw Bottom	D 6cm x L 20cm	1
G022-975		D 6cm x L 24cm	1
G022-985		D 6cm x L 37.5cm	1
G022-905	Standard Diathermy Quiver	D 4.8cm x L 24.5cm	1
G022-910		D 6cm x L 23cm	1
G022-950		D 6cm x L 37.5cm	1
G022-915	Julian Taylor Diathermy Quiver with Drainage Holes	D 4.4cm x L 25cm	1

REUSABLE INSTRUMENTATION



Aural Speculum


GRUBER AURAL SPECULUM

ROUND END

- E034-013** - No 1, 3mm
- E034-014** - No 2, 4mm
- E034-015** - No 3, 5mm
- E034-017** - No 4, 7mm
- E034-010** - Nos 1-4, Set of 4
- E034-011** - Nos 2-4, Set of 3

OVAL END

- E034-023** - No 1, 2.5mm x 3.5mm
- E034-024** - No 2, 3.5mm x 4.5mm
- E034-025** - No 3, 4mm x 6.0mm
- E034-026** - No 4, 6mm x 8.0mm
- E034-020** - Nos 1-4, Set of 4
- E034-021** - Nos 2-4, Set of 3



HARTMAN AURAL SPECULUM

Round End, Stainless Steel

- E034-612** - No 1, 2mm
- E034-613** - No 2, 3mm
- E034-614** - No 3, 4mm
- E034-615** - No 4, 5.5mm
- E034-616** - No 5, 6.5mm
- E034-617** - No 6, 7.5mm
- E034-600** - Set of 3, Nos 2-4
- E034-602** - Set of 3, Nos 3-5
- E034-604** - Set of 4, Nos 3-6

Black Finish

- E034-633** - No 1, 3mm
- E034-634** - No 2, 4mm
- E034-635** - No 3, 5mm
- E034-637** - No 4, 7mm
- E034-639** - No 5, 7.5mm
- E034-640** - No 6, 9mm
- E034-630** - Set of 4, Nos 1-4
- E034-631** - Set of 6, Nos 1-6



ROSEN AURAL SPECULUM

Black Finish, Unslotted

Round End

- E034-073** - No 1, 3mm
- E034-074** - No 2, 4mm
- E034-075** - No 3, 5mm
- E034-076** - No 4, 6mm
- E034-077** - No 5, 7mm
- E034-071** - Set of 4, Nos 2-5
- E034-070** - Set of 5, Nos 1-5

Oval End


- E034-092** - 2.5mm x 3.5mm, No 1 Oval
- E034-094** - 3.5mm x 4.5mm, No 2 Oval
- E034-095** - 4mm x 6mm, No 3 Oval
- E034-096** - 5mm x 7mm, No 4 Oval
- E034-097** - 6mm x 8mm, No 5 Oval
- E034-090** - Set of 4, Nos 2-5



TUMARKIN TRANSMEATAL AURAL SPECULUM

Black Finish, Round, Slotted

- E034-345** - No 0, 3.5mm
- E034-355** - No 0.5, 4.5mm
- E034-365** - No 1, 5.5mm
- E034-370** - No 2, 6mm
- E034-375** - No 3, 6.5mm
- E034-380** - No 4, 7mm
- E034-385** - No 5, 7.5mm
- E034-390** - No 6, 8mm
- E034-391** - No 7, 10mm
- E034-310** - Set of 9, Nos 0-7



SHEA AURAL SPECULUM

Oblique End, Black, 38mm

- E034-245** - No 1
- E034-255** - No 2
- E034-275** - No 3
- E034-285** - No 4
- E034-295** - No 5
- E034-200** - Set of 4, Nos 2-5
- E034-205** - Set of 5, Nos 1-5




REUSABLE INSTRUMENTATION



Nasal Speculum / Nasal Forceps

ST CLAIR THOMSON NASAL SPECULUM

- E182-438** - 38mm Blade
- E182-451** - 50mm Blade
- E182-464** - 65mm Blade
- E182-476** - 76mm Blade
- E182-400** - Set of 4



THUDICUM NASAL SPECULUM

- E182-501** - No 1, 12mm Blade
- E182-502** - No 2, 14mm Blade
- E182-503** - No 3, 16mm Blade
- E182-504** - No 4, 18mm Blade
- E182-505** - No 5, 20mm Blade
- E182-506** - No 6, 22mm Blade
- E182-507** - No 7, 24mm Blade



KILLIAN NASAL SPECULUM


- E182-150** - 5cm Blade
- E182-165** - 6.5cm Blade
- E182-175** - 7.5cm Blade
- E182-190** - 9cm Blade
- E182-100** - Set of 4



COTTLE NASAL SPECULUM


Satin Finish

- E182-010** - No 1, 20mm Blade
- E182-015** - No 2, 32mm Blade
- E182-020** - No 3, 41mm Blade
- E182-025** - No 4, 60mm Blade
- E182-050** - Set of 4



LEMPERT ENDAURAL SPECULUM

- E182-293** - Small, 13.5cm
- E182-295** - Medium, 13.5cm
- E182-299** - Large, 13.5cm



SHAH NASAL SPECULUM

- E182-305** - Left
- E182-310** - Right
- E182-315** - Insulated, Left
- E182-320** - Insulated, Right



LUC NASAL TURBINATE FORCEP

Oval Shaped Jaw, 19.5cm

- E066-405** - Small
- E066-410** - Medium
- E066-415** - Large



LUC NASAL TURBINATE FORCEP


Heart Shaped Jaw, 19.5cm

- E066-305** - Small
- E066-310** - Medium
- E066-315** - Large



GRUNWALD NASAL TURBINATE PUNCH FORCEP 11.5cm

- E066-105** - Pattern A
- E066-110** - Pattern B
- E066-115** - Pattern C



IRWIN MOORE NASAL TURBINATE FORCEP

- E066-205** - 19.5cm




REUSABLE INSTRUMENTATION



Aural Forceps

HARTMAN AURAL FORCEP


Heavy Pattern,
Screw Joint



E038-105 - 14.5cm

HENCKEL (STRUMPEL) AURAL FORCEP


Crocodile Action,
Fine Fenestrated Blade



E038-205 - Effective Length 7.5cm

TILLEY AURAL FORCEP


Serrated Jaw,
Box Joint



E038-405 - Working Length 6.7cm
Overall Length 14.5cm

PRITCHARD AURAL FORCEP

Slender Jaw,
Square Serrated Ends



E038-305 - Effective Length 5cm
Overall Length 14.5cm

CUMBERBATCH AURAL FORCEP

Spring Pattern,
Serrated



E040-105 - 14cm

TONSIL HOLDING FORCEP

Denis Browne,
Box Joint, Angled



E082-105 - 20cm

FAGGE AURAL FORCEP

Bayonet,
Serrated Jaw



E040-205 - 14cm

HUNTER TODD AURAL FORCEP

Bulbous,
Serrated Jaw



E040-305 - 12.5cm



HARTMAN AURAL FORCEP

Crocodile Action, Satin Finish

E036-415 - Standard Serrated Jaw, 7mm x 1.25mm

E036-420 - Extra Fine Serrated Jaw, 4mm x 0.8mm

E036-410 - Round Cup Jaw, 2mm Diameter



REUSABLE INSTRUMENTATION




Aural Scissors / Micro Aural Scissors

CAWTHORNE AURAL SCISSORS

Extra Fine,
Satin Finish,
4mm x 1.5mm


- E044-105** - Straight
- E044-110** - Curved Left
- E044-115** - Curved Right



CAWTHORNE MICRO AURAL SCISSORS

Satin Finish,
3mm x 0.8mm

- E048-105** - Straight
- E048-110** - Curved Left
- E048-115** - Curved Right



ORMEROD AURAL SCISSORS

Extra Fine,
Satin Finish,
4mm x 1.5mm


- E044-205** - Curved Upwards



ORMEROD MICRO AURAL FORCEP

Crocodile Action,
Circular Cup,
0.5mm Dia. Cup

- E046-105** - Straight Jaw



HOUGH MICRO AURAL FORCEP

Serrated Jaw,
Satin Finish,
4mm x 0.5mm


- E046-310** - Straight
- E046-315** - Curved Left
- E046-320** - Curved Right



HOUGH MICRO AURAL FORCEP

Oval Cup Jaw,
Satin Finish,
4mm x 0.5mm


- E046-305** - Straight



BELLUCCI (SHEA) AURAL SCISSORS

Crocodile Action,
Satin Finish,
8mm x 1.5mm


- E044-305** - Curved Left
- E044-310** - Curved Right



CHEVALIER JACKSON FOREIGN BODY FORCEP

Fenestrated Jaw,
Lightweight

- E059-030** - 30cm
- E059-035** - 35cm
- E059-040** - 40cm
- E059-045** - 45cm
- E059-050** - 50cm



REUSABLE INSTRUMENTATION



Aural Syringe / Percussors

AURAL SYRINGE

E184-105 - 2oz capacity, chrome plated, with 2 pipes

E184-110 - 3oz capacity, chrome plated, with 2 pipes

E184-115 - 4oz capacity, chrome plated, with 2 pipes

E184-205 - Bulbous metal pipe, with wings

E184-305 - Conical metal pipe, with wings



NATIONAL HOSPITAL PERCUSSOR (ROSSIER)

Nylon Handle,
5cm Rubber Head



G162-130 - 35cm

TAYLOR PERCUSSOR

G162-135 - 19.5cm



Other Percussors available:

BABINSKI / BUCKS / DEJERINE/ HURST



Snares / Curettes

GLEGG NASAL SNARE

E178-205 - Open Ended, 16.5cm

E178-210 - Self Forming, 16.5cm

E178-215 - Open Ended, 23cm


E178-220 - Self Forming, 23cm



BALLANCE AURAL SNARE

Fine

E178-105 - 15cm



EVE TONSIL SNARE

Slide Action,
Stainless Steel



E180-105 - 28cm

BECKMANN ADENOID CURETTE

21.5cm



Sizes ranging from 10mm - 17mm

***Contact our Customer Service Team for Product Codes**

Other Curettes available:

BARNHILL / BUCK / WILLIGER / ANTRUM / VOLKMANN / HOUSE / MACKIE / ST CLAIR THOMSON

REUSABLE INSTRUMENTATION




Forceps / Needle Holders

TILLEY NASAL DRESSING FORCEP

Serrated Jaw

E054-105 - Working Length 8.2cm,
Overall Length 16.5cm



TROUSSEAU TRACHEAL DILATING FORCEP

E056-205 - Child 10.5cm
E056-210 - Adult 14cm



MAYO NEEDLE HOLDER

Serrated Jaw,
Box Joint


G124-140 - 12.5cm

G124-142 - 14cm **G124-152** - 20cm

G124-144 - 15cm **G124-154** - 23cm

G124-146 - 16cm **G124-156** - 25cm

G124-148 - 18cm **G124-158** - 30cm



KILNER NEEDLE HOLDER

Serrated Jaw,
Box Joint

G124-112 - 13.5cm

G124-114 - 16.5cm



WATSON WILLIAMS NASAL POLYPUS FORCEP

E062-205 - Serrated Jaw,
Effective Length 12cm

E062-207 - Plain, 6mm

E062-209 - Plain, 8mm

E062-211 - Plain, 10mm



KALUSKAR FISHBONE FORCEP

Satin Finish

E061-750 - 12.5cm



LISTER SINUS FORCEP

Serrated Jaw,
Straight

E078-105 - 12.5cm

E078-110 - 15cm

E078-115 - 18cm

E078-120 - 20cm



CRILE ARTERY FORCEP

Serrated Jaw,
Box Joint

STRAIGHT

G052-040 - 14cm

G052-042 - 16cm


G052-044 - 18cm

CURVED

G052-041 - 14cm

G052-043 - 16cm

G052-045 - 18cm




MAGILL INTRODUCING FORCEP

Endotracheal,
Screw Joint

A090-050 - Adult, 25cm

A090-100 - Child, 20cm

A090-150 - Infant, 15cm



LEMPERT NIBBLING FORCEP

Straight,
2.5mm Bite

E068-300 - 16cm



REUSABLE INSTRUMENTATION



Forceps / Scissors

GILLIES DISSECTING FORCEP

15cm



G062-266 - Serrated

G062-268 - 1:2 Teeth

TREVES DISSECTING FORCEP

1:2 Teeth

G062-484 - 12.5cm

G062-486 - 15cm

G062-488 - 18cm

G062-490 - 20cm

G062-492 - 25cm



MCINDOE DISSECTING FORCEP

15cm



G062-340 - Serrated

G062-345 - 1:2 Teeth

UNIVERSAL WIRE CUTTING SCISSORS

Angled Serrated Blade



E172-805 - 12cm

COTTLE DRESSING SCISSORS

Curved



E172-205 - 14cm

FOMAN SCISSORS

E172-305 - BL/BL,

Angled, 14cm

E172-310 - Upper Lateral,

Curved, 14cm

E172-315 - TC Blade,

BL/BL, Angled, 14.5cm



KNAPP SCISSORS

Satin Finish,
Blunt Points,
11.5cm



E172-505 - Curved

E172-510 - Straight

MULTI-CUT UTILITY SCISSORS

S/S Serrated Blade,
Plastic Handles,
Autoclavable,
Pks of 10



JMR/MC2 - 18cm

Colours Available:

Yellow, Blue,
Green & Red



REUSABLE INSTRUMENTATION



Middle Ear Instruments

CAWTHORNE MIDDLE EAR HOOK

Satin Finish, 16.5cm



E108-205 - No 1, Slight Curve

E108-210 - No 2, Medium Curve

E108-215 - No 3, Full Curve

E108-220 - No 4, Double Curve

E108-225 - No 5, Angled 90°, 0.5mm

E108-230 - No 6, Angled 90°, 1mm

E108-235 - No 7, Retro Curve

E108-240 - No 8, Angled 45°, 0.5mm

E108-245 - No 9, Straight

E108-200 - Set of 9

SCHUKNECHT MIDDLE EAR HOOK

Satin Finish, 16cm



E108-805 - Elevator, Oval Window

E108-810 - Hook, Downward, 0.4mm

E108-815 - Hook, Upward, 0.4mm

E108-820 - Hook, Downward, 0.7mm

E108-825 - Hook, Upward, 0.7mm

E108-830 - Hook, Angled 90°, Upward

E108-835 - Pick, Straight

E108-840 - Spatula

GARDINER BROWN TUNING FORKS



E083-105 - Pitch C0, 128dv

E083-115 - Pitch C2, 512dv

E083-125 - Pitch C4, 2048dv

E083-135 - Set of 6, Complete with Box

E083-110 - Pitch C1, 256dv

E083-120 - Pitch C3, 1024dv

E083-130 - Pitch C5, 4096dv

ROSEN MIDDLE EAR HOOK

Satin Finish



E108-705 - No 1, Knife, Meatal Incision

E108-710 - No 2, Elevator, Broad

E108-715 - No 3, Elevator, Narrow

E108-720 - No 4, Elevator, Drum Lifting

E108-725 - No 5, Curette, Tympanic Ring, Small

E108-730 - No 6, Curette, Tympanic Ring, Large

E108-735 - No 7, Probe, Fine

E108-740 - No 8, Mobilizer, Slight Curve

E108-745 - No 9, Mobilizer, Full Curve

E108-760 - No 11, Scoop, Oval, Size 000

E108-765 - No 11, Scoop, Oval, Size 00

E108-770 - No 11, Scoop, Oval, Size 0

E108-775 - No 11, Scoop, Oval, Size 1

E108-780 - Fenestrator, Oval Window

REUSABLE INSTRUMENTATION



Wax Hooks

ST BARTHOLEMEWS WAX HOOK

Blunt

E096-505 - 17cm



HARTMAN WAX HOOK

Sharp

E096-405 - 15cm



FORMBY WAX HOOK


Hook & Scoop, Double Ended

E096-205 - 18cm



CAWTHORNE WAX HOOK

E096-105 - 14.5cm



GARDINER BROWN WAX HOOK

Hook & Imray Loop,
Double Ended

E096-305 - 12.5cm



REUSABLE INSTRUMENTATION



Probes

JOBSON HORNE PROBE

Double Ended

E150-305 - 14cm

E150-310 - 18cm



HEATH PROBE

Nasal


E150-205 - 18cm



HEATH PROBE

Wool Carrier

E150-210 - 17cm



LOGAN TURNER PROBE

E150-505 - 13.5cm



KINGS COLLEGE HOSPITAL

Cocaine Wire

E150-405 - 9cm


E150-410 - 13.5cm



TUMARKIN

Twisted Wire, Pk Size 10

E150-705 - 15cm



MUSGROVE PROBE

Double Ended

E150-605 - 18cm



PROBES WITH EYE

Single Ended, Silver or Nickel Silver

Sizes ranging from 12.5cm - 30cm

*Contact our Customer Service Team for Product Codes



REUSABLE INSTRUMENTATION



Knives / Gouges

LEMPERT SICKLE KNIFE

Meatal Flap Incision

E126-105 - 16cm



LUMSDEN KNIFE

First Incision

E120-305 - 16.5cm



COTTLE KNIFE

Single Ended,
Straight

E118-135 - 13.5cm



BALLENGER NASAL CARTILAGE KNIFE

Straight Stem, Revolving Blade

E118-105 - 3mm Wide



JENKINS GOUGE

14.5cm

Sizes ranging from 2mm - 14mm

*Contact our Customer Service Team for Product Codes




COUPLAND GOUGE

P070-003 - 3mm

P070-004 - 4mm

P070-005 - 5mm



KILLIAN NASAL GOUGE

Bayonet

E092-205 - 16cm



TILLEY NASAL GOUGE

Bayonet Shape

E092-305 - 16cm



REUSABLE INSTRUMENTATION



Mouth Gags

BOYLE DAVIS MOUTH GAG

15mm Wide Fitting

- E084-105** - Child, 13.5cm x 7cm
- E084-110** - Adult, 15cm x 9cm
- E084-150** - Replaceable Tooth Guard, Silicone Rubber, In Pairs



BOYLE DAVIS TONGUE PLATE

Plain

- E084-205** - 6.5cm
- E084-210** - 7.5cm
- E084-215** - 9cm
- E084-220** - 10cm
- E084-225** - 11.5cm



BOYLE DAVIS TONGUE PLATE

With Anaesthetic Tube

- E084-305** - 6.5cm
- E084-310** - 7.5cm
- E084-315** - 9cm
- E084-320** - 10cm
- E084-325** - 11.5cm



DOUGHTY TONGUE PLATE

Slotted


- E084-405** - 6.5cm
- E084-410** - 7.5cm
- E084-415** - 9cm
- E084-420** - 10cm
- E084-425** - 11.5cm



DOYEN MOUTH GAG

With Cam Rack

- E086-105** - Child, 11.5cm
- E086-110** - Adult, 12.5cm



FERGUSSON MOUTH GAG


Ackland Jaw

- E086-205** - Child, 17cm
- E086-210** - Adult, 20cm



KILNER DOTT TONGUE PLATE


- E084-605** - Mouth Gag, Complete Set, Frame, with 3 East Grinstead Tongue Plates, S, M & L
- E084-610** - Small
- E084-615** - Medium
- E084-620** - Large



KILNER DOTT TONGUE PLATE

Insulated

- E084-705** - Mouth Gag, Complete Set, Insulated Frame, with 3 East Grinstead Tongue Plates, S, M & L
- E084-710** - Small
- E084-715** - Medium
- E084-720** - Large



REUSABLE INSTRUMENTATION



Suction Tubes / Trocars

ZOELLNER SUCTION TUBE

Original Pattern

E196-975 - 16cm



FORMBY SUCTION TUBE

E196-275 - Angled

Tapered Tube to
10fg (3.3mm)



LEMPERT SUCTION TUBE

E196-515 - 6fg

E196-520 - 9fg



BELLUCCI SUCTION TUBE

Chrome Finish

E196-205 -

Tube & 3 Ends,
20, 22 & 24 swg

E196-210 -

Tube only, 14cm



MAGILL SUCTION TUBE

Original Pattern

E196-605 - 3mm Dia, S

E196-615 - 4mm Dia, L

E196-610 - 3.5mm Dia, M

E196-620 - 5mm Dia, XL



BARRON SUCTION TUBE

Diameters ranging from 1mm (3fg) - 2.3mm (7fg)

Length Options: 13cm, 17cm, 23cm


***Contact our Customer Service Team for Product Codes**



TILLEY TROCAR

Harpoon

E190-100 - 17cm



TILLEY LICHTWITZ TROCAR

Luer Lock, with Tubing Mount


E190-105 - 5fg

E190-108 - 8fg

E190-106 - 6fg

E190-109 - 9fg

E190-107 - 7fg



REUSABLE INSTRUMENTATION




Supplementary ENT Instrumentation

LEMPERT SCOOP

Oval Scoop, 20cm

Sizes ranging from 2mm x 3.5mm - 3mm x 5.5mm

*Contact our Customer Service Team for Product Codes



LEMPERT RUGINE

Narrow Curved Blade, 3mm

E168-205 - 19cm



WARWICK JAMES ELEVATOR

P050-140 - Right

P050-145 - Left

P050-150 - Straight



CRYERS ELEVATOR

P050-040 - Left

P050-045 - Right



FREER ELEVATOR

E032-305 - Single Ended, 19cm, Sharp

E032-310 - Single Ended, 19cm, Blunt

E032-315 - Double Ended, 18cm, Sharp/Blunt



FREER ELEVATOR

E032-320 - Double Ended, with Knife
Curved to Side, 19cm

E032-350 - Separator & Elevator,
Double Ended, 19cm



REUSABLE INSTRUMENTATION




Supplementary ENT Instrumentation

MITCHELL TRIMMER

Double Ended

P065-101 - 16cm



COLLEGE DRESSING TWEEZERS

15cm

P062-300 - Angled, Serrated Jaw

P062-302 - Angled, Plain Smooth Jaw



TILLEY HENCKEL PUNCH FORCEP


Ethmoid, Graduated Upper Shaft, 11.5cm

E074-405 - 4.5mm Bite

E074-415 - 6.5mm Bite

E074-410 - 5.5mm Bite

E074-420 - 8mm Bite



CITELLI ANTRUM PUNCH FORCEP

7.5cm

E033-102 - 2mm Bite

E033-103 - 3mm Bite

E033-104 - 4mm Bite



LACK TONGUE DEPRESSOR

E024-010 - Child, 13mm Wide

E024-020 - Adult, 22mm Wide


E024-015 - Adolescent, 19mm Wide



CARTRIDGE SYRINGE

Dental

P320-100 - 2.2cc



REUSABLE INSTRUMENTATION




Accessories

MAYO INSTRUMENT PINS

Ball End


G128-090 - 11.5cm



BAG CLIPS

Designed for multiple uses in Theatres & Associated Areas. Manufactured with Robust Spring Construction.


G090-100 - 2.5cm Wide



SLING CLIPS

Designed for holding drapes, dressings, slings, etc. Manufactured with Robust Spring Construction. Serrated Jaw

G090-150 - 6.5cm Wide



AGNEW POLITZER MYRINGOTOME

Angled, Round Handle

E138-105 - 19cm



LARYNGEAL MIRROR

Mirror with Fixed Handle

Sizes ranging from 6mm - 30mm


***Contact our Customer Service Team for Product Codes**



LARYNGEAL MIRROR

E134-205 - Mirror Handle Only

E134-210 - Mirror Handle, Screw Adjustment



LARYNGEAL MIRROR

Mirror Only, No Handle

Sizes ranging from 6mm - 30mm

***Contact our Customer Service Team for Product Codes**



ST CLAIR THOMSON MIRROR

Rhinoscopic

Sizes ranging from 6mm - 16mm

***Contact our Customer Service Team for Product Codes**



Flexible & easy to insert, these **Single-Use Oesophageal Dilating Bougies** are perfect for the **dilation & expansion of the oesophagus and cardia** during Endoscopic procedures.

- Individually **STERILE** Packed.
- Anti-shedding design offering **safety and reliability**.
- **Two radiopaque marks** on both ends for ease of identification under X-ray.
- **Printed markers** on both sides, allowing the clinicians to identify the insertion.
- Length: 850mm | Guide Wire: 0.035".
- Made with **malleable medical polymer** materials.
- Various **sizes ranging from 16fg – 60fg**, to accommodate your requirements.



Product Code	Supply Chain Code	Product Description	Pk Size
BDS-05	DIRECT ONLY	Oesophageal Bougie, 16fg	5
BDS-06	DIRECT ONLY	Oesophageal Bougie, 18/20fg	5
BDS-07	DIRECT ONLY	Oesophageal Bougie, 22fg	5
BDS-08	DIRECT ONLY	Oesophageal Bougie, 24/26fg	5
BDS-09	DIRECT ONLY	Oesophageal Bougie, 28fg	5
BDS-10	DIRECT ONLY	Oesophageal Bougie, 30/32fg	5
BDS-11	DIRECT ONLY	Oesophageal Bougie, 34fg	5
BDS-12	DIRECT ONLY	Oesophageal Bougie, 36/38fg	5
BDS-13	DIRECT ONLY	Oesophageal Bougie, 40fg	5
BDS-14	DIRECT ONLY	Oesophageal Bougie, 42/44fg	5
BDS-15	DIRECT ONLY	Oesophageal Bougie, 46fg	5
BDS-16	DIRECT ONLY	Oesophageal Bougie, 48/50fg	5
BDS-17	DIRECT ONLY	Oesophageal Bougie, 52fg	5
BDS-18	DIRECT ONLY	Oesophageal Bougie, 54/56fg	5
BDS-19	DIRECT ONLY	Oesophageal Bougie, 58fg	5
BDS-20	DIRECT ONLY	Oesophageal Bougie, 60fg	5



TIE WIRES

STERILE


These **Single-Use, STERILE Dental / Oral Tie Wires** are commonly used in **Orthodontic procedures** and are suitable for **long term invasive use** within the oral cavity.

Features & Benefits

- STERILE Packed
- Reduce Costs & Save Time
- Improve Patient Safety
- Remove the Need to Sterilise

Specifications

- Pre-Stretched 10"
- High Quality 316 Annealed Stainless Steel
- Manufactured in the UK



Product Code	Supply Chain Code	Product Description	Pk Size
JMO030S	DIRECT ONLY	0.30mm Diameter	20*
JMO035S	DIRECT ONLY	0.35mm Diameter	20*
JMO040S	DIRECT ONLY	0.40mm Diameter	20*
JMO045S	DIRECT ONLY	0.45mm Diameter	20*
JMO050S	DIRECT ONLY	0.50mm Diameter	20*
JMO055S	DIRECT ONLY	0.55mm Diameter	20*

Pack Size 20* - 2 wires per pouch, 20 pouches per box

SNARE WIRES

STERILE


These **Single-Use, STERILE pre-formed Nasal & Tonsil Snare Wires** are intended to be used to **remove tonsils, polyps and other tissues** in **Tonsillectomy & Nasal** procedures.

Features & Benefits

- STERILE Packed - For ease of use
 - For use with a variety of Nasal & Tonsil snares
- *Contact our Customer Service Team for further information**

Specifications

- Pre-formed
- High Quality 316 hard-drawn Stainless Steel
- 2 wires per pouch, 20 pouches per box



Code: JMN628S | Nasal | 28swg | 0.350mm
Code: JMT725S | Tonsil | 25swg | 0.500mm

HEADLIGHTS & LIGHT SOURCES

LLS-9500 Battery Powered LED Headlight System


CUDA® Surgical's LLS-9500 LED headlight produces an **unparalleled bright, white light with superior battery life** so you can work without interruption. The revolutionary **Rear Cranial Support design** will create a new wearing experience which is **lightweight, comfortable, flexible and durable!**



CUDA® Surgical Battery and Charger:

Each battery lasts **2 hours 30 minutes at full intensity, 5 hours at half power (output)**. The charger will fully charge an empty four cell battery in 3 ½ hours.

The **LLS-9500 system comes with two batteries for easy replacement** when one battery becomes low during procedures and is now available with a single bay or a dual bay charger.



LED Light Source – LED4000

The **LED4000** is a state-of-the-art LED Light Source from **CUDA® Surgical** designed to provide illumination for **Headlights & Endoscopes**.

Expected LED Life of 30,000 hours eliminates the need to purchase replacement lamps, thus saving money on operational costs.



- Compact, efficient, durable source of brilliant “daylight” quality light for headlights & endoscopes.
- LED automatically shuts off when cable is removed from the active port.
- Costly periodic lamp changes eliminated.
- Lumen output comparable to many Xenon light sources but more stable and consistent.
- **DISCOUNTS** available for **Bundle Options**:
(Bundle consists of: Light Source, Headlight, Headlight Cable, Floorstand)

Xenon Light Sources also available

ENDOSCOPIC CABLES

CUDA®
SURGICAL

CUDA® Surgical's range of Endoscopic cables are designed to help deliver optimal light for surgical procedures.

High quality, **versatility**, performance & reliability go hand-in-hand with these **Universal** Fibre-Optic cables.



Compatible with:

- KeyMed®
 - Stryker®
 - Eschmann®
 - Wolf®
 - Dyonics®
 - Karl Storz®
 - Olympus® & GU®
- Light Sources**



BLUE

PINK

RED

GREEN

Features & Benefits:

- Bundle Diameter: 5.0mm
- Cable Length: 7.5ft
- Latex Free Silicone Sheathing
- Fully Autoclavable
- Interchangeable Adaptors to suit your requirements
- Available in different colours for ease of identification

**Other colours & size options
available upon request**

FIBRE LIGHT CABLES & ADAPTORS

Interchangeable Cable Adaptors for use with our range of **CUDA® Surgical Universal Endoscopic Cables**.

A Light Source Adaptor and an Instrument Adaptor are required with each cable. Simply choose your required light source and instrument adaptor and screw to the cable ends.

Cable Adaptors (Light Source End)

Karl Storz®
Code C0100



Wolf®/Dyonics®
Code C0101



ACMI®
Code C0102



Olympus®
Code C0103



Fujinon®
Code C0104



Machida®
Code C0105



Pentax®
Code C0106



ACMI®/Olympus®
Code C0107

Not illustrated

Wolf® Extended
Code C0108

Not illustrated



Cable Adaptors (Instrument End)

**Storz®/
Olympus®**
Code C0200



**Wolf®/Dyonics®
Female**
Code C0201



**ACMI®
Female**
Code C0202



**ACMI®
Male**
Code C0204



**ACMI®
Female Clip**
Code C0205



**Wolf®/Dyonics®
External Clip**
Code EX-0203



**Wolf®/Dyonics®
Male**
Code C0203

Not illustrated

**Olympus®
Retractable**
C0206

Not illustrated

Our range of adaptors provide users with the **versatility** of creating a custom Fibre Light Cable, in a **cost effective** manner, by simply changing the adaptors on each end.











CATALOGUE RESOURCES & PRODUCT VIDEOS

Our **product catalogues** are available for **download directly from our website**, providing you with instant access to our comprehensive range of products at your convenience. Alternatively, if you prefer a **physical copy** you can easily **request one** and we'll ensure it is promptly **delivered** to you.

Explore our extensive selection of products in the format that suits you best!



SCAN TO VIEW
OUR CATALOGUES



SCAN TO VIEW OUR
VIDEO LIBRARY



If you're looking for some informative content, be sure to check out our **YouTube channel**. Our videos cover a wide range of specialities that provide product features & benefits through to useful tips.

Dive into our videos and join our growing community today!



QUALITY & ACCREDITATIONS

Judd Medical is certified to **ISO 13485: 2016 and ISO 14001: 2015** across our two locations with all relevant policies and procedures being robustly implemented.

These certifications collectively highlight our commitment to excellence and our leadership in both quality assurance and environmental sustainability.



ISO 13485: 2016

Demonstrating our commitment to the safety & quality of our medical devices.



ISO 14001: 2015

Demonstrating our commitment to identifying, managing, monitoring and controlling our environmental processes.

LSAS: Level 2

In addition, we have achieved **Labour Standards Assurance System (LSAS) Level 2.**

NOTES



HOW TO REACH US



01527 559010



sales@judd-medical.co.uk



www.judd-medical.co.uk



**Unity House, Buntsford Park Road,
Bromsgrove, Worcestershire, B60 3DX**



Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk



We are an Approved Supplier to the
NHS Supply Chain



Supply Chain