



Judd
medical

Gynae Theatre & Outpatients Catalogue

www.judd-medical.co.uk



NHS

Supply Chain

APPROVED SUPPLIER

TABLE OF CONTENTS

3	Who We Are / How Can We Help
4	What We Do / Training & Education
5	How To Order / Delivery Options
6	Diathermy Pencils (Disposable & Reusable)
7	Diathermy Smoke Evacuation Pencil (Disposable)
8	Bipolar Forceps (Reusable)
9-10	Bipolar Forceps (Disposable)
11-12	Monopolar Forceps (Reusable)
13	Monopolar Forceps (Disposable)
14	Abbey Needles & Laparoscopic Hooks
15	Monopolar & Bipolar Cables (Reusable)
16-18	Diathermy Electrodes (Disposable & Reusable)
19	Colposcopy Loop Electrodes (Disposable)
20	Electrocauteries (Disposable)
21	Diathermy Quivers (Disposable)
22	Diathermy Quivers (Reusable)
23-36	Instrumentation (Reusable)
37-39	Headlights, Light Sources, Cables & Adaptors
40	Catalogue Resources & Product Videos
41	Quality & Accreditations
42	Notes
43	How To Reach Us

Disclaimer: Product specifications in this catalogue are correct at time of publication but may change from time-to-time without notice. It is the customer's responsibility to validate that a particular product with the properties described in the product specification meets the requirements of their application.

WHO WE ARE



Judd Medical Limited is recognised as being amongst one of the most progressive medical companies in the UK and during the last 25 years have seen numerous developments that have enabled us to grow into a leading medical device supplier in the industry.

With the expansion of our product portfolio to include many innovative products and services, it has allowed Judd Medical to make a significant impact on the healthcare market to NHS hospitals, private groups and independent clinics throughout the UK.

HOW CAN WE HELP



CUSTOMER SERVICE LEVELS

We ensure that should you need any help and advice or encounter any problems we have a fully trained team ready to support you. Our specifically trained Customer Service Team offer friendly, knowledgeable advice and support for every requirement.

GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any products proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts, standing and call off orders dependent on your requirements.

JUDD MEDICAL

“YOUR PARTNER IN PRECISION HEALTHCARE”

WHAT WE DO.....



Our company is proud to offer a comprehensive product portfolio of
Single Use & Reusable Surgical Instruments and Medical Devices.

-  Arthrodesis & Kirschner Wires, Twist Drills, Steinmann Pins & Denham Pins
-  Electrosurgery: Diathermy Products – Single Use & Reusable
-  Light Sources: Xenon & LED, Surgical Headlights, Fibre Light Cables & Adaptors
-  Pulse Lavage Systems
-  Replacement Surgical Saw Blades & Burs
-  Single-Use & Reusable Surgical Instrumentation

At Judd Medical we are committed to our partnership with the NHS Supply Chain, NHS Scotland and NHS Wales, ensuring optimal value for our NHS Customers.

We are an approved supplier to the NHS Supply Chain.



TRAINING & EDUCATION



We are dedicated to providing comprehensive training and education on regulatory compliance, organising international seminars and hosting healthcare roundtable events. We facilitate peer-to-peer clinical discussions and offer training on the care and handling of surgical instrumentation and diathermy safety.

Through these initiatives, we aim to enhance professional knowledge, ensure adherence to industry standards and promote the highest levels of safety and efficacy in healthcare practices.

HOW TO ORDER:

You can place an order with us via:



01527 559010



sales@judd-medical.co.uk



NHS Framework Agreements



SCAN TO VIEW OUR
NHS FRAMEWORK
AGREEMENTS

DELIVERY OPTIONS:



Shipping & Despatch – Where stock is available this will be despatched on the same day.
(Order must be processed before 3pm)

Delivery options available:

- Standard Next Day Delivery
- Pre 9am Next Day Delivery*
- Pre 10.30am Next Day Delivery*
- Pre 12 Noon Next Day Delivery*
- Same Day Delivery – Must be processed before 12 noon on date of request.*


*Price on Application

We deliver to:

- Mainland UK
- Northern Ireland
- Channel Islands
- Scottish Highlands & Islands



**Centrally located in Bromsgrove,
Worcestershire**



DIATHERMY PENCILS

DISPOSABLE PENCIL with PUSH-BUTTON CONTROLS

Designed to fulgurate, coagulate, cut or dissect biological tissue during electrosurgery.

- Push-button controls for **Cut & Coagulation** functions
- Available with either a **standard** blade or **non-stick** blade electrode
- Individually Sterile Packed

Disposable Diathermy Pencil - Standard

Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4797	FDJ2684	320mm	2.38mm	Valleylab	Blade F4046	25



Disposable Diathermy Pencil - Non-Stick

Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4797TEF	Direct Only	320cm	2.38mm	Valleylab	Blade F4046/TEF	25

REUSABLE PENCIL with PUSH-BUTTON CONTROLS

These Reusable Diathermy Pencils, with **Push-Button Controls for Cut & Coagulation functions**, are designed for optimal performance during electrosurgery and are complete with **standard blade electrodes**.

REUSABLE (can be sterilised up to **100 times**)




Product Code	Supply Chain Code	Cable Length	Shaft	Connector	Electrode	Pk Size
F4244	FDJ2835	300cm	2.38mm	Valleylab	Blade F4060	1

Disposable Hand Switch Pencil housing 2.38mm (3/32") with cut-coag double button and Smoke Evacuator System.

This Diathermy Pencil with Smoke Evacuation System is designed to **fulgurate, coagulate, cut or dissect biological tissue and also provide the effective removal of surgical smoke** produced during electrosurgery.

- Available with either **2 standard blade or 2 teflon (non-stick) blade electrodes**
- Generator Connector: Universal tripolar moulded **Valleylab**



- 2 x Cannula with screw connections. (Lengths 37mm & 122mm)
- Pencil Length - 170mm
- Material - High resistance plastic, Latex Free
- Corrugated tubing, highly flexible, 300cm, + -5%

Product Code	Supply Chain Code	Product Description	Pk Size
F4797SE	FDJ2771	Smoke Evacuation Pencil (includes 2 Standard Electrodes)	10
F4797SE/TEF	FDJ2736	Smoke Evacuation Pencil (includes 2 Teflon Electrodes)	10

BIPOLAR FORCEPS

Reusable

Precision Reusable Bipolar Forceps designed for grasping & coagulation of tissue during electrosurgery.



- Types available: **Adson, Bayonet, Bayonet (angled down), Jeweller#5, McPherson & Straight options**
- Tip size options: **0.5mm or 1.0mm**

50 times Autoclavable



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-103/R1.0	Direct Only	Straight Forcep	12cm	1.0mm	1
PH786-105/R1.0	FDJ2209	Straight Forcep	15cm	1.0mm	1
PH786-109/R1.0	Direct Only	Straight Forcep	18cm	1.0mm	1
PH786-111/R1.0	Direct Only	Straight Forcep	20cm	1.0mm	1
PH786-113/R1.0	Direct Only	Straight Forcep	22cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-127/R1.0	Direct Only	Bayonet Forcep	18cm	1.0mm	1
PH786-130/R1.0	Direct Only	Bayonet Forcep	20cm	1.0mm	1
PH786-131/R1.0	Direct Only	Bayonet Forcep (angled down)	20cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-141/R1.0	Direct Only	Adson Forcep	12cm	1.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-151/R0.5	Direct Only	Jeweller Forcep #5 (straight point)	11.5cm	0.5mm	1
PH786-155/R0.5	Direct Only	McPherson Forcep	10.2cm	0.5mm	1

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

BIPOLAR FORCEPS

(with Attached Cable)

Disposable

STANDARD

- Individually sterile packaged & boxed in 25's
- Tip Sizes: **0.5mm & 1.0mm**
- Various Lengths
- Flexible bridge connector with 2 x 4mm plugs
- Suitable for the most popular diathermy machines including:
Valleylab, Conmed, ERBE & more
- Pre-attached 3M cable
- All STERILE bipolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-105/DWC0.5	FDJ2195	Straight Forcep	15cm	0.5mm	STERILE	25
PH786-105/DWC1.0	FDJ2196	Straight Forcep	15cm	1.0mm	STERILE	25
PH786-107/DWC0.5	FDJ2680	Straight Forcep	16cm	0.5mm	STERILE	25
PH786-107/DWC1.0	FDJ2724	Straight Forcep	16cm	1.0mm	STERILE	25
PH786-109/DWC1.0	FDJ2197	Straight Forcep	18cm	1.0mm	STERILE	25
PH786-111/DWC1.0	FDJ2198	Straight Forcep	20cm	1.0mm	STERILE	25
PH786-112/DWC1.0	FDJ2199	Straight Forcep (angled tip)	20cm	1.0mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-124/DWC1.0	FDJ2200	Bayonet Forcep	16cm	1.0mm	STERILE	25
PH786-130/DWC1.0	FDJ2201	Bayonet Forcep	20cm	1.0mm	STERILE	25
PH786-133/DWC1.0	FDJ2202	Bayonet Forcep	22cm	1.0mm	STERILE	25
PH786-136/DWC1.0	FDJ2203	Bayonet Forcep	24cm	1.0mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-141/DWC1.0	FDJ2194	Adson Forcep	12cm	1.0mm	STERILE	25
PH786-151/DWC0.5	FDJ2193	Jeweller Forcep	11.5cm	0.5mm	STERILE	25
PH786-155/DWC0.5	FDJ2192	McPherson Forcep	10.2cm	0.5mm	STERILE	25

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

BIPOLAR FORCEPS

(with Attached Cable)

Disposable

NON-STICK

- Individually sterile packaged & boxed in 25's
- Tip Sizes: **0.5mm & 1.0mm**
- Various Lengths
- Gold Coloured Tips
- Flexible bridge connector with 2 x 4mm plugs
- Suitable for the most popular diathermy machines including:

Valleylab, Conmed, ERBE & more

- Pre-attached 3M cable
- All STERILE bipolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue



Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-105/NS0.5	FDJ2733	Non-Stick Straight Bipolar Forcep	15cm	0.5mm	STERILE	25
PH786-105/NS1.0	FDJ2754	Non-Stick Straight Bipolar Forcep	15cm	1mm	STERILE	25
PH786-111/NS1.0	FDJ2700	Non-Stick Straight Bipolar Forcep	20cm	1mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-124/NS1.0	FDJ2755	Non-Stick Bayonet Bipolar Forcep	16cm	1mm	STERILE	25
PH786-130/NS1.0	FDJ2653	Non-Stick Bayonet Bipolar Forcep	20cm	1mm	STERILE	25

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-141/NS1.0	FDJ2751	Non-Stick Adson Bipolar Forcep	12cm	1mm	STERILE	25
PH786-151/NS0.5	FDJ2648	Non-Stick Jeweller Bipolar Forcep	11.5cm	0.5mm	STERILE	25

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

MONOPOLAR FORCEPS

Reusable

Reusable Monopolar Forceps for grasping & coagulation of tissue during electrosurgical diathermy procedures.



- Types available: **Adson, Gillie (Serrated), McIndoe, Waugh, Gillies 1:2 Toothed, DeBakey & Straight options**
- Tip Size: **2.0mm**
- Various Lengths available



McIndoe Forcep



Straight Forcep



Adson Forcep



Gillie (Serrated) Forcep

50 times Autoclavable

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-301/R2.0	FDJ2666	Straight Standard Forcep	12cm	2.0mm	1
PH786-303/R2.0	FDJ2738	Straight Standard Forcep	15cm	2.0mm	1
PH786-305/R2.0	FDJ2732	Straight Standard Forcep	18cm	2.0mm	1
PH786-307/R2.0	FDJ2670	Straight Standard Forcep	20cm	2.0mm	1
PH786-309/R2.0	FDJ2656	Straight Standard Forcep	22cm	2.0mm	1
PH786-311/R2.0	FDJ2725	Straight Standard Forcep	25cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-313/R2.0	FDJ2645	McIndoe Forcep	15cm	2.0mm	1
PH786-315/R2.0	FDJ2762	Gillie (Serrated) Forcep	15cm	2.0mm	1
PH786-317/R2.0	FDJ2675	Adson Forcep	12cm	2.0mm	1
PH786-319/R2.0	FDJ2187	Adson Forcep	15cm	2.0mm	1


Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

MONOPOLAR FORCEPS



Reusable

- Tip Size: **2.0mm**
- Various Lengths available



Waugh Forcep



De Bakey Forcep



Gillies 1:2 Forcep

50 times Autoclavable

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-321/R2.0	FDJ2756	Waugh Forcep	15cm	2.0mm	1
PH786-323/R2.0	FDJ2763	Waugh Forcep	18cm	2.0mm	1
PH786-325/R2.0	FDJ2721	Waugh Forcep	20cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-335/R2.0	FDJ2635	Gillies 1:2 Toothed Forcep	15cm	2.0mm	1
PH786-337/R2.0	FDJ2765	Gillies 1:2 Toothed Forcep	18cm	2.0mm	1
PH786-339/R2.0	FDJ2746	Gillies 1:2 Toothed Forcep	20cm	2.0mm	1

Product Code	Supply Chain Code	Product Description	Size	Tip Size	Pk Size
PH786-345/R2.0	FDJ2695	De Bakey Forcep	15cm	2.0mm	1
PH786-347/R2.0	FDJ2848	De Bakey Forcep	18cm	2.0mm	1
PH786-349/R2.0	FDJ2649	De Bakey Forcep	20cm	2.0mm	1

Listings shown do not reflect the entire range of Forceps available, additional sizes also available.


MONOPOLAR FORCEPS

(with Attached Cable)

Disposable

- Suitable for the most popular diathermy machines including:
 - **Valleylab**
 - **Conmed**
 - **ERBE & more**
- Pre-attached 3M cable with **8mm generator pin**
- All STERILE monopolar forceps are supplied **DOUBLE PACKED**
- Lightweight to prevent fatigue

Straight
Forcep



McIndoe
Forcep



1:2 Teeth
Gillies Forcep



Waugh
Forcep





Product Code	Supply Chain Code	Product Description	Size	Tip Size	Single Use	Pk Size
PH786-305/DWC2.0	FDJ2191	Serrated Straight Dissecting Monopolar Forcep	18cm	2.0mm	STERILE	25
PH786-313/DWC2.0	FDJ2189	Serrated McIndoe Monopolar Forcep	15cm	2.0mm	STERILE	25
PH786-339/DWC2.0	FDJ2184	1:2 Teeth Gillies Monopolar Forcep	20cm	2.0mm	STERILE	25
PH786-523/DWC2.0	FDJ2190	Serrated Waugh Monopolar Forcep	18cm	2.0mm	STERILE	25



Listings shown do not reflect the entire range of Forceps available, additional sizes also available.

ABBAY NEEDLES

Disposable



Product Code	Supply Chain Code	Product Description	Size	Pk Size
PH786-389/DWC	FDJ2786	Abbay Needle	17.5cm	25

LAPAROSCOPIC HOOKS

Disposable

These **Non-Stick (PTFE) coated** Laparoscopic Hooks resist eschar build-up which is a common problem during electrosurgery.

Features & Benefits:

- Minimises Eschar Accumulation
- Minimises Risk
- Latex Free
- 8mm Pin Connector
- Pre-attached - 3M Cable
- Sold in Packs of 10

Laparoscopic Hook Electrodes are:

- Sterilised by Ethylene Oxide
- Individually double pouched in tough Tyvek® backed polyethylene pouches



Product Code	Product Description	Length	Pin Connector	Pk Size
PH786-660/DWC8mm	Laparoscopic Hook	330mm	8mm	10

MONOPOLAR CABLES




Reusable

These Monopolar Cables are compatible with most standard electrosurgical devices and compatible with the following generators and forceps:

Valleylab, Conmed, Megadyne, Bovie, Medtronic, Dremed and all others that are similar to the aforementioned brands.

50 times Autoclavable

Product Code	Product Description	No. of Uses	Machine Fitting	Instrument Fitting	Pk Size
PHF685-3M	Monopolar Cable 4mm pin to 4mm	50			1
PHF690-3M	Monopolar Cable 4mm pin to 4.8mm	50			1
PHF675-3M	Monopolar Cable 8mm pin to 4mm	50			1
PHF680-3M	Monopolar Cable 8mm pin to 4.8mm	50			1

BIPOLAR CABLES



Reusable

Product Code	Product Description	No. of Uses	Machine Fitting	Instrument Fitting	Pk Size
OBS581F	Reusable Bipolar Cable Blue with 4mm Pins & Flexible Bridge Connector European Block End Connector 3M Length	50			1
OBS580TF	Reusable Bipolar Cable Blue with 4mm Pins & Flexible Bridge Connector US Two Pin Forceps Connector 3M Length	50			1

Universal Bipolar Cables available upon request.

ELECTRODES

Disposable

The **Teflon (PTFE) coating** makes these electrodes “Non-Stick” and is only present on the active part of the electrode. Suitable for general diathermy.

STAINLESS STEEL - Sterile

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044	FDJ2652	70mm	BALL 4mm	24
F4046	FDJ2715	70mm	BLADE	10
F4046/S	Direct Only	50mm (44mm sleeve)	BLADE	24
F4048	FDJ2661	70mm	NEEDLE	10
F4048/S	Direct Only	69mm (63mm sleeve)	NEEDLE	24
F4049	FDJ2636	54mm	BALL 2mm	24
F4050	Direct Only	153mm	BLADE	24
F4050/S	Direct Only	153mm (147mm sleeve)	BLADE	24
F4052	Direct Only	152mm	NEEDLE	24
F4052/S	Direct Only	152mm (146mm sleeve)	NEEDLE	24
F4053	Direct Only	100mm	BLADE	24
F4067	FDJ2787	156mm	BALL 4mm	24
F4067/S	FDJ2792	157mm	BALL 5mm	24



TEFLON COATED - Sterile - NON STICK

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044/TEF	FDJ2836	70mm	BALL 4mm	24
F4046/TEF	FDJ2647	70mm	BLADE	10
F4046S/TEF	FDJ2658	70mm (65mm sleeve)	BLADE	10
F4048/TEF	FDJ2663	70mm	NEEDLE	10
F4048S/TEF	FDJ2729	70mm (65mm sleeve)	NEEDLE	10
F4050/TEF	FDJ2757	150mm	BLADE	10
F4050S/TEF	FDJ2639	150mm (145mm sleeve)	BLADE	10
F4052/TEF	Direct Only	150mm	NEEDLE	24
F4053/TEF	FDJ2846	100mm	BLADE	10
F4053S/TEF	FDJ2791	100mm (95mm sleeve)	BLADE	10
F4067/TEF	FDJ2828	150mm	BALL 4mm	10
OBS0312C64	Direct Only	64mm (shielded tip)	BLADE	10
OBS0324C120	Direct Only	120mm	BALL 3mm	10
OBS0323C120	Direct Only	120mm	BALL 5mm	10
OBS0314C130	Direct Only	130mm	BLADE	10




These devices are Latex Free.

Images shown do not reflect the size or entire range of Electrodes available.

The **Titanium coating** makes these electrodes “Non-Stick” and is only present on the active part of the electrode.

TITANIUM COATED - Sterile - NON STICK

Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4044T	FDJ2651	70mm	BALL 4mm	10
F4046T	Direct Only	70mm	BLADE	24
F4048T	Direct Only	70mm	NEEDLE	24
F4050T	FDJ2671	153mm	BLADE	24
F4052T	Direct Only	152mm	NEEDLE	24
F4067T	FDJ2840	156mm	BALL 4mm	24




Disposable **Micro Dissection Needles** suitable for a wide variety of diathermy procedures and allowing for greater accuracy. The **TUNGSTEN Dissection Needles** offer **precise tissue dissection** with improved cutting & coagulation performance.

MICRO DISSECTION NEEDLES - Sterile



Product Code	Supply Chain Code	Product Description	Shape	Pk Size
F4051/TU108	FDJ2837	108mm (104mm sleeve)	NEEDLE	10
F4051/TU75	Direct Only	75mm (71mm sleeve)	NEEDLE	24
F4048/TU	Direct Only	65mm (61mm sleeve)	NEEDLE	24
F4051/TU	FDJ2808	55mm (51mm sleeve)	NEEDLE	10
F4051/45TU	FDJ2854	55mm (51mm sleeve)	NEEDLE 45° Angled	10
F4042/TU	FDJ2676	55mm (51mm sleeve)	NEEDLE 30° Angled	10
F4047/TU	Direct Only	45mm (41mm sleeve)	NEEDLE	24



Images shown do not reflect the size or entire range of Electrodes available.

ELECTRODES

Reusable

STAINLESS STEEL 2.38MM SHAFT

Product Code	Supply Chain Code	Length	Shape	Pk Size
F4045	FDJ2816	156mm	SPHERE 4mm	5
F4045/10	FDJ2655	156mm	SPHERE 4mm, 10° angled	5
F4045/20	FDJ2735	156mm	SPHERE 4mm, 20° angled	5
F4045/30	FDJ2654	156mm	SPHERE 4mm, 30° angled	5
F4045/40	FDJ2728	156mm	SPHERE 4mm, 40° angled	5
F4053R	FDJ2822	100mm	BLADE	5
F4060	FDJ2826	69mm	BLADE	5
F4061	FDJ2820	100mm	ROUND NEEDLE, 20° angled	5
F4062	FDJ2827	150mm	CURVED BLADE	5
F4063	FDJ2779	152mm	NEEDLE 0.78mm	5
F4064	FDJ2857	150mm	SPHERE 2.5mm	5
F4065	FDJ2782	153mm	BLADE	5
F4066	FDJ2781	152mm	NEEDLE 1.30mm	5
F4068	FDJ2780	70mm	SPHERE 4mm	5
F4084	FDJ2803	69mm	NEEDLE 0.75mm	5
F4086	FDJ2818	69mm	NEEDLE 1.30mm	5



STAINLESS STEEL 4MM SHAFT

Product Code	Supply Chain Code	Length	Shape	Pk Size
F4064/4	Direct Only	55mm (51mm sleeve)	BLADE	5
F4065/4	FDJ2793	160mm	BLADE	5
F4068/4	FDJ2847	50mm	SPHERE 4mm	5
F4084/4	FDJ2807	55mm	NEEDLE 0.8mm	5


EXTENSIONS

Product Code	Length	Pk Size
EXTENSIONS DISPOSABLE		
F4902	102mm	5
F4903	150mm	5
F4904	182mm	5
EXTENSIONS REUSABLE		
F4952	102mm	5
F4953	150mm	5
F4954	182mm	5



Product Code	Supply Chain Code	Product Description	Pk Size
OB0353C	Direct Only	Colposcopy Loop Electrode, 20mm x 12mm	25
OB0354C	Direct Only	Colposcopy Loop Electrode, 20mm x 15mm	25
OBS0344	Direct Only	Colposcopy Loop Electrode, 20mm x 20mm	25
F4620	FDJ2823	Disposable Round Loop Electrode 10mm x 3mm	24
F4621	FDJ2817	Disposable Round Loop Electrode 10mm x 5mm	24
F4622	FDJ2804	Disposable Round Loop Electrode 10mm x 7mm	24
F4623	FDJ2166	Disposable Round Loop Electrode 10mm x 10mm	24
F4624	FDJ2707	Disposable Round Loop Electrode 15mm x 10mm	24
F4625	FDJ2168	Disposable Round Loop Electrode 20mm x 10mm	24
F4626	FDJ2167	Disposable Round Loop Electrode 15mm x 12mm	24
F4630	FDJ2712	Disposable Square Loop Electrode 10mm x 4mm	24
F4631	FDJ2170	Disposable Square Loop Electrode 10mm x 8mm	24
F4632	Direct Only	Disposable Square Loop Electrode 10mm x 10mm	24
F4640	Direct Only	Disposable Trapezoidal Loop Electrode 30mm x 5mm	24
F4641	Direct Only	Disposable Reversed Curved Loop Electrode 23mm x 13mm	24
F4642	Direct Only	Disposable Triangular Loop Electrode 17mm x 12mm	24

Individually STERILE Packaged



FIAB Sterile Cautery pens are used in haemostasis, granting a **fast and precise cauterisation of small blood vessels**. They are extremely **easy to handle**, thanks to an **internal long life battery**; several models are available, different from each other for **working temperatures (1200c, 800c and 600c)**, **shape, thickness and length** of the cautery tip.

The available range allows the user to select a cautery pen that perfectly matches the clinical demands of the procedure.



Supplied in Packs of 10



Product Code	Supply Chain Code	Tip Shape	Tip Length (mm)	Temp at activation (°C/°F)	Power Supply (V)	Diameter (mm)	Weight (GR)	Series
F7244	FDJ2694	FINE	28	1200/2200	3,0	16	65	HTC*
F7266	FDJ2770	LARGE	28	1200/2200	3,0	16	65	HTC*
F7234	FDJ2679	FINE	75	1200/2200	3,0	16	78	HTC*
F7277	FDJ2764	FINE	125	1200/2200	3,0	16	78	HTC*
F7288	FDJ2683	THICK	28	800/1400	3,0	16	65	MTC*
F7255	FDJ2711	FINE	28	600/1100	1,5	16	45	LTC*
F7255UF	FDJ2772	ULTRA FINE	28	600/1100	1,5	16	45	LTC*

HTC* High Temperature Cautery
MTC* Med. Temperature Cautery
LTC* Low Temperature Cautery



DIATHERMY QUIVER with Laparoscopic Extension

This Disposable Diathermy Quiver is a **cost effective** means of offering protection to patients and medical staff by **safely securing instrumentation during surgery**.

An additional benefit of this product is the availability of a **Laparoscopic Extension** which is designed to be used in conjunction with the aforementioned quiver. This open ended extension fits into the quiver providing the end user with an **extended overall working length of 400mm**.



Features & Benefits:

- Single Use, STERILE Packaged
- Lightweight, Strong & Durable
- Integral Towel Clip ring
- Clear in colour
- Double Wrapped
- COMPETITIVELY PRICED
- Ease of accessibility



Product Code	Supply Chain Code	Product Description	Pk Size
AVM4047	Direct Only	Single Use Diathermy Quiver - 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	65
AVM4085	Direct Only	Single Use Laparoscopic Extension - 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	40
AVM4086	Direct Only	Laparoscopic Quiver Combined - 20cm + 20cm Material: Polypropylene, Colour: Clear Sterile - Double wrapped	40

With **4 different options available**, these **Reusable Diathermy Quivers** provide a secure and hygienic method of storing diathermy instruments while not in use during surgical procedures.

Options Available:

- **Standard Diathermy Quiver**
- **Diathermy Quiver with Removable Screw Bottom**
- **Julian Taylor Diathermy Quiver with Drainage Holes**



Standard Diathermy Quiver



Diathermy Quiver with Removable Screw Bottom



Julian Taylor Diathermy Quiver with Drainage Holes



Features & Benefits:

- REUSABLE & SUSTAINABLE
- Robust & Long Lasting
- Easy to Clean
- Available in different sizes
- UK Manufactured
- Material: Polypropylene
- DISCOUNTS available for quantity purchases
- COMPETITIVELY PRICED, excellent value for money

Product Code	Product Description	Size	Pk Size
G022-987	Diathermy Quiver with Removable Screw Bottom	D 6cm x L 20cm	1
G022-975		D 6cm x L 24cm	1
G022-985		D 6cm x L 37.5cm	1
G022-905	Standard Diathermy Quiver	D 4.8cm x L 24.5cm	1
G022-910		D 6cm x L 23cm	1
G022-950		D 6cm x L 37.5cm	1
G022-915	Julian Taylor Diathermy Quiver with Drainage Holes	D 4.4cm x L 25cm	1

REUSABLE INSTRUMENTATION



Cannula & Catheters

Leech Wilkinson Uterine Cannula

F022-123 (Small)

F022-125 (Medium)

F022-127 (Large)



Spackmann Uterine Cannula with Permanent Metal Acorn

F022-144 (14fg)

F022-140 (18fg)



Female Catheter Round

Diameters ranging from

2mm - 6mm

*Contact our Customer Service Team for Product Codes



Rotunda Catheter 20cm

F022-309 (8eg)

F022-310 (10eg)



REUSABLE INSTRUMENTATION




Hysterectomy Clamps

Bonney Berkley Hysterectomy Clamp

1:2 Teeth, Curved


F034-280 (19cm)



Bonney Berkley Hysterectomy Clamp

1:2 Teeth, Angled


F034-290 (19cm)



Faure Hysterectomy Clamp

1:2 Teeth, Longitudinal Serrations

F034-317 (20cm)



Gwilliam Hysterectomy Clamp

1:2 Teeth, Curved

F034-325 (20cm)



REUSABLE INSTRUMENTATION




Hysterectomy Clamps

Gwilliam Hysterectomy Clamp

1:2 Teeth, Straight

F034-320 (20cm)




Heaney Hysterectomy Forcep

Curved

F034-445 (Double Toothed Jaw, 20cm)

F034-440 (Single Toothed Jaw, 20cm)




Maingot Hysterectomy Forcep

1:2 Teeth

F034-344 (Curved, 20cm)

F034-342 (Straight, 20cm)




DeBakey (Rogers) Hysterectomy Forcep

Atraumatic Jaw

Options:

- Full Curve, Medium Curve, Light Curve, Straight
- 22cm or 24cm

*Contact our Customer Service Team for Product Codes



REUSABLE INSTRUMENTATION



Curettes

Sims Uterine Curette (UK Pattern)

Single Ended, Blunt

F026-305 (6mm, Small, 25cm)

F026-310 (10mm, Medium, 25cm)

F026-315 (13mm, Large, 25cm)



Sims Uterine Curette (UK Pattern)

Single Ended, Sharp

F026-320 (6mm, Small, 25cm)

F026-325 (10mm, Medium, 25cm)

F026-330 (13mm, Large, 25cm)



Sims Uterine Curette (UK Pattern)

Double Ended, Sharp/Blunt Combined

F026-503 (4.5mm Wide, Extra Small, 25cm)

F026-505 (6mm Wide, Small, 25cm)

F026-510 (10mm Wide, Medium, 25cm)

F026-515 (13mm Wide, Large, 25cm)



Recamier Uterine Curette

Rigid Shaft, Sharp, 31cm

Sizes ranging from 5mm - 40mm

*Contact our Customer Service Team for Product Codes



Recamier Uterine Curette

Rigid Shaft, Blunt, 31cm

Sizes ranging from 5mm - 20mm

*Contact our Customer Service Team for Product Codes



REUSABLE INSTRUMENTATION




Dilators

Bonney Barker Dilator

S/S

Sizes ranging from 3/5mm - 18/20mm

*Contact our Customer Service Team for Product Codes



Hawkin Ambler Dilator

Single Ended

Sizes ranging from 3/6mm - 18/21mm

*Contact our Customer Service Team for Product Codes



Hegar Uterine Dilator

S/S, Double Ended

Sizes ranging from 1/2mm - 25/26mm

*Contact our Customer Service Team for Product Codes




Hegar Uterine Dilator

S/E, Stainless Steel

Sizes ranging from 1mm - 26mm

*Contact our Customer Service Team for Product Codes



REUSABLE INSTRUMENTATION




Dressing Forceps & Myoma Screws

Bozemann Uterine

Dressing Forcep

Grooved Jaw, Curved

F033-326 (26.5cm)




Bozemann Uterine

Dressing Forcep

Serrated Jaw, Curved


F023-328 (26.5cm)



Doyen Myoma Screw

Abdominal

F047-115 (15cm)



Doyen Myoma Screw

Vaginal

F047-120 (18cm)



REUSABLE INSTRUMENTATION



Obstetric Forceps

**Anderson (Simpson-Braun)
Obstetric Forcep**

F036-074 (35.5cm)



**Barnes Neville
Obstetric Forcep**

F036-116 (35.5cm)



**Kielland
Obstetric Forcep**

F036-080 (41cm)



**Simpson
Obstetric Forcep**

F036-120 (30cm)



**Wrigley
Obstetric Forcep**

F036-090 (25cm)



REUSABLE INSTRUMENTATION



Ovum & Peritoneum Forceps


McClintock Ovum Forcep

F038-727 (24cm)



Greenhalgh Ovum Forcep

F038-624 (24cm)




Heywood Smith Ovum Forcep

F038-726 (26.5cm)



Goodwin Amniotomy Forcep Double Curve

F031-424 (24cm)



Green Armitage Forcep Screw Joint

F032-220 (Straight, 20cm)
F032-225 (Angled, 20cm)



REUSABLE INSTRUMENTATION



Punch Forceps

Eppendorf

Cervical Biopsy Forcep

4.5mm x 8mm Round Basket Jaw

F032-405 (20cm Effective Length)




Leech Wilkinson

Cervical Biopsy Forcep

Round Basket Jaw

F032-420 (20cm Effective Length)



Schumacher

Cervical Biopsy Forcep

4mm Hooked Triangular Jaw

F032-430 (24cm Effective Length)



Wittner

Cervical Biopsy Forcep

4.5mm x 8mm Toothed Round Basket Jaw

F032-730 (Straight, 20cm)

F032-740 (Angled, 20cm)



REUSABLE INSTRUMENTATION



Retractors

Doyen Pozzi Retractor

F134-513 (Small, 4.8cm x 6cm Blade)

F134-515 (Medium, 4.8cm x 9cm Blade)

F134-517 (Large, 4.8cm x 12.5cm Blade)



Doyen Retractor

Flat Blade

F134-500 (3cm x 10cm Blade)



Landon Vaginal Retractor

2.5cm x 9cm Blade

F134-260 (19.5cm)



Wertheim Vaginal Retractor

F134-800 (3.8cm x 10cm Blade)



Breisky (Vienna) Vaginal Retractor

Blade Sizes ranging from
2.0cm x 13cm - 4.0cm x 16cm

*Contact our Customer Service Team for Product Codes



REUSABLE INSTRUMENTATION



Scissors

Barnes Episiotomy Scissor

Sharp/Blunt, Angled

F046-075 (12.5cm)



Braun Stadler Episiotomy Scissor

Angled

G140-065 (15cm)

G140-067 (22cm)



Sims Uterine Scissor

Blunt/Blunt

F046-310 (Straight, 20cm)

F046-312 (Straight, 23cm)

F046-322 (Curved, 20cm)

F046-324 (Curved, 23cm)



USA Pattern Umbilical Scissor

Heavy Curve

G140-910 (10cm)



Vant Episiotomy Scissor

Blunt/Blunt, Straight

F046-124 (24cm)



REUSABLE INSTRUMENTATION



Sounds & Vaginal Speculum

Sims Uterine Sound

Malleable, Graduated in cms

F048-042 (30cm)



Horrocks Uterine Sound

Graduated in cms, Silver Plated

F048-038 (29cm)



Cusco Speculum

UK Pattern, Centre Screw

Blade sizes ranging from:
17mm x 70mm - 35mm x 95mm

*Contact our Customer Service Team for Product Codes



Sims Speculum

Double Ended,
Sizes: Small, Medium, Large

Sims Vaginal Speculum

Double Ended, Insulated
Sizes: Small, Medium, Large

*Contact our Customer Service Team for Product Codes



Grave (Brewer) Vaginal Speculum

F050-095 (Small, 2cm x 7.5cm Blade)

F050-100 (Medium, 3cm x 9.5cm Blade)

F050-105 (Large, 3.5cm x 11.5cm Blade)



REUSABLE INSTRUMENTATION



Vaginal Speculum

Sims Vaginal Speculum

Single Ended

F050-410 (Small, 2cm x 7.5cm Blade)

F050-415 (Medium, 2.8cm x 8.5cm Blade)

F050-420 (Large, 4cm x 9cm Blade)



Auvard Vaginal Speculum

Stainless Steel, Fixed Weight, Seamless

F050-055 (Blade 6.5cm long, 4.5cm wide tapering to 4cm)



Kristeller Vaginal Speculum only

Standard Pattern

F051-213 (70mm x 27mm)

F051-215 (75mm x 30mm)

F051-217 (90mm x 36mm)



Kristeller Vaginal Speculum only

Long Bladed Pattern

F051-311 (110mm x 18mm)

F051-313 (110mm x 21mm)

F051-315 (110mm x 27mm)

F051-317 (110mm x 30mm)

F051-319 (110mm x 36mm)



REUSABLE INSTRUMENTATION



Vulsellum Forceps

Bonney Uterine Polypus Forcep

Straight

F040-020 (24cm)



Vant Uterine Polypus Forcep

Curved


F040-028 (24cm)



Lewer Vulsellum Forcep

2:2 Teeth


F040-045 (23cm)



Luer Vulsellum Forcep

1:1 Teeth

F040-050 (23cm)



Schroeder Vulsellum Forcep

1:1 Teeth

F040-100 (25cm)



Teale Vulsellum Forcep

Curved, 3:4 Teeth

F040-118 (18cm)

F040-123 (23cm)



Museux Vulsellum Forcep

2:2 Teeth, 24cm

Straight or Curved options

Widths ranging from 6mm - 10mm

*Contact our Customer Service Team for Product Codes



HEADLIGHTS & LIGHT SOURCES

LLS-9500 Battery Powered LED Headlight System


CUDA® Surgical's LLS-9500 LED headlight produces an **unparalleled bright, white light with superior battery life** so you can work without interruption. The revolutionary **Rear Cranial Support design** will create a new wearing experience which is **lightweight, comfortable, flexible and durable!**



CUDA® Surgical Battery and Charger:

Each battery lasts **2 hours 30 minutes at full intensity, 5 hours at half power (output)**. The charger will fully charge an empty four cell battery in 3 ½ hours.

The **LLS-9500 system comes with two batteries for easy replacement** when one battery becomes low during procedures and is now available with a single bay or a dual bay charger.



LED Light Source – LED4000

The **LED4000** is a state-of-the-art LED Light Source from **CUDA® Surgical** designed to provide illumination for **Headlights & Endoscopes**.

Expected LED Life of 30,000 hours eliminates the need to purchase replacement lamps, thus saving money on operational costs.



- Compact, efficient, durable source of brilliant “daylight” quality light for headlights & endoscopes.
- LED automatically shuts off when cable is removed from the active port.
- Costly periodic lamp changes eliminated.
- Lumen output comparable to many Xenon light sources but more stable and consistent.
- **DISCOUNTS** available for **Bundle Options**:
(Bundle consists of: Light Source, Headlight, Headlight Cable, Floorstand)

Xenon Light Sources also available

ENDOSCOPIC CABLES

CUDA®
SURGICAL

CUDA® Surgical's range of Endoscopic cables are designed to help deliver optimal light for surgical procedures.

High quality, **versatility**, performance & reliability go hand-in-hand with these **Universal** Fibre-Optic cables.



Compatible with:

- KeyMed®
 - Stryker®
 - Eschmann®
 - Wolf®
 - Dyonics®
 - Karl Storz®
 - Olympus® & GU®
- Light Sources**



BLUE

PINK

RED

GREEN

Features & Benefits:

- Bundle Diameter: 5.0mm
- Cable Length: 7.5ft
- Latex Free Silicone Sheathing
- Fully Autoclavable
- Interchangeable Adaptors to suit your requirements
- Available in different colours for ease of identification

**Other colours & size options
available upon request**

FIBRE LIGHT CABLES & ADAPTORS

Interchangeable Cable Adaptors for use with our range of **CUDA® Surgical Universal Endoscopic Cables**.

A Light Source Adaptor and an Instrument Adaptor are required with each cable. Simply choose your required light source and instrument adaptor and screw to the cable ends.

Cable Adaptors (Light Source End)

Karl Storz®
Code C0100



Wolf®/Dyonics®
Code C0101



ACMI®
Code C0102



Olympus®
Code C0103



Fujinon®
Code C0104



Machida®
Code C0105



Pentax®
Code C0106



ACMI®/Olympus®
Code C0107

Not illustrated

Wolf® Extended
Code C0108

Not illustrated



Cable Adaptors (Instrument End)

**Storz®/
Olympus®**
Code C0200



**Wolf®/Dyonics®
Female**
Code C0201



**ACMI®
Female**
Code C0202



**ACMI®
Male**
Code C0204



**ACMI®
Female Clip**
Code C0205



**Wolf®/Dyonics®
External Clip**
Code EX-0203



**Wolf®/Dyonics®
Male**
Code C0203

Not illustrated

**Olympus®
Retractable**
C0206

Not illustrated

Our range of adaptors provide users with the **versatility** of creating a custom Fibre Light Cable, in a **cost effective** manner, by simply changing the adaptors on each end.



CATALOGUE RESOURCES & PRODUCT VIDEOS

Our **product catalogues** are available for **download directly from our website**, providing you with instant access to our comprehensive range of products at your convenience. Alternatively, if you prefer a **physical copy** you can easily **request one** and we'll ensure it is promptly **delivered** to you.

Explore our extensive selection of products in the format that suits you best!



SCAN TO VIEW
OUR CATALOGUES



SCAN TO VIEW OUR
VIDEO LIBRARY



If you're looking for some informative content, be sure to check out our **YouTube channel**. Our videos cover a wide range of specialities that provide product features & benefits through to useful tips.

Dive into our videos and join our growing community today!



QUALITY & ACCREDITATIONS

Judd Medical is certified to **ISO 13485: 2016 and ISO 14001: 2015** across our two locations with all relevant policies and procedures being robustly implemented.

These certifications collectively highlight our commitment to excellence and our leadership in both quality assurance and environmental sustainability.



ISO 13485: 2016

Demonstrating our commitment to the safety & quality of our medical devices.



ISO 14001: 2015

Demonstrating our commitment to identifying, managing, monitoring and controlling our environmental processes.

LSAS: Level 2

In addition, we have achieved **Labour Standards Assurance System (LSAS) Level 2.**

NOTES



HOW TO REACH US



01527 559010



sales@judd-medical.co.uk



www.judd-medical.co.uk



**Unity House, Buntsford Park Road,
Bromsgrove, Worcestershire, B60 3DX**



Judd Medical Limited
Unity House
Buntsford Park Road
Bromsgrove
Worcestershire, B60 3DX

01527 559010
sales@judd-medical.co.uk
www.judd-medical.co.uk



We are an Approved Supplier to the
NHS Supply Chain



Supply Chain